

Plants Research

Tackling key challenges in the life sciences
with advanced biological tools



Introduction

Plants offer a cost-effective means to produce therapeutically useful monoclonal antibodies. However, antibodies produced in plants differ from those produced in mammalian cells. Notably, the glycan composition of plant antibodies can affect their binding activity.

“One in ten plant species contains anticancer substances of variable potency, but relatively few have been bioassayed.”

Edward O. Wilson (born 1929)

Contents

Antibody.....	1
Biochemicals.....	3
Kits	5
Proteins	7
ChIP protocols with plant samples.....	9



Target	Modifications	Tested application	Cat. no.
Actin		WB	ab197345
		IHC-P, WB	ab230169
atpF		WB	ab65382
beta Tubulin		ICC/IF, IHC-P, IP, WB	ab15568
		ICC/IF, IHC-P, WB	ab231082
CAMKIV		IHC-Fr, WB	ab3557
CSN5		WB	ab195635
EIF2S1		Dot, Flow Cyt, IHC-Fr, IHC-P, WB	ab32157
		Dot, Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	ab214434
gamma Tubulin		Flow Cyt, ICC/IF	ab191114
Glucocorticoid Receptor alpha		ICC/IF, IHC-P, WB	ab3580*
Histone H3		ChIP, EM, IHC-P, IP, WB	ab1791*
	acetyl K9	ChIP, ICC/IF, IHC-P, IP, WB	ab10812*
	acetyl K9	ChIP, ICC/IF, IHC-P, WB	ab12179*
	acetyl K27	ChIP, ICC/IF, IHC-P, PepArr, WB	ab4729*
	di methyl K4	ChIP, CHIPseq, Flow Cyt, ICC/IF, IHC-Fr, IHC-P, IP, WB	ab173324
	di methyl K4	ChIP, Flow Cyt, ICC/IF, IHC-Fr, IHC-P, IP, WB	ab32356
	di methyl K4, tri methyl K4	ChIP, Flow Cyt, ICC, ICC/IF, WB	ab237973
	di methyl K4, tri methyl K4	ChIP, Flow Cyt, ICC/IF, WB	ab6000*
	di methyl K9	ChIP, ELISA, IHC-P, IP, WB	ab1220*
	di methyl K27	ChIP, ICC/IF, IHC-P, IP, WB	ab24684*
	tri methyl K4	ChIP, ELISA, Flow Cyt, ICC, ICC/IF, IHC-Fr, PepArr, WB	ab1012*
	tri methyl K4	ChIP, ELISA, Flow Cyt, ICC, ICC/IF, PepArr, WB	ab185637
	tri methyl K4	ChIP, ICC/IF, IHC-Fr, IP, PepArr, WB	ab8580*
	tri methyl K27	ChIP, ELISA, ICC/IF, IHC - Wmt, IHC-P, WB	ab6002*
	tri methyl K36	ChIP, ICC/IF, IHC-P, WB	ab9050*
	Histone H4	acetyl K5	ChIP, ELISA, ICC/IF, IHC-P, IP, WB
Hsp70		Flow Cyt, ICC/IF, IHC-P, IP, WB	ab5439
KDEL Receptor		IHC-P, WB	ab96720
thylakoid membrane proteins		WB	ab65574
xyloglucan		ICC/IF	ab190146

*ChIP-grade antibody

Highlighted

Target	Tested application	Cat. no.
Cytochrome f	WB	ab243204
Sucrose-Phosphate Synthase	WB	ab243200
Sulfite Reductase	WB	ab243201

Anti-Histone H3 antibody - Nuclear Loading Control and ChIP Grade (ab1791)

★★★★★ Reviews (212) Specific References (3065)

Description: Rabbit polyclonal to Histone H3 - Nuclear Loading Control and ChIP Grade
 Reactivity: Mouse, Rat, Chicken, Dog, Human, and more.

Call +44 (0)1223 696000 or contact us

abcam





Biochemicals

Product name	Purity	CAS number	Cat. no.
1-Triacontanol, Naturally occurring plant-growth regulator	= 98%	593-50-0	ab145129
3-Indole-acetic acid sodium salt, plant growth regulator	> 98%	6505-45-9	ab146403
3-Indolebutyric acid (Indolebutyric acid), Plant growth auxin	> 99%	133-32-4	ab143118
Astilbin, flavonoid found in plants	> 98%	29838-67-3	ab143575
Azure B, plant tissue stain		531-55-5	ab146367
Difenoconazole, Plant fungicide	> 98%	119446-68-3	ab142287
Indole-3-butyric acid potassium salt, auxin-family plant hormone	> 98%	60096-23-3	ab146396
Puerarin, Plant derived isoflavone	> 98%	3681-99-0	ab142939
trans-Zeatin, Plant growth hormone	> 98%	1637-39-4	ab145640

3-Indole-acetic acid sodium salt, plant growth regulator
(ab146403)

[References \(2\)](#)

Purity: > 98%

Call +44 (0)1223 696000 or contact us

abcam

ab146403





Assay Kit	Assay type	Sample type	Cat. no.
Cellulase	Enzyme activity (quantitative)	Tissue Lysate	ab189817
D-Mannitol	Quantitative - Colorimetric	Cells/tissues	ab155890
DNA demethylase	Quantitative	Cell culture extracts, Tissue Extracts, Nuclear Extracts	ab156908
DNMT	Enzyme activity (quantitative)		ab113467
	Enzyme activity		ab113468
Glucose-1-Phosphate	Quantitative	Tissue, Adherent cells, Suspension cells, Tissue Extracts	ab155892
Hexokinase	Enzyme activity (quantitative)	Cell Lysate, Tissue Lysate	ab211103
Histone Demethylase	Enzyme activity	Tissue Extracts, Cell Lysate	ab113455
Histone H3	Quantitative	Cell culture extracts, Tissue Extracts, Nuclear Extracts	ab156914
Histone H3 modifications	Quantitative	Tissue, Adherent cells, Suspension cells	ab185910
KDM1 / LSD1	Enzyme activity		ab113459
	Enzyme activity		ab113460
KDM6A/KDM6B	Enzyme activity	Cell culture extracts, Tissue Extracts, Cell Lysate, Nuclear Extracts, Purified protein	ab156910
	Enzyme activity	Nuclear Extracts, Purified protein	ab156911
NAD/NADH	Quantitative	Cell Lysate, Tissue Lysate	ab176723
Oxalate	Quantitative	Urine, Serum, Plasma, Tissue	ab196990
Oxalate decarboxylase	Enzyme activity	Tissue	ab196993
SIRT1	Sandwich	Nuclear Extracts, Purified protein	ab156915
TET1	Enzyme activity	Cell culture supernatant, Tissue, Nuclear Extracts, Purified protein	ab156912
	Enzyme activity	Nuclear Extracts, Purified protein	ab156913

Highlighted

Extraction kit	Component	Size	Cat. no.
ChIP Kit	100X Protease Inhibitor Cocktail, 5X Lysis Buffer I, 8-Well Assay Strips (with Frame) , 8-Well Strip Caps, Antibody Buffer, Anti-H3K9me2 (1 mg/mL), Binding Buffer, ChIP Dilution Buffer, DNA Release Buffer, Elution Buffer, F-Collection Tube	24 tests/ 48 tests	ab117137
Chromatin Extraction	1000X Protease Inhibitor Cocktail, 5X Lysis Buffer, Chromatin Buffer, Extraction Buffer A, Extraction Buffer B, Extraction Buffer C	50 tests	ab156906
Tissue extraction	Homogenization Pestles, Homogenization Tubes, Plant Extraction Buffer, Protease Inhibitor Cocktail (200X)	50 tests	ab206999

DNMT Activity Assay Kit (Colorimetric) (ab113467)

★★★★★ Reviews (2) Specific References (8)

Reactivity: Mouse, Rat, Human, Drosophila melanogaster, and more.

Call +44 (0)1223 696000 or contact us

abcam





Target	Species	Tag	Protein length	Application	Cat. no.
Autolysin protein	Chlamydomonas reinhardtii	His tag	Full length protein	SDS-PAGE	ab226279
BAS1	Arabidopsis thaliana	His tag	Full length protein	SDS-PAGE	ab241226
FLS1	Arabidopsis thaliana	His tag N-Terminus	Full length protein	SDS-PAGE	ab239438
Gliadin protein	Wheat	Histidine tag C-terminal	Protein fragment	ELISA, SDS-PAGE, Western blot	ab124981
HDT2	Arabidopsis thaliana	His tag N-Terminus	Full length protein	SDS-PAGE	ab241230
NADPH-dependent oxidoreductase 2-alkenal reductase protein	Arabidopsis thaliana		Full length protein	SDS-PAGE	ab225591
OSM34	Arabidopsis thaliana	His tag N-Terminus	Full length protein	SDS-PAGE	ab225564
Osmotin	Tobacco	His tag N-Terminus	Full length protein	SDS-PAGE	ab224842
Polyadenylate-binding protein RBP47A	Arabidopsis thaliana	His tag N-Terminus	Full length protein	SDS-PAGE	ab238280
Trypsin/alpha-amylase inhibitor CMx1/CMx3 protein	Wheat	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	Full length protein	SDS-PAGE	ab238371
Trypsin/alpha-amylase inhibitor CMX2 protein	Wheat	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	Full length protein	SDS-PAGE	ab238385

Recombinant A. thaliana HDT2 protein (His tag) (ab241230)

Expression system: Yeast
Application: SDS-PAGE
Species: Arabidopsis thaliana
Protein length: Full length protein

Call +44 (0)1223 696000 or contact us

abcam





OP

017850

ChIP using plant samples: Arabidopsis

Detailed procedure and tips for ChIP using plant samples. Edited from a protocol kindly provided by Werner Aufsatz.

Introduction

Eukaryotic chromatin is a complex of DNA and associated histone proteins that are involved in the higher order packaging of DNA into chromosomes. The chromatin state of a given DNA sequence influences transcriptional activity and replication timing and is regulated by potentially reversible covalent modifications of DNA and histones.

Histone modifications at conserved lysine and arginine residues within the flexible N-terminal tails, such as phosphorylation, acetylation and methylation, specify a code that serves as an interaction platform with specific domains of chromatin-associated proteins.

The immunoprecipitation (IP) of cross-linked chromatin with antibodies specific for certain histone modifications (chromatin immunoprecipitation, ChIP), followed by PCR to detect a potential enrichment or depletion of a DNA sequence of interest within IP fractions, constitutes an elegant and direct method to query specific chromatin states of individual genes and is already routinely used in many labs.

In contrast to animal cells, however, plant cells have a rigid cell wall which poses limitations to the simple utilization of protocols established for animals. In this protocol, the method described is used to study histone modifications in the plant model organism *Arabidopsis thaliana*. This protocol is an adapted version of the original procedure published by Lawrence and co-workers (Lawrence et al., 2004).

This protocol describes how chromatin is prepared from *Arabidopsis*, which can subsequently be used for chromatin immunoprecipitation (ChIP). The exact chromatin concentration should be determined before starting the X-ChIP assay. View our cross-linking chromatin immunoprecipitation (X-ChIP) protocol, which should be used after the chromatin preparation detailed below.

Procedure

Arabidopsis seeds are stratified for 48 hr in 0.1% Phytoblend w/v at 4°C and then sown onto soil. 1.5 g of whole, 3-4week-old seedlings, are used per chromatin preparation. It is imperative to avoid contamination with soil as much as possible during harvest.

Day 1

Chromatin cross-linking

1. Harvest 1.5 g seedlings and place them into a 50 ml tube.
2. Rinse seedlings twice with 40 ml double distilled water. Remove as much water as possible after second wash.
3. Add 37 ml 1% w/v formaldehyde solution. Gently submerge seedlings at the bottom of the tube by stuffing the tube with nylon mesh. Screw on cap and poke cap with needle holes. Put in exsiccator and draw vacuum for 10 min.
4. Release vacuum slowly and shake exsiccator slightly to remove air bubbles. Seedlings should appear translucent.
5. Add 2.5 ml 2 M glycine to quench cross-linking. Draw vacuum for 5 min.
6. Again, release vacuum slowly and shake exsiccator slightly to remove air bubbles.
7. Remove nylon mesh, decant supernatant and wash seedlings twice with 40 ml of double distilled water. After second wash, remove as much water as possible and put seedlings between two layers of kitchen paper. Roll up paper layers carefully to remove as much liquid as possible.

At this step, plant material can be snap-frozen in liquid nitrogen and stored at -80°C.

Chromatin preparation

1. Pre-cool mortar with liquid nitrogen. Add two small spoons of white quartz sand and plant material. Grind plant material to a fine powder.
2. Use cooled spoon to add powder to 30 ml of extraction buffer 1 stored on ice. Vortex to mix and keep at 4°C until solution is homogenous.
3. Rotate for 30 min at 4°C on a turning wheel or equivalent.
4. Filter extract through into a new, ice-cold 50 ml conical tube. Press to recover extract from solid material.
5. Repeat step 4.
6. Centrifuge extract at 4000 rpm for 20 min at 4°C.
7. Gently pour off supernatant and resuspend pellet in 1 ml of extraction buffer 2 by pipetting up and down. Transfer solution to Eppendorf tube.
8. Spin in cooled benchtop centrifuge at 13000 rpm for 10 min.
9. Remove supernatant and resuspend pellet in 300 µl of extraction buffer 2 by pipetting up and down.
10. Add 300 µl of extraction buffer 3 to fresh Eppendorf tube. Use pipette to carefully later solution from step 9 onto it.
11. Spin in cooled benchtop centrifuge at 13000 rpm for 10 min. In meantime, prepare 10 ml nuclei lysis buffer and 20 ml ChIP dilution buffer. Put buffers in cold room.
12. Remove supernatant and resuspend pellet in 300-500 µl of cold nuclei lysis buffer. Resuspend by pipetting up and down and by vortexing. Keep solution cold between vortexing. Incubate for 20 min on ice.

Chromatin preparation (continued)

13. Remove 10 μ l to run on an agarose gel.
14. Sonicate for 10 min at 4°C with sonicator: Setting "HIGH", 10 sec "ON cycle", 45 sec "OFF cycle". Make sure that the solution does not foam during sonication, e.g., by cooling the tubes with a mix of 100% w/v ethanol in ice during the sonication step.
15. Spin in cooled benchtop centrifuge at 13000 rpm for 10 min. Add supernatant to new Eppendorf tube.
16. Repeat step 14. Remove 10 μ l to run on an agarose gel.
17. Separate aliquots from steps 12 and 15 on 1.5% w/v agarose gel. In the sonicated samples, DNA should be shifted and more intense compared to untreated samples and range between 200-2000 bp, centering around 500 bp. Following step 16, the chromatin samples can be "snap-frozen" in liquid nitrogen and stored at -80°C. Repeated freezing/thawing cycles, however, should be avoided.

Pre-clearing and immunoprecipitation (IP)

1. Transfer chromatin to 15 ml Falcon tubes. Dilute 1/10 with fresh, ice-cold ChIP dilution buffer.
2. Prepare protein A agarose beads pre-adsorbed with sheared salmon sperm DNA by rinsing the required amount of beads 3 times with 1 ml ChIP dilution buffer in an Eppendorf tube. Spin in cooled benchtop centrifuge for 30 sec at 13000 rpm between the washes to pellet the beads. After the last wash, resuspend the beads with ChIP dilution buffer to yield a 25% slurry (for pipetting reasons).
3. Pre-clear each chromatin sample by adding 140 μ l of washed beads. Rotate for 1 hr at 4°C.
4. Spin Falcon tubes in a cooled centrifuge for 3 min at 3000 rpm to pellet the beads and transfer the supernatant to a new Falcon tube. Be careful not to carry over beads.
5. Store at 60 μ l aliquot of pooled chromatin at -20°C. This will serve as input control later on.
6. Add 600 μ l of chromatin solution per IP to an Eppendorf tube with an appropriate antibody. The optimal chromatin/antibody ratio has to be determined empirically for each antibody used. As a "rule of thumb", use about 10 μ g of antibody per IP for polyclonal antibodies. For monoclonal antibodies, usually a 1.5x to 4x higher concentration has to be used. Also set up 600 μ l of chromatin solution with an unrelated antibody as a mock IP.
7. Add 80 μ l of washed Protein A agarose beads to the chromatin IPs (25% slurry; for the preparation see step 2). Rotate overnight at 4°C.

Day 2

Collection, washes and elution of immune complexes

1. Prepare fresh elution buffer and place it at 65°C.
2. Spin IPs in cooled benchtop centrifuge at 5000 rpm for 30 sec to collect beads and discard the supernatant.
3. Add 1 ml of low salt wash buffer per tube. Rotate for 5 min at 4°C.
4. Spin in cooled benchtop centrifuge at 5000 rpm for 30 sec to collect beads and discard supernatant.
5. Add 1 ml of high salt wash buffer per tube. Rotate for 5 min at 4°C.
6. Spin in cooled benchtop centrifuge at 5000 rpm for 30 sec to collect beads and discard supernatant.
7. Add 1 ml of LiCl wash buffer per tube. Rotate for 5 min at 4°C.
8. Spin in cooled benchtop centrifuge at 5000 rpm for 30 sec to collect beads and discard supernatant.
9. Add 1 ml of TE buffer per tube. Rotate for 5 min at 4°C.
10. Spin in cooled benchtop centrifuge at 5000 rpm for 30 sec to collect beads and discard supernatant.
11. Repeat TE wash. Spin in cooled benchtop centrifuge at 5000 rpm for 30 sec to collect beads and discard supernatant.
12. Elute immune complexes by adding 250 µl of elution buffer. Vortex briefly to mix and incubate at 65°C for 15 min. Spin in benchtop centrifuge at 13000 rpm for 30 sec and transfer supernatant to a fresh Eppendorf tube.
13. Repeat elution and finally combine the two elutes.

Reverse cross-linking

1. Add 20 µl of 5M NaCl to samples. Incubate overnight at 65°C.
2. Add 109 µl of TE buffer, 7/1 µl of 5 M NaCl and 8.7 µl 20% SDS to the 60 µl input control aliquoted on Day 1. Incubate overnight at 65°C.

Day 3

DNA cleanup

1. Add 10 µl of 0.5M EDTA, 20 µl 1M Tris-HCl pH 6.5 and 1 µl of 20 mg/ml proteinase K to the IP samples. Add 1.2 µl of 0.5M EDTA, 2.4 µl 1M Tris-HCl pH 6.5 and 1 µl of 20 mg/ml proteinase K to the input control samples. Incubate for 1-3 hr at 45°C. The samples can be gently shaken through this incubation.
2. Purify the DNAs using a silica-gel membrane (e.g. PCR purification kit). Elute DNAs twice with 50 µl 10 mM Tris-HCl pH 8.0 and pool elutes. Proceed to PCR reactions. An excellent guide to quantification of ChIP by real-time PCR is provided by Haring et al. (2007).

Materials and reagents

Extraction buffer 1

- 0.4 M Sucrose
- 10 mM Tris-HCl, pH 8.0
- 10 mM MgCl₂
- 5 mM β-mercaptoethanol
- Protease inhibitors

Extraction buffer 2

- 0.25 M Sucrose
- 10 mM Tris-HCl, pH 8.0
- 10 mM MgCl₂
- 1% w/v Triton X-100
- 5 mM β-mercaptoethanol
- Protease inhibitors

Extraction buffer 3

- 1.7 M Sucrose
- 10 mM Tris-HCl, pH 8.0
- 2 mM MgCl₂
- 0.15% w/v Triton X-100
- 5 mM β-mercaptoethanol
- Protease inhibitors

Nuclei lysis buffer

- 30 mM Tris-HCl, pH 8.0
- 10 mM EDTA
- 1% w/v SDS
- Protease inhibitors

Protease inhibitor

- 1 mM PMSF (final concentration)
- Protease inhibitor cocktail (follow manufacturer's instructions)

ChIP dilution buffer

- 1.1% Triton X-100
- 1.2 mM EDTA
- 16.7 mM Tris-HCl, pH 8.0
- 167 mM NaCl

Elution buffer

- 1% SDS
- 0.1 M NaHCO₃

Low salt wash buffer

- 150 mM NaCl
- 0.1% SDS
- 1% Triton X-100
- 2 mM EDTA
- 20 mM Tris-HCl, pH 8.0

High salt wash buffer

- 500 mM NaCl
- 0.1% SDS
- 1% Triton X-100
- 2 mM EDTA
- 20 mM Tris-HCl, pH 8.0

LiCl wash buffer

- 0.25 M LiCl
- 1% Nonident P-40
- 1% sodium deoxycholate
- 1 mM EDTA
- 10 mM Tris-HCl, pH 8.0

TE buffer

- 10 mM Tris-HCl, pH 8.0
- 1 mM EDTA

References:

- Haring M, Offerman S, Danker T, Horst I, Peterhansel C, Stam M (2007). Chromatin immunoprecipitation: optimization, quantitative analysis and data normalization. *Plant methods*, 3:11.
- Lawrence RJ, Earley K, Pontes O, Silva M, Chen ZJ, Neves N, Viegas W and Pikaard CS (2004). A concerted DNA methylation/histone methylation switch regulates rRNA gene dosage control and nucleolar dominance. *Mol Cell*, 13: 599-609.
- Lippman Z, Gendrel AV, Black M, Vaughn MW, Dedhia N, McCombie WR, Lavine K, Mittal V, May B, Kasschau KD, Carrington JC, Doerge RW, Colot V and Martienssen R (2004). Role of transposable elements in heterochromatin and epigenetic control. *Nature*, 430: 471-476.
- Xie Z, Johansen LK, Gustafson AM, Kasschau KD, Lellis AD, Zilberman D, Jacobsen SE and Carrington JC (2004). Genetic and functional diversification of small RNA pathways in plants. *PLoS Biol.*, 2: 642-652.

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

Новосибирск

пр. Академика
Лаврентьева, д. 6/1
тел.
(383) 328-0048
nsk@dia-m.ru

Казань

ул. Парижской
Коммуны, д. 6
тел.
(843) 210-2080
kazan@dia-m.ru

С.-Петербург

ул. Профессора
Попова, д. 23
тел.
(812) 372-6040
spb@dia-m.ru

**Ростов-
на-Дону**

пер. Семашко, д. 114
тел.
(863) 303-5500
rnd@dia-m.ru

Пермь

Представитель
тел.
(342) 202-2239
perm@dia-m.ru

Воронеж

Представитель
тел.
(473) 232-4412
voronezh@dia-m.ru

Армения

Представитель
тел.
(094) 01-0173
armenia@dia-m.ru

Узбекистан

Представитель
тел.
(90) 354-8569
uz@dia-m.ru