

BIO-PRINT

GEL DOCUMENTATION & FLUORESCENCE IMAGING

ДИА•М
современная лаборатория

www.dia-m.ru
заказ on-line



GEL DOCUMENTATION IMAGING

VILBER
Smart Imaging



The Bio-Print protocol driven image acquisition is as quick as intuitive: adjust your exposure, save, print or quantify.



The Bio-Print delivers more significant quantitative data compared to other imagers.



THE LABORATORY STANDARD SO EASY

The Bio-Print is the laboratory standard for DNA and protein gel imaging. The system is based on an advanced set of performance ideal to achieve high fidelity and quantitative scientific measurement. Its excellent imaging properties guarantee the best performance for nucleic acid and protein samples.

The images generated allow for critical image analysis thanks to their high depth of data. Moreover, they can easily be printed out in an attractive way for the laboratory book. The Super-Bright UV illumination and filter technology enhances the image quality especially for DNA and RNA gels. The compact darkroom is made entirely of stainless steel and aluminium.

DNA, RNA & PROTEIN GELS UNRIVALED SENSITIVITY

Fluorescence is the main method used for gene expression and protein detection. The fluorophore absorbs the excitation light, reaching a higher energy state. By returning to its former state, it emits fluorescent light. The aim of the imaging system is to separate the emitted light from the excitation light in order to obtain an optimum sample image.

Our systems can accommodate up to 6 excitation channels in the UV and visible area. This is ideal for a large array of applications such as DNA or RNA gels, 1D protein fluorescent gels, stain free gels... A large number of dyes could be used such as Ethidium Bromide, Sybr-Safe, Sybr Green, Gel-Red, Gel-Green, Sybr-gold, GFP, Pro-Q Emerald 300, Sypro-Ruby, FITC, DAPI...

Up to 6 excitation and emission channels in the visible and UV area.

COLORIMETRIC STAINS THE APPLICATION MASTER

A colorimetric stain can be used for protein detection after the bands have been separated by electrophoresis. The colorimetric stains require visible light illumination for sample visualization and documentation. Our systems offer several reflective and transmitting illumination options for colorimetric stains.

Our white light screen provides unmatched white light illumination. In our darkroom, LED-based reflective epi-white light excitation provides even white light illumination, compared to other systems with standard bulbs. A large number of stains could be used such as Coomassie blue, Silver stain, Ponceau S Red, Copper stain, Zinc stain... Our systems are ideal for colorimetric stained protein gels, X-ray film imaging, autorads, SSCP gels, colony dish and flask imaging...

For fluorescence, photobleaching and phototoxicity are reduced.

Bio-Print



ONE CLICK TO IMAGE

Automatic control of the camera, lens and lighting for an unrivaled ease of use



APPS STUDIO

A complete library of gel imaging applications to ensure reproducibility.




SUPER RESOLUTION

Very high level of detail and more quantitative information.



SCIENTIFIC USB-3 CAMERA

Interchangeable light pad for UV, blue, green & red fluorescence.



SIMPLY PRECISE ONE CLICK TO THE IMAGE

The Bio-Print has been designed for maximum ease of use. From its simple installation to its intuitive user interface, this system is plug-and-play. The Bio-Print software is the easiest software to take an image. Place your gel on the table, select your application, click on Start and automatically the system auto-exposes your

gel image, your marker image and combines the two together.

The Bio-Print includes our unique Apps Studio approach to imaging. The Apps Studio is a library which contains 40 different protocols for your blot and gel. The protocol oriented Apps Studio ensures reproducibility and one click acquisition for the best ease of use.

The Bio-Print has rich features and guides you into the advanced functions in a very ergonomically designed user interface. The simple and self-explanatory menu is ideal in a multi-user environment.

The Bio-Print includes our unique Dynamic 3D scan technology. The 3D Dynamic Scan images your sample in real time and reconstructs the data to create live three dimensional models. The 3D reconstruction provides direct information regarding the image dynamic, background level and protein or DNA quantity. A little change of exposure time will refresh the 3D view automatically.

BIO-PRINT

Upgradeable

- Add a white light conversion plate for colorimetric stained protein gels, X-ray film imaging, autorads, SSCP gels, colony dish and flask imaging
- Add a blue light conversion plate for blue light DNA/ RNA/ Protein fluorescence

Performance At A Budget

- Scientific camera made in Germany
- Motorized zoom lens
- USB-3 interface
- Camera passive cooling to maximize the signal to noise ratio
- Auto-exposure and auto-focus
- Automatic light illumination control




Ideal For Quantification

- Reproducible and comparable quantification data
- ImageMaster™ technology to obtain the optimum image for quantification
- Scientific TIFF file or proprietary file format
- Clarity™ technology for razor sharp band appearance without affecting data integrity

Easy & Intuitive

- One click to get the image
- Auto-exposure and auto-focus
- Automatic light illumination according to your application
- Protocol driven image acquisition



CAMERA & OPTICS

- CX4 camera:
- Scientific grade camera - Made in Germany
- 2 megapixels resolution extendable to 7,6 megapixels
- Passive cooling
- 16-bit – 65 536 grey levels
- USB-3 connection
- Motorized zoom lens with feedback
- Field of view: 20x20cm

EASE OF USE

One-Click-to-the-Image™

Auto-exposure

Auto-focus

Auto-lighting

HARDWARE

Smart Darkroom technology:

Software control of the lighting

White light led panels with automatic intensity adjustment

UV cut-off filter

Fix position transilluminator with optional protection screen

Steel and stainless steel darkroom for long lasting robustness. Wide access door with UV safety shut-off

SOFTWARE

Free Bio-Vision software for image acquisition with full GLP compliance. Molecular weight calculation, band quantification, colony counting, distance calculation, text annotation and image enhancement included.

CFR21 Part 11 ready

APPLICATIONS

DNA, RNA, protein gels. Fluorescence and colorimetric stain imaging:

- Ethidium bromide, Sybr-Safe, Sybr-Green, Gel-Red, Gel-Green, Sybr-Gold, GFP, Pro-Q Emerald, Sypro ruby, FITC, DAPI
- Colorimetric stained protein gels, X-Ray film, autorads, SSCP gels, colony dish and flask imaging and other EPI white light applications
- Coomassie blue, Silver stain, Ponceau S Red, Copper stain...

Sky-Pad or Blue light conversion screen for DNA/RNA detection (avoid «nicking» DNA):

- Sybr Safe, Sybr Green, eGFP

TECHNOLOGY & INNOVATION

- Apps Studio™
- 3D Dynamics Scan™
- SuperResolution™
- ImageMaster™ assistant

OPTIONS

White light or blue light conversion screen

Advanced Bio-1D quantification software

CFR21 Part 11 module

HEADQUARTER

Vilber Lourmat
ZAC de Lamirault
Collegien
F-77601 Marne-la-Vallee cedex 3
France
Phone: + 33 (0) 1 60 06 07 71
info@vilber.com

GERMANY

Vilber Lourmat
Deutschland GmbH
Wielandstrasse 2
D-88436 Eberhardzell
Deutschland
Phone : + 49 (0) 7355 931 380
info@vilber.de

CHINA

Vilber China
Room 127 Building A
N° 111 Yuquangying
Fengtai District – Beijing
China
Phone : + 86 1361 1131 545
info@vilber.cn

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krrsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemmerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

