




# RWD



In Vivo Imaging Solution



Neural Signal Regulation Solution



Pathological Section Solution



Cellular and Molecular Biology Research Solution




## PRODUCT CATALOGUE

[www.rwdstco.com](http://www.rwdstco.com)  
[rwd@rwdstco.com](mailto:rwd@rwdstco.com)

# Speed up Life Science Research Achieve Better Healthcare for Human & Animals



Learn more about us at:  
[www.rwdstco.com](http://www.rwdstco.com)





# CONTENTS

|  |           |   |           |
|--|-----------|---|-----------|
| <b>Animal Surgery and Modeling Solution</b> .....              | <b>01</b> | <b>Cellular and Molecular Biology Research Solution</b> ..... | <b>57</b> |
| Stereotaxic Instruments .....                                  | 1         | Single Cell Suspension Dissociator .....                      | 57        |
| Precise Impactor .....   | 12        | Cell Separation Kit .....                                     | 58        |
| Animal Anesthesia Machine .....                                | 14        | CO <sub>2</sub> Incubator .....                               | 59        |
| Animal Ventilator .....  | 23        | Automated Cell Counters .....                                 | 61        |
| Surgical Auxiliary Equipment .....                             | 25        | Microcentrifuges .....  | 63        |
| Veterinary Laser Therapy .....                                 | 30        | <b>Surgical Instruments &amp; Matrix</b> .....                | <b>66</b> |
| Smart Animal ICU .....   | 31        |   |           |
| Surgical Instruments & Matrix .....                            | 32        |   |           |
| <b>Neural Signal Regulation Solution</b> .....                 | <b>33</b> |   |           |
| Optogenetics System .....                                      | 33        |   |           |
| Fiber Photometry System .....                                  | 36        |   |           |
| Micropipette Puller .....                                      | 39        |   |           |
| Motorized Micromanipulator .....                               | 40        |   |           |
| <b>In Vivo Imaging Solution</b> .....                          | <b>41</b> |   |           |
| Laser Speckle Imaging System .....                             | 41        |   |           |
| <b>Animal Administration and Blood Sampling Solution</b> ..... | <b>43</b> |   |           |
| Drug Delivery Cannula .....                                    | 43        |   |           |
| Osmotic Pumps .....  | 48        |   |           |
| Blood Sampling and Administration System .....                 | 49        |   |           |
| Syringe Pumps .....  | 51        |   |           |
| <b>Pathological Section Solution</b> .....                     | <b>53</b> |   |           |
| Rotary Microtome .....   | 53        |   |           |
| Cryostat .....   | 55        |   |           |

## Introduction

Established in 2002, RWD Life Science has been dedicated to being a leading manufacturer in research and veterinary equipments worldwide. We specialize in anesthesia, physiological and surgical equipments, which are used by pre-clinical researchers and veterinarians. Our mission is to provide high quality products and integrated solutions to our customers with prompt and professional support.

For more information about our products, such as Laser Speckle Perfusion Imager, Anesthesia Machine, Ventilator, Stereotaxic Instrument, Syringe Pump, Cannula System and Surgical tools etc., please check our website:  
[www.rwdstco.com](http://www.rwdstco.com)

### RWD Life Science Co., Ltd

**Add:** 19-20/F, Building 7A, Shenzhen International Innovation Valley, Dashi 1st Road, Nanshan District, Shenzhen, Guangdong, China

**Web:** [www.rwdstco.com](http://www.rwdstco.com)

**E-mail:** [rwd@rwdstco.com](mailto:rwd@rwdstco.com)

### RWD Life Science Inc.

**Add:** 6540 Lusk Blvd, Suite C161 San Diego, CA 92121

**Tel:** 1-858-900-5879

**Tech Support:** [service@rwdstco.com](mailto:service@rwdstco.com)

## Animal Surgery and Modeling Solution


Stereotaxic instrument is an important research equipment widely used in many fields such as neuroanatomy, neurophysiology, neuropharmacology and neurosurgery. Up to now, RWD's stereotaxic instrument has helped global neuroscience researchers to publish more than 2300 articles in Cell, Nature, Science and other academic journals, and has been widely recognized in the industry.

- Meet the needs of rats, mice, guinea pigs, tree shrews, cats, rabbits, dogs, monkeys and other different animal experiments.
- Three kinds of accuracy 100 $\mu$ m, 10 $\mu$ m, 1 $\mu$ m are available.
- Except for the automated model 71000, other models can be upgraded to digital or dual manipulator models.



Automated Stereotaxic Instrument

Operation accuracy is 1 $\mu$ m;  
Built-in rodent brain atlas and three automatic programs.



# Small Animal Stereotaxic Instruments

## ▶ 71000 Automated Stereotaxic Instrument

The 71000 is an automated, high precise (1 $\mu$ m) and intelligent stereotaxic instrument integrated with rodent atlases and built-in software which allows for three automated procedures (skull window, tissue removal and multi-point programs) and It provides intuitive observation of the probe position relative to the brain area, which is much convenient, efficient and easy to operate. It greatly reduces error and damage caused by human operation. In addition, the anti-collision function would prevent animal from damage caused by wrong operation, which is much more safer.

### Features & Technical Parameters:

● **Precise**

1. The stereotaxic instrument uses a stepping motor to control the manipulator, with a resolution of 1 $\mu$ m, which meets the needs of higher precision experiments.
2. Anti-backlash function can greatly eliminate the error caused by the gap between the motor wheels when the motor moves in the reverse direction.
3. Calibration function, which can calibrate the actual scale value on the manipulator and the reading on the software.



● **Convenient and efficient**



1. The software integrates mouse and rat atlas, which can display the trajectory and position of the probe in real time.
2. Multiple control movement methods: micro-manipulation panel, software and computer keyboard to control the movement of the three-dimensional manipulators, the movement speed is adjustable.
3. When reaching the Bregma/Lambda point, you can set the Bregma/Lambda point as the origin with one key, and you can return to the origin with one key at any coordinate.
4. Skull window, tissue removal and multi-point programs, which can meet a variety of experimental needs, reduce injuries and errors caused by manual operations, and are easy to operate.

● **Safe**

1. Z-axis retraction and motional collision function to prevent animal damage during moving.
2. The emergency stop button can immediately stop the running program to avoid animal damage caused by operating errors.

### Tissue removal module (optional)

It is used for the removal of brain tissue.



Animal Surgery and Modeling Solution

**Features :**

- Configure a vacuum pump, the software can set the suction depth, sequence, and number of repetitions to automatically and accurately remove brain tissue;
- The software can set the depth and speed of lens implantation, one-click implantation;
- Ensure the flatness of the wound end surface and reduce damage;
- Improve the survival of neurons per unit area and enhance imaging effect;
- Reduce the error caused by traditional manual operation and enhance the repeatability of experiment.

**▶ Rotational Digital Stereotaxic Instrument**

The 69100 Rotational Digital Stereotaxic allows for quick and precise adjustment of the animal's skull level and alignment of the sagittal suture without the need of repeated setting Bregma or repeated fixing animal, which can greatly reduce the time required for skull level adjustment. The Alignment Indicator is used to indicate the level and also to improve the perpendicularity of the syringe or cannula to the brain, increasing of positioning

**Features:**

- The rotation adaptor can adjust head position in 3 dimensional space.
- The displacement distance of the head rotating adaptor X, Y and Z is 30mm, the left and right angles are rotated  $\pm 15^\circ$ , the front and rear swing angle is  $\pm 10^\circ$ , and the bottom plate is rotated  $\pm 10^\circ$ .
- 0.01mm resolution in digital display, supply by battery, no electronic interference.
- The UP indicator on z-axis prevents z-axis manipulate from misoperation.
- The manipulator drives by electrode more than 100,000 times with minimum deviation.
- Double head screw design makes the 3-axis more stable and accurate.



**Recommended Accessories:**



**40X Center Magnifier(69130)**

1. Used to assist in coinciding the Bregma point with the crossover point of the center height gauge, the crossover point is the crossover point of the rotation axis, the coordinates of point remain unchanged when rotating, so there is no need to repeat set up the Bregma point when adjusting the level.
2. Help to find Bregma and Lambda points.



**Alignment Indicator(69131)**


When the two probes touch the skull, the difference in the deflection angle of the micrometer pointer can feedback whether the height of the animal's head is horizontal from left to right, front to back. The resolution of the Alignment Indicator is 0.01mm.

**Order Information:**

| Model | Specification   |
|-------|---|
| 69100 | Ultra Precise Rotational Stereotaxic Instrument for Mouse, Digital, SGL M, 69110 Adaptor Incl.  |
| 69101 | Ultra Precise Rotational Stereotaxic Instrument for Mouse, Digital, Dual M, 69110 Adaptor Incl. |
| 69102 | Ultra Precise Rotational Stereotaxic Instrument for Rat, Digital, SGL M, 69112 Adaptor Incl.    |
| 69103 | Ultra Precise Rotational Stereotaxic Instrument for Rat, Digital, Dual M, 69112 Adaptor Incl.   |
| 69130 | Centering Scope 40X.  |
| 69131 | Alignment Indicator.  |

## ► Standard Stereotaxic Instrument

Classical U-type base design, professional services in the laboratory for 17 years!



Composition:


68001 and 68002 are complete sets of stereotaxic (see figure), both containing 68861N base plate, 68401 three-dimensional manipulator - left, 68021 rat adaptor, 68301 rat 18 ° ear bar, 68201 standard holder, and 68002 includes an extra 68402 three-dimensional manipulator - right.

### Order Information:

| Model  | Product Description                                  |
|--------|--|
| 68001  | Stereotaxic for Rat, SGL M, 18 Deg                   |
| 68002  | Stereotaxic for Rat, Dual M, 18 Deg                  |
| 68005  | Stereotaxic for Rat, SGL M, 45 Deg                   |
| 68006  | Stereotaxic for Rat, Dual M, 45 Deg                  |
| 68037  | Stereotaxic for Mouse, SGL M, 68030 Adaptor Incl     |
| 68038  | Stereotaxic for Mouse, Dual M, 68030 Adaptor Incl    |
| 68043  | Stereotaxic for Mouse, SGL M, 68055 Adaptor, 60 Deg  |
| 68044  | Stereotaxic for Mouse, Dual M, 68055 Adaptor, 60 Deg |
| 68535  | Stereotaxic for Mouse, SGL M, 68077 Adaptor, 60 Deg  |
| 68536  | Stereotaxic for Rat, SGL M, 68078 Adaptor, 18 Deg    |
| 68861N | Stereotaxic Base Plate, Standard                     |

## ▶ Desktop Digital Stereotaxic Instrument

Digital displacement display, higher accuracy, higher efficiency!




### Composition:

68025 and 68026 are complete sets of stereotaxic (see figure), both containing 68861N base plate, 68409 3-Axis Digital Manipulator, Left-hand, 68021 rat adaptor, 68301 rat 18 ° ear bar, 68201 standard holder, 68026 includes an extra 68410 3-Axis Digital Manipulator, Right-hand.

### Order Information:

| Model  | Product Description   |
|--------|---|
| 68025  | Stereotaxic for Rat, SGL M, Digital, 18 Deg                   |
| 68026  | Stereotaxic for Rat, Dual M, Digital, 18 Deg                  |
| 68027  | Stereotaxic for Rat, SGL M, Digital, 45 Deg                   |
| 68028  | Stereotaxic for Rat, Dual M, Digital, 45 Deg                  |
| 68018  | Stereotaxic for Mouse, SGL M, Digital, 68030 Adaptor Incl     |
| 68019  | Stereotaxic for Mouse, Dual M, Digital, 68030 Adaptor Incl    |
| 68045  | Stereotaxic for Mouse, SGL M, Digital, 68055 Adaptor, 60 Deg  |
| 68046  | Stereotaxic for Mouse, Dual M, Digital, 68055 Adaptor, 60 Deg |
| 68537  | Stereotaxic for Mouse, SGL M, Digital, 68077 Adaptor, 60 Deg  |
| 68538  | Stereotaxic for Rat, SGL M, Digital, 68078 Adaptor, 18 Deg    |
| 68861N | Stereotaxic Base Plate, Standard                              |


## ► Compact Mouse Stereotaxic Instrument



### Order Information:

| Model | Product Description                                 |
|-------|---|
| 68507 | Stereotaxic for Mouse, SGL M, Portable              |
| 68510 | Stereotaxic for Mouse, Dual M, Portable             |
| 68513 | Stereotaxic for Mouse, Digital and Portable, SGL M  |
| 68516 | Stereotaxic for Mouse, Digital and Portable, Dual M |

## ▶ Portable Rat & Mouse Stereotaxic Instrument



**Composition:**

68526 is a complete set stereotaxic configurations, containing 68867S base plate, 68401 three-dimensional manipulator - left, 68087 mouse adaptor and 68088 rat adaptor, 68301 rat 18° ear bars and 68306 mouse 60 ° ear bars, 68201 standard holder; Digital version 68528 adds digital display module.

**Order Information:**

| Model  | Product Description   |
|--------|---|
| 68526  | Stereotaxic for Rat and Mouse, SGL M, Portable              |
| 68527  | Stereotaxic for Rat and Mouse, Dual M, Portable             |
| 68528  | Stereotaxic for Rat and Mouse, SGL M, Digital and Portable  |
| 68529  | Stereotaxic for Rat and Mouse, Dual M, Digital and Portable |
| 68867S | Stereotaxic Base Plate for Rat and Mouse, Portable, SGL M   |
| 68867D | Stereotaxic Base Plate for Rat and Mouse, Portable, Dual M  |
| 68087  | Mouse adaptor (Stereotaxic for Rat and Mouse )              |
| 68088  | Rat adaptor (Stereotaxic for Rat and Mouse )                |

## ► Large Animal Stereotaxic Instrument

Adaptors are available for dogs, monkeys, pigs, and other large animals.

Unique jaw plate design better stabilizes the heads of large animals of varying weights.

Y-axis rail design, longer moving range:  $\pm 100\text{mm}$ , 1mm resolution

Combined with microdriver, the resolution of vertical reading can be promoted to  $50\mu\text{m}$

Holds up to six manipulator arms



### Composition:

68902 is a complete set of stereotaxic, containing 68868N base plate, 68941 two-dimensional manipulator - left, 68942 two-dimensional manipulator - right, 68081 dog/monkey adaptor, 68303 cat/ monkey 18 ° ear bar, 68201 standard holder; Digital version 68916 adds digital display module.

### Order Information:

| Model  | Product Description   |
|--------|---|
| 68901  | Stereotaxic for Dog/Monkey, SGL M, 68081 Adaptor Incl           |
| 68902  | Stereotaxic for Dog/Monkey, Dual M, 68081 Adaptor Incl          |
| 68911  | Stereotaxic for Cat/Monkey, SGL M, 68041 Adaptor Incl           |
| 68912  | Stereotaxic for Cat/Monkey, Dual M, 68041 Adaptor Incl          |
| 68916  | Stereotaxic for Dog/Monkey, Digital, SGL M, 68081 Adaptor Incl  |
| 68917  | Stereotaxic for Dog/Monkey, Digital, Dual M, 68081 Adaptor Incl |
| 68920  | Stereotaxic for Cat/Monkey, Digital, SGL M, 68041 Adaptor Incl  |
| 68921  | Stereotaxic for Cat/Monkey, Digital, Dual M, 68041 Adaptor Incl |
| 68868N | Stereotaxic Base Plate for Dog, Pig, Monkey etc                 |
| 68941  | 2-axis Manip for Dog/Monkey Stereotaxic, Left Hand              |
| 68942  | 2-axis Manip for Dog/Monkey Stereotaxic, Right Hand             |
| 68944  | 2-axis Digital Manip for Dog/Monkey Stereotaxic, Left Hand      |
| 68945  | 2-axis Digital Manip for Dog/Monkey Stereotaxic, Right Hand     |

## Stereotaxic Instruments Accessories

### Electrode Holder

68201 Standard Probe Holder-Corner



V-groove design, widely used for syringe needles or electrodes with diameter 0.3mm-1.5mm.

68211 Cannula Holder



Compatible with stereotaxic instruments with dual manipulators. Holds electrodes or cannulas of 2.1mm to 5.0mm in diameter.


68210 Electrode Holder



Holder is specifically designed for electrophysiological experiments. Two electrodes, controlled individually, can be located very close to each other. Holds thin electrodes with diameter from 0.7mm to 2.5mm.


### Cannula Holders

68205 Cannula Holder




It's made of metal. Suitable for holding the single cannula or double cannula with the outer diameter of 3.5mm plastic.

68217 Cannula Holder



Suitable for holding the single cannula or double cannula with the outer diameter of 3.5mm plastic base.

68214/68215 Ceramic Ferrule Holder



Hold cannula or Ceramic Ferrule during implantation.

- Versions for Ø1.25 mm — 68214 Cannula or Ceramic Ferrule.
- Versions for Ø2.5 mm — 68215 Cannula or Ceramic Ferrule.

### Multi-function Holders

68202 Large Probe Holder




Lateral fixation, holds micro syringe needles, large diameter electrodes and flexible drill handles etc, V-groove design with 6mm-12mm diameters.

68218 Syringe Holder



Holder is applicable for micro-injection, which holds syringe and needle simultaneously, and the clamping of the syringe is more stable, the syringe holding range is 3mm to 12mm and needle's is 0.3mm to 1.5mm.

68605 Microdrill Holder




Holder secures the microdrill to stereotaxic instruments. Controls the depth of drilling by the manipulator. Easy and accurate operation protects brain tissue from damage caused by excessive drilling. Holds the microdrill with diameter about 14.5mm.

## Adaptors


### Mouse Adaptors

68030 Mouse/Neonatal Rat Adaptor  
(with 18° ear bars 68034 and 45° ear bars 68313)



This adaptor can not only fix the mouse and neonatal rat head, but also avoid damage to the skull. Four kinds of earstick tips are available, the height of the earstick and earstick are adjustable, and the adjustment range is 20mm.

68055 Mouse Adaptor




It holds the mouse head by a palate and incisor, and offers vertical regulation to make the head level or at certain angle. The Vertical adjustable range from +10mm to -20mm, with 100µm resolution; Horizontal adjustable range of 43.5mm; Horizontal rotatable range of 35°.

### 68057 Mouse Gas Anesthesia Head Holder


This model offers maximal adjustment for positioning the head. The adaptor can rotate the head to multi angle without blocking the anesthetic gas delivery.

- Vertical adjustable range: 30mm (+10mm~-20mm), with 100µm accuracy.
- Vertical rotatable range: 70° (-35°~ +35°), with 5° accuracy.
- Horizontal rotatable range: 70° (-35°~ +35°), with 5° accuracy.
- Horizontal adjustable range: 43.5mm
- Diameter of inlet/outlet of the mask: 2mm



## Rat Adaptors

### 68021 Rat Adaptor



It holds rat head firmly through incisor, nose clip and ear bars. This adaptor can be adjusted 30 mm vertically along dovetail slide with 100µm resolution and 50mm horizontally, which makes it applicable for rats of different weight.

### 68053 Rat Anesthesia Adaptor



This adaptor integrates the gas anesthesia mask function while hold the rat head firmly. The adaptor can rotate the head to multi angle without blocking the anesthetic gas delivery.

- Vertical adjustment range is 30mm (-20mm ~ + 10mm), with 100µm accuracy.
- Vertically rotatable range is 70°(-35° ~ + 35°), with 5°accuracy.
- Horizontal rotatable range is 70° (-35°~ + 35°), with 5°accuracy.
- Horizontal adjustment range is 43.5mm.

## Multi-function Adaptor

### 68091 Spinal Cord Adaptor for Rat & Mouse

The adaptor is ideal for spinal securing of rat and mouse. The base plate is 210mm(L) x 165mm(W) x 15mm(H), which can be attached to RWD 68000 series small animal stereotaxic instruments. There are 3 vertical bars and locking blocks. which can be easily disassembled through the T slots without wrench. The adjustment range of vertical bar is 56mm. The height of locking block can be easily adjusted by screw. This adaptor can be adjusted in Vertical/Lateral, Anterior/Posterior and Ventral/Dorsal directions. Standard accessories contain: Base plate, 3 vertical bars and locking blocks and spinal cord clamp (3 pairs).

#### Specifications:

- Accurate and reliable positioning.
- Wide various components and accessories can be used in module.
- Simple, fast and effective clamping and locking.
- Compatible with small animal stereotaxic instruments.



## Ear Bars

### RWD ear bar with the scale of 35mm, 1mm resolution, square 6.15mm

68301




68302




68306




68314





## Stereotaxic Frame Nosecone Masks

- This type of anesthetic masks are specifically designed to fit most stereotaxic frame and deliver anesthetic gases with minimal exposure during surgery. They are locked onto the adaptor with the cone opening toward the open "U". The animal's teeth is placed over the incisor bar.
- Two kinds of Stereotaxic Frame Nosecone Masks are available: one is "passive" mask with collinear inlet and outlet, and the other one is "active" mask with separate (non-collinear) inlet and outlet on both sides. the "passive" mask is used with charcoal canister only, and the "active" mask can be used with evacuation apparatus.



### Order Information:

| Model | Product Description   | Remarks                |
|-------|---|------------------------|
| 68601 | Stereotaxic Anesthesia Mask Kit -mice or neonatal rats (<70g, Canister & Tubing included) | For Passive Scavenging |
| 68602 | Stereotaxic Anesthesia Mask Kit - rats (<300g, Canister & Tubing included)                |                        |
| 68663 | Stereotaxic Anesthesia/Evacuation mask- mice or neonatal rats (<30g)                      | For Active Scavenging  |
| 68665 | Stereotaxic Anesthesia/Evacuation mask- mice or neonatal rats (30~70g)                    |                        |
| 68666 | Stereotaxic Anesthesia/Evacuation mask- rats (<300g)                                      |                        |
| 68667 | Stereotaxic Anesthesia/Evacuation mask- rats (>300g)                                      |                        |

## Precise Impactor

### 68099 II Precise Impactor

The 68099 II Precise Impactor is an instrument used for traumatic brain injury and spinal cord injury model. The system adopts pneumatic-electric control, which can precisely adjust the speed, depth and dwell time to achieve precise impact. Touch screen design, convenient and friendly operation. The zero method uses a sensor contact mechanism to automatically detect the zero interface, which is sensitive and efficient.



#### Features:

- Pneumatic-electric design, stable blow without shaking, to ensure the repeatability of the experiment;
- It is suitable for small animal brain and spinal cord injury models such as mice, rats, and guinea pigs;
- LCD touch screen, wide-angle visibility, convenient and friendly operation;
- Sound alarms and text message prompts improve the human-computer interaction between the equipment and the user, and reduce human errors in the operation process;
- The zero-contact method can automatically determine the zero interface to reduce human error;
- A variety of cylindrical flat-head head hammers with outer diameters of 1, 2, 3, 4 and 5mm, which can meet the needs of traumatic brain injury and spinal cord injury model;
- Speed calibration function to ensure the stability and repeatability of the experiment;
- Speed range is 0.5-5.6m/s, resolution 0.1m/s; depth range is 0.00-5.00mm, resolution 0.01mm; dwell time range is 0.00-5.00s, resolution 0.01s;
- The rotatable design ensures that the brain and spinal cord are hit vertically and improves the repeatability of the experiment.

#### Traumatic Brain Injury(TBI)

- Two kinds of brain adaptors for rats and mice are optional. The rotational adaptor can rotate left and right  $\pm 30^\circ$ , the height adjustment range of  $+10\sim -20\text{mm}$ , and the precise rotational adaptor can rotate left and right  $\pm 30^\circ$ , swing up and down  $\pm 10^\circ$ , three-axis displacement 30mm;
- The base is movable, 0~116mm in X direction and 0~98mm in Y direction.



Precise Rotational Adaptor




Rotational Adaptor



Head hammers

### Spinal Cord Injury(SCI)

- The rat and mouse spinal cord adaptor adopts the method of supporting the spinal cord from the lower sides of the spine to avoid the depth error caused by the collapse of the animal when breathing or impacting;
- The three-axis micro-manipulated rotatable platform includes a universal adjustment device, which can adjust the level of the spinal cord to ensure vertical strikes and improve experiment repeatability.



Spinal Cord Adaptor



Three-axis Micro-manipulated Rotatable Platform

#### Order Information:

| Model         | Product Description                          | Remarks   |
|---------------|--|---|
| 68099 II      | Precise Impactor-Brain                       | Standard configuration (without adaptor)  |
| 68099 II -S-M | Precise Spinal Impactor for Mice             | Standard configuration (including adaptor)  |
| 68099 II -S-R | Precise Spinal Impactor for Rat              | Standard configuration (including adaptor)  |
| 68063         | Rotational Rat Adaptor, 60°Ear Bars Incl.    | Optional  |
| 68065         | Rotational Mouse Adaptor, 60°Ear Bars Incl.  | Optional  |
| 69100-03      | Precise Rotational adaptor for Rat           | Optional, can fix the head of mouse collaborately use mouse anesthesia adaptor(69110 or 69111) and ear bars(69114). |
| 69100-04      | Precise Rotational adaptor for Mouse         | Optional, can fix the head of rat collaborately use rat anesthesia adaptor(69112) and ear bars(69115).              |
| 69110         | Mouse Anesthesia Adaptor with Tubing(<30g)   | Optional (69100-04 including)   |
| 69111         | Mouse Anesthesia Adaptor with Tubing(30~70g) | Optional  |
| 69112         | Rat Anesthesia Adaptor with Tubing           | Optional (69100-03 including)   |
| 69114         | 60°Ear Bars for Mouse, One pair              | Optional (69100-04 including)   |
| 69115         | 18°Ear Bars for Rat, One pair                | Optional (69100-03 including)   |


# Animal Anesthesia Machine

## Small Animal Anesthesia Machines ▶

### R500 Compact Anesthesia Machine

Regular assembly (R500-A), vaporizer(R580), Inlet/Outlet End Caps(R580-IO)and tube accessory kit are included in the R500 anesthesia machine. Different air supply, induction chambers, small animal masks, filter canisters and gas evacuation apparatuses can be selected for anesthesia recovery of different animals.

Suitable for rats, mice, hamsters, guinea pigs, rabbits, cats and other animals within 7kg



Note: Oxygen flush is commonly used to flush out anesthetic gas to the filter canister before opening the induction chamber to prevent the operators from anesthetic gas.

#### Order Information:


| Model  | Product Description   | Configuration | Remarks   |
|--------|---|---------------|---|
| R500IE | Compact Small Animal Anesthesia Device-Isoflurane, Easy Fill, 0-4L  | ●             | Contains connection tubing for mask&induction chamber |
| R500IP | Compact Small Animal Anesthesia Device-Isoflurane, Pour Fill, 0-4L  | ●             |   |
| R500IK | Compact Small Animal Anesthesia Device-Isoflurane, Key Fill, 0-4L   | ●             |   |
| R500SE | Compact Small Animal Anesthesia Device-Sevoflurane, Easy Fill, 0-4L | ●             |   |
| R500SP | Compact Small Animal Anesthesia Device-Sevoflurane, Pour Fill, 0-4L | ●             |   |
| R500SK | Compact Small Animal Anesthesia Device-Sevoflurane, Key Fill, 0-4L  | ●             |   |
| R510-S | Portable Stand for R510 Series Anesthesia Machine                   | ○             | Upgrade to portable type                              |

Note: ●Recommended configuration ○ Optional

## R550 Multi-output Anesthesia Machine

Regular assembly (R550-A), vaporizer (R580), Inlet/Outlet End Caps (R580-IO) and tube accessory kit are included in the R550 anesthesia machine. Different air supply, induction chambers, small animal masks, filter canisters and gas evacuation apparatuses can be selected to form the experimental platform for anesthesia recovery of different animals.

With the latest design, it can be applies to simultaneous anesthesia of multiple mice, rats, cats, rabbits and other animals (1-5 small animals), and every channel can be controlled independently.



### Order Information:


| Model  | Product Description  | Configuration | Remarks  |
|--------|--|---------------|--|
| R550IE | Enhanced Small Animal Anesthesia Device-Isoflurane, Easy Fill, 0-4L  | ●             | Contains connection tubing for mask & induction chamber. |
| R550IP | Enhanced Small Animal Anesthesia Device-Isoflurane, Pour Fill, 0-4L  | ●             |  |
| R550IK | Enhanced Small Animal Anesthesia Device-Isoflurane, Key Fill, 0-4L   | ●             |  |
| R550SE | Enhanced Small Animal Anesthesia Device-Sevoflurane, Easy Fill, 0-4L | ●             |  |
| R550SP | Enhanced Small Animal Anesthesia Device-Sevoflurane, Pour Fill, 0-4L | ●             |  |
| R550SK | Compact Small Animal Anesthesia Device-Sevoflurane, Key Fill, 0-4L   | ●             |  |
| R510-S | Portable Stand for R510 Series Anesthesia Machine                    | ○             | Upgrade to portable type                                 |

Note: ● Recommended configuration ○ Optional

## R530 Portable Anesthesia Machine

Regular assembly (R540-A), vaporizer (R580), Inlet/Outlet End Caps (R580-IO), induction chamber, fixed bracket for gas evacuation apparatus, air pump, and portable assembly (R510-S), tubing are included in the R520 anesthesia machine.

Different air supply, induction chambers, small animal masks, filter canisters and gas evacuation apparatus should be purchased separately. The difference between R520 & R530 is only the anesthesia assembly.



★ The circled item should be purchased separately

### Order Information:

| Model  | Product Description                                       | Configuration | Remarks   |
|--------|---|---------------|---|
| R530IE | Portable Anesthesia Machine, Isoflurane, Easy Fill, 0-4L  | ●             | Including fixed bracket for induction chamber, evacuation apparatus, filter canister. |
| R530SE | Portable Anesthesia Machine, Sevoflurane, Easy Fill, 0-4L | ●             |   |

Note: ● Recommended configuration ○ Optional

## Veterinary Anesthesia Machines ▶

### R620-S1 Universal Veterinary Anesthesia Machine

R620-S1 universal veterinary anesthesia machine is applicable to clinical anesthesia of dogs, cats, pigs, monkeys, rodent, reptiles, birds etc. With "one-click" operation mode, it would be easy, safe and reliable to perform surgery on pets.

#### User-friendly design

Blocking bar design for breathing circuit, wheel brake for Universal wheel, and tray shelf for placing monitor

#### Rail-mounted CO<sub>2</sub> absorber

One step to assemble and disassemble, save time and effort to a great extent


#### Full inspection policy

Inspection would be made before shipment to make sure high accuracy

#### Minimal operation design

Just press to close APL valve; turn the switch to shift between rebreathing and non-rebreathing circuit; with only one action to install and uninstall the vaporizer; With one switch to turn, the smart weighing canister start to work

| Item                     | Specification  |
|--------------------------|--|
| Vaporizer                | Output concentration: 0~5%(Isoflurane);<br>0~8%(Sevoflurane)<br>Accuracy: ±0.1(0-1%); ±0.15(>1%)<br>work temperature: 10~35°<br>Volume>120ml |
| Flowmeter                | 0.1-4L/min   |
| CO <sub>2</sub> Absorber | 1.8L   |
| Filter canister          | Self-weight is 800g, Maximum<br>absorption capacity is 200g  |
| Weighing device          | Alarm range>1000g<br>(R510-31 filter canister);<br>Resolution:1g   |



Optimal: MRI Conditional Animal Anesthesia Machine are available for purchase.

#### Order Information:

| Model        | Product Discription   | Configuration | Remarks  |
|--------------|---|---------------|--|
| R620-S1-IECS | Veterinary Anesthesia Machine-Pole Mount, Isoflurane, Easy Fill, Cagemount  | ●             | Super package. Comparedwith the basic package, an additional oxygen valve,calcium lime,a trachea cannula hanger ,Mounting Braket for Weighing Meter, and an automatic Weighing Apparatus are included. |
| R620-S1-SECS | Veterinary Anesthesia Machine-Pole Mount, Sevoflurane, Easy Fill, Cagemount | ●             |  |
| R620-S1-IECP | Veterinary Anesthesia Machine-Pole Mount, Isoflurane, Easy Fill, Cagemount  | ●             |  |
| R620-S1-SECP | Veterinary Anesthesia Machine-Pole Mount, Sevoflurane, Easy Fill, Cagemount | ●             |  |
| R620-MRI     | MRI Conditional Animal Anesthesia Machine                                   | ●             |  |

Note: ● Recommended configuration ○ Optional


## Complete Anesthesia and Respiration Solutions ▶

### Multi-function Anesthesia Solutions



The experimental system is equipped with a gas recovery system to provide a safe experimental environment to maximize the protection of the operator; by replacing different types of anesthesia machines and masks to meet the needs of animal experiments, such as cerebral ischemia and myocardial ischemia model, osteoporosis model, tail vein injection, abdominal aorta blood sampling, heart blood sampling, ventricular injection, tissue / organ removal and other animal surgery experiments.

### Stereotaxic Anesthesia Solutions



By changing the different types of stereotaxic instruments and adaptors to meet the needs of different animal experiments, such as the neurological disease models (such as Parkinson, Alzheimer's disease, craniocerebral injury, spinal cord injury, etc.), brain drug injection, Cannula Implantation (long term repeated administration), physiology experiment (such as nerve stimulation, signal recording, wireless telemetry, etc.), microdialysis experiments (probe implantation and short-term dialysis) and other experiments.

## Combination Solutions of Ventilator to Anesthesia System

### Small animals



This program is the combination of anesthesia machine and ventilator. While keeping the animal under anesthesia, it provides a certain volume (tidal volume) of oxygen and respiratory rate to support the normal breathing of animals.

It can meet the requirements of long surgeries or surgeries on animals suffering from certain diseases (such as cardiovascular and cerebrovascular diseases, hypertension, cerebral infarction, thrombosis), it's suitable for mice, rats, rabbits, cats and other animals within 1kg.

### Large animals



Circular anaesthetic mask

Laryngoscopes

Oxygen Concentrator

Endotracheal Tubes

Calcium lime

This program is suitable for the anesthesia and auxiliary respiration of animals within 100kg, such as pigs, monkeys, dogs and cats.

## Anesthesia Machine Accessories ▶

### (Anesthesia Machine) Vaporizer

Adopted advanced design, R580 series vaporizer can be perfectly collaborated with small animal anesthesia machine and large animal anesthesia machine to make sure accurate transmission of anesthetic gas. RWD offers Pour Fill, Easy Fill and Key Fill Cagemont and Selectatec.

Two kinds of vaporizers  
Three dosing methods  
Two installation



Isoflurane vaporizer



Sevoflurane vaporizer

#### Features:

- The built-in temperature compensator ensures the stable concentration of anesthetic gas at different temperatures and flow rates, the flow range is 0.2-10L/min.
- Mixed anesthetic gas output pressure range  $\Delta P \leq 2.5kPa$ .
- Applicable temperature range: 10-35 °C.
- Good tightness, internal 50kPa pressure can be maintained, zero leakage, safe and reliable.
- Internal capacity is not less than 120ml, and the range of isoflurane vaporizer concentration is 0-5% (Sevoflurane: 0-8%).
- Anti-accident open lock structure and closed state security protection structure to ensure safety.
- The average filling rate is greater than 2ml/s, and the dosing adapter can withstand the pressure of 0.1mpa without leakage (Key Fill).



Pour-Fill



Easy-Fill



Key-Fill

#### Note:

When purchasing a Cagemont vaporizer, vaporizer endcaps (model: R580-IO) should be purchased separately.  
Vaporizer of Easy Fill includes adapter R582-01 or R583-01.  
Vaporizer of Key Fill includes adapter R584-01 or R585-01.

#### Order Information:

| Cat No. | Product Description                        |
|---------|--|
| R580S   | Vaporizer - Isoflurane/Pour Fill/Cagemont  |
| R581S   | Vaporizer - Sevoflurane/Pour Fill/Cagemont |
| R582S   | Vaporizer - Sevoflurane/Easy Fill/Cagemont |
| R583S   | Vaporizer - Isoflurane/Easy Fill/Cagemont  |

| Cat No. | Product Description                     |
|---------|---|
| R584S   | Vaporizer-Isoflurane/Key Fill/Cagemont  |
| R585S   | Vaporizer-Sevoflurane/Key Fill/Cagemont |
| R580-IO | Vaporizer Endcaps (1 pair)              |

### Isoflurane & Sevoflurane

#### Order Information:

| Model      | Product Description               |
|------------|-----------------------------------|
| R510-22-16 | Isoflurane-100mL (16 bottles/box) |
| R510-22-8  | Isoflurane-100mL (8 bottles/box)  |



Isoflurane

## Gas Evacuation Apparatus

The gas evacuation apparatus (GEA) was designed to safely and actively remove waste gas associated with everyday inhalation anesthesia. The waste gas from cone mask, concentric tubing mask and stereotaxic frame nosecone masks can be pulled through activated charcoal canister so that the air can be re-circulated into the room. The GEA will accommodate waste gas supply from at least 5 masks. The new one has been upgraded with a wider range of exhaust flow, Low noise, and Exquisite exterior design.

### Features:

- New design make it more exquisite and easy to carry.
- A wide range of flow rate-8~60LPM, real-time display through LED screen.
- Small size 232mm (W)\*200mm (D)\*184mm (H).
- Noise reduction system within 50dB. One click set mute to prevent noise.
- Two level overweight alarm: 990g and 1010g.
- Power: 100-240V, 50/60Hz.
- Automatic temperature compensation system can help stably run even in 5°C-40°C.
- Integrated negative pressure recovery system and weighing alarm functions, can weigh and display the weight of the gas filter tank at anytime to confirm whether the adsorption is saturated.



R546Pro

### Order Information:

| Model   | Product Description                                       |
|---------|---|
| R546Pro | Gas Evacuation Apparatus with Weight Monitoring, 110-240V |

## Gas Filter Canister



Used to absorb isoflurane, sevoflurane and enflurane and other gases, disposable, large and small size are available, large canister (net weight: 800g) absorb 200-240g gas; small canister (net weight: 250g) absorb 50g gas.

### Order Information:

| Model       | Product Description                      |
|-------------|--|
| R510-31-6   | Gas Filter Canister, Large, pkg of 6pcs  |
| R510-31S-6  | Gas Filter Canister, Small, pkg of 6pcs  |
| R510-31S-16 | Gas Filter Canister, Small, pkg of 16pcs |

## Induction Chambers

The anesthetic gas is inducted into the chamber, and exhausted to the connected charcoal canister. With following features:

- Transparent plexiglass material, durable, easy to clean, convenient to observe the anesthetic state of animals.
- Positive pressure gaskets at lids.
- Both ports on one side to reduce tubing in procedure area.



### Order Information:

| Model | Product Description                             | Remarks  |
|-------|---|--|
| V100  | Induction Chamber-mouse (15cm*10cm*10cm)        | Contains 1 plexiglass box, 1.2m bellows and 1 gas filter canister (R510-31S) |
| V101  | Induction Chamber-rat (24cm*12cm*18cm)          |  |
| V102  | Induction Chamber-rabbit/cat (40cm*18.5cm*25cm) |  |
| V103  | Induction Chamber-large animal (50cm*30cm*30cm) |  |

## Small Animal Anesthesia Mask

### • Concentric Tubing Masks

Using Concentric Tubing design, sturdy and durable, can be placed arbitrarily on the desktop. The entire mask (including tubing) is made of non-metallic material suitable for MRI environments.

- The mask can be secured by fixed bracket (68625), height of the mask is adjustable, and it can also be fixed and placed in MRI, PET and SPECT imaging systems.



### Order Information:

| Model   | Product Description                                | Model   | Product Description   |
|---------|--|---------|---|
| R510-86 | SGL MRI Compatible Mask for Mouse                  | R510-87 | SGL MRI Compatible Mask for Rat                                   |
| R510-88 | Dual MRI Compatible Masks for Mouse                | R510-89 | Dual MRI Compatible Masks for Rat                                 |
| R510-54 | Five-channel Tubing Anesthesia Mask Kits-rat/mouse | 68625   | MRI Compatible Mounting Bracket for Tubing or Cone Mask, optional |

### • Cone Masks

- Special designed with inlet tubing inside make the animal anesthesia exhaust gas can be absorbed by gas filter canister.
- Allows access to eyes, ears and top of head.
- The cone masks can be installed onto 68630 operation bed for single cone mask (300\*210\*75mm) or 68657 manifold for five cone masks.



### Order Information:

| Model | Product Description                                     | O.D. (size) | I.D.(size) | D(size) | Nosecone Model |
|-------|---|-------------|------------|---------|----------------|
| 68680 | Cone Mask with Tubing for Neonatal Mouse(<15g)          | 13mm        | 5mm        | 24mm    | 68683          |
| 68681 | Cone Mask with Tubing for Mouse or Neonatal Rat(15-40g) | 20mm        | 10mm       | 28.5mm  | 68684          |
| 68682 | Cone Mask with Tubing for Rat(200-350g)                 | 25mm        | 15mm       | 38mm    | 68685          |

| Model     | Product Description                            | Model | Product Description           |
|-----------|--|-------|-------------------------------|
| R-CT-6/15 | Concentric Tubing, 6mm Inlet/15mm Outlet, 45cm | 68657 | Manifold for Five Cone Masks  |
| 68656     | Manifold for Single Cone Mask                  | 68630 | Anesthesia Operation Platform |

# Animal Ventilator

## VentStar Small Animal Ventilator

Animal ventilator is the most commonly used equipment to assist animal respiration. It is widely used in basic medicine, clinical medicine, animal medicine and other scientific research experiments (such as myocardial ischemia model making, lung function evaluation), respiratory management, animal first aid and respiratory treatment in clinical surgery. Small animal ventilator suitable for 10g-1kg mice, rats and guinea pigs and other small animal.

Suitable for mice, rats and guinea pigs and other animals weighing 10g-1kg

Input power range: 100 ~ 240V, 50/60Hz, 24V, 40W, 1.67A adapter DC power supply, safe and effective

Small size, compact layout, flexible operation

Provide safe and effective intermittent positive pressure ventilation (IPPV), volume and pressure modes are available.

7 inches LCD touch screen, wide-angle visual

The pressure-time graph can be displayed in real time in volume and pressure modes

Sound, intelligent text, error code and other fault alarm alerts

10 groups of different respiratory parameters can be stored and called, which is convenient and efficient

R415-AP

Set the respiratory parameters such as PIP/PEEP/Sigh/INSP.Hold/EXP. Hold

After sumitting the animal weight, and the device automatically sets the reference breathing parameters

The adjustable range of I:E ratio is from 20% to 80%; the sigh increases the tidal volume from 0 to 20%, and can set auto-sigh once after every 10 to 999 breaths or manual control

The breathing frequency can be switched between 10 ~ 300bpm

Tidal volume range: 0.05ml ~ 5ml

Airway pressure range: 1~50cmH<sub>2</sub>O, the accuracy range: ±0.7 cmH<sub>2</sub>O; positive end expiratory pressure (PEEP) range is 0~10cmH<sub>2</sub>O

| Model   | Product Description                             |
|---------|---|
| R415    | Small Animal Ventilator                         |
| R415-AP | Accessory package for Anesthesia and Ventilator |

## RE20 Animal Endotracheal Intubation Kit

Noninvasive endotracheal intubation kit for animals guided by cordless optic fiber to support mouse, rat and other similar small rodents in a stable and comfortable position for experimenters to perform endotracheal intubation, drug infusion and other operations, effectively reducing animal damage, and enhancing the experimental efficiency and precision of the experiment results.

### Features:

- PU material, moderate soft and hard, comfortable and not harmful to animals;
- The cordless portable optic fiber lighting system, which has been patented for an invention, improves the visibility of the airway and facilitates users to quickly and accurately finish endotracheal intubation;
- The stand angle can be adjusted, including: 45°, 60°, 90°, to meet the requirements of users with different operation habits;
- An endotracheal magnifier with auxiliary positioning function and a tracheoscope mirror that helps to confirm intubation result are more user-friendly.




### Order Information:

| Model  | Product Description                                  |
|--------|--|
| RE20-R | Endotracheal intubation tube for rat (length 50mm)   |
| RE20-M | Endotracheal intubation tube for mouse (length 25mm) |

## R419 Intelligent Animal Ventilator

R419 intelligent animal ventilator as mechanical ventilation device, which is completely according with veterinary clinical practice. There are 3 modes of operation, including VCV, PIPCV and APNEA, to meet multiple clinical applications.



- Adjustable power supply (110 / 220V)
- Light weight, portable, user-friendly design, multi-parameter fixed platform, easy to operate
- Provides safe and effective intermittent positive pressure ventilation (IPPV)
- Provide internal battery power, which can work more than 4 hours

### Order Information:

| Model   | Product Description                  |
|---------|--------------------------------------|
| R419    | R419 Intelligent animal ventilator   |
| R409-MB | Mounting Bracket for R409 Ventilator |

## Surgical Auxiliary Equipment

### DOM-1001 Digital Microscope

DOM-1001 is a digital microscope without eyepiece, with dual light sources and a cross arm bracket. High-end optical components and 16 million pixel CMOS camera can be used to capture real-time images and video, realizing multi-person sharing. Mainly used in life science, university teaching, industrial manufacturing, material science and other fields.



#### Features:

- **Novel and unique integrated design**

High-end optical components and a 16 million pixels digital camera, combined with LED ring lights and gooseneck light control components, and optional brackets for various scenarios, together create a powerful DOM-1001 digital surgical microscope system.

- **Excellent ergonomic design**

The comfortable electric adjustment knob on the side is combined with the foot focus adjustment, and with the 11.6-inch display, the operator can maintain a correct and comfortable working posture, reducing fatigue and damage caused by long-term use.

- **Real-time ultra high-definition projection screen**

Built-in digital camera with zoom optical components, high depth of field (maximum depth of field 6.7mm), adjustable magnification 8-44X(11.6-inch display), 18-100X(27-inch display), support real-time display with high-definition monitor or notebook computer connection.

- **Convenient and practical remote control**

With the infrared remote control, you can directly perform remote control, including many functions such as taking pictures, recording videos, flipping and zooming, adjusting brightness and contrast of images, and white balance.

- **Homogeneous and coaxial optical design**


Always displayed in the center of the field of view in a clear and optimal state-no need to refocus even when changing the magnification.

- **Greater operating space**

The working distance of 180 mm, the lens body can be adjusted at multiple angles, which greatly improves the operating space and ensures that observation and operation are performed simultaneously.

## Microscope

77001S Desktop Binocular Stereo Microscope equips with universal support; researchers can observe in high resolution and three-dimension.



77001S

### Features:

| NO. | Products Name       | Products Detail  |
|-----|---------------------|--|
| 1   | Eyeiece             | 10X Eyepieces, View field diameter $\phi$ 22mm   |
| 2   | Objective Lens      | 0.67X-4.5X Objective lens ratio 1:6.7  |
| 3   | Binocular Drawtube  | Pupil Distance 55-77mm Diopter adjustment $\pm$ 5 ,Diopter 45°Incline 360°Spin               |
| 4   | View field Range    | 31.2 mm -5.1mm   |
| 5   | Total Magnification | 6.7X-45X(10X Eyepiece),3.35X-22.5X(10X Eyepiece +0.5X Assist lens)                           |
| 6   | Working Distance    | Working distance 105mm(10X Eyepiece),Working distance 177mm (10X Eyepiece +0.5X Assist lens) |
| 7   | Focusing Mechanism  | Focusing hand wheel is adjustable  |
| 8   | Universal Support   | 360°Spin, Retractable  |
| 9   | Lighting            | TY-B-2 Ring LED Illumination   |

On the basis of the binocular microscope, 77008S Trinocular Zoom Stereo microscope increase the interface, so that the camera, video camera can be installed in the microscope to achieve real time observing of the objects.



77008S

## 77002 Vertical Animal Operating Microscope

The vertical animal operating microscope is a simple operating microscope with single person and binocular coaxial illumination. The body is small and light with high flexibility, which can meet the requirements of general animal microsurgery.

### Features:

- Eyepiece magnification: 12.5x, objective lens focal length: 200 mm.
- Working distance: 190mm, main mirror magnification: 5.3x, 8x, 12x.
- Field of view diameter: 38mm, 25mm, 17mm.
- The adjustment range of diopter:  $\pm$  5D, the adjustment range of pupil distance: 50mm-70mm.
- Lighting source: 12V / 100W, cold reflection medical halogen bulb.
- Lighting type: 6° + 0° cold source coaxial lighting.
- Illumination of coaxial lighting surface:  $\geq$  20000 LX, transverse arm extension radius: 870 mm.
- Vertical adjustment range (ground to large objective): 700mm-1100mm, fine focus stroke: 30mm.
- Voltage: AC220 V  $\pm$  10%, 50 Hz  $\pm$  1 Hz, power: 120 VA.
- Fuse: AC250V T1.25A weight: 41kg.



77002

## Animal Restrainer

### Tubular type Restrainer

Tubular type restrainer is used to restrict the animal from turning or moving, by fixation on rear-end and flexibly regulated. Unique design of tapered nose restrainer provide easy-access for animals, meanwhile makes animals feel safety.

Taper restrainer is a new fixed tool design for animal injection, suitable for fast and efficient batch injection.



Restrainer



Taper restrainer

### Ordering Information:

| Model | Products Description                                    |
|-------|---|
| 9131  | Small restrainer for up to 15g Mouse/L110xI.D.20mm      |
| 9132  | Medium restrainer for 15-40g Mouse/L135xI.D.30mm        |
| 9133  | Large restrainer for 50-75g Mouse/L160xI.D.40mm         |
| 9134  | Small restrainer for 75-200g Rat/L200xI.D.50mm          |
| 9135  | Medium restrainer for 200-300g Rat/L210xI.D.60mm        |
| 9136  | Large restrainer for 300-500g Rat/L220xI.D.70mm         |
| 9137  | Extra Large restrainer for above 500g Rat/L230xI.D.80mm |
| 9138  | Taper restrainer for Mice/L150xH130xO.D.100mm           |

## Rats & Mice Shaver

### Features:


- Applicable to shave off hair on animals' cephalosome, face, legs and ear.
- AA batteries use.
- Low vibration and noise as well as ergonomic design.
- Suitable for small animals.

CP-5200



## Manual Syringe Pumps(68606)

Compatible with RWD stereotaxic instruments, simple installation, with 10µm displacement accuracy and up to 25mm distance, it is an economical and practical microinjection pump!



### Features:

- Precisely-controlled infusion/withdrawal of picoliter volume.
- Provides an economical means of injecting or withdrawing a calibrated volume with simple, precise manual control.
- The injector has a universal adaptor for all RWD and Hamilton 700, 7000 or Gas-Tight Series syringers, ranging from 1 to 100µl sizes and enables a wide range of volumes to be delivered-both infusion and withdrawal operation.
- The Manual Stereotaxic injector is graduated in 0.01mm increments. One complete revolution advances the plunger of the syringe 0.5mm.

## Cordless Microdrill

It is battery-powered, easy to carry and space-saving.

### Features:

- Cordless design, simple & controllable operation.
- The double-bearing design is used at the joint of the drill bit, make drilling and grinding more stable.
- 750mA lithium battery, can be used continuously for 45 mins after fully charged.
- Wide range of applications, can be used with a 2.3mm diameter drill bit for drilling and polishing.
- It can be equipped with a cranial drill holder (68228) fixed to Stereotaxic, free hands operation.
- Speed: 12000rpm.
- Stall torque: 1.25 oz./in.
- Power: 100-240V.
- Dimension: length 20cm, body diameter 19.5mm .



### Order Information:

| Model | Product Description  |
|-------|--|
| 87001 | Cordless Microdrill (Bits Not Included, some details see below)              |
| 68228 | Cordless Microdrill Holder(Holds the microdrill with diameter about 19.5mm.) |

## Microdrill

It takes only one second to drill the skull, easy to operate and control. It can also be used for grinding the skull.



### Features:

- Speed up to 38,000rpm.
- Forward or reverse available.
- Controlled manually or by foot pedal.
- Can be installed onto the stereotaxic instrument by the microdrill holder (Model: 68605) and controlled by the manipulator.
- Various drill bits available, diameter from 0.5-2.3mm.

### Order Information:

| Cat No. | Product Description                          | Cat No. | Product Description                          |
|---------|--|---------|--|
| 78001   | Microdrill (bits not included)               | 78050   | Drill Bits HM1027 2.7mm, Round Tip, pkg of 2 |
| 78040   | Drill Bits HM1005 0.5mm, Round Tip, pkg of 5 | 78060   | Drill Bits HM31008 0.8mm, Flat Tip, pkg of 5 |
| 78041   | Drill Bits HM1006 0.6mm, Round Tip, pkg of 5 | 78061   | Drill Bits HM31010 1.0mm, Flat Tip, pkg of 5 |
| 78042   | Drill Bits HM1008 0.8mm, Round Tip, pkg of 5 | 78062   | Drill Bits HM31012 1.2mm, Flat Tip, pkg of 5 |
| 78043   | Drill Bits HM1010 1.0mm, Round Tip, pkg of 5 | 78063   | Drill Bits HM31014 1.4mm, Flat Tip, pkg of 5 |
| 78044   | Drill Bits HM1012 1.2mm, Round Tip, pkg of 5 | 78064   | Drill Bits HM31016 1.6mm, Flat Tip, pkg of 5 |
| 78045   | Drill Bits HM1014 1.4mm, Round Tip, pkg of 5 | 78065   | Drill Bits HM31018 1.8mm, Flat Tip, pkg of 5 |
| 78046   | Drill Bits HM1016 1.6mm, Round Tip, pkg of 5 | 78066   | Drill Bits HM31021 2.1mm, Flat Tip, pkg of 5 |
| 78047   | Drill Bits HM1018 1.8mm, Round Tip, pkg of 5 | 78067   | Drill Bits HM31023 2.3mm, Flat Tip, pkg of 5 |
| 78048   | Drill Bits HM1021 2.1mm, Round Tip, pkg of 5 | 68605   | Microdrill Holder                            |
| 78049   | Drill Bits HM1023 2.3mm, Round Tip, pkg of 5 |         |  |

## Fiber-lite Illuminator

The halogen fiber-lite illuminator adopts dual-pipe pipe optic coaxial lighting, provide better field of view, with good brightness; stable output voltage, safe and reliable light bulb; size is about 280 x, 100 x 165mm, and the light intensity can be adjusted continuously , with color temperature 3200K.



Halogen fiber-lite illuminator

### Order Information:

| Model | Product Description                                      | Remarks                                     |
|-------|--|---|
| 76301 | Fiber-Lite Illuminator(150W)                             | Dual-Pipe Fiber Optic, Bulb and Main Engine |
| 76304 | Main Engine  |   |
| 76310 | Dual-Pipe Fiber Optic(Included in illuminator)           |   |
| 76306 | Replacement Bulb for 76301 Illuminator (21V 150W) -OSRAM |   |

## ThermoStar Homeothermic Monitoring System

During the process of animal surgery or acute anesthetized animal experiment, maintaining animals' body temperature is very important for the accuracy and reliability of results, furthermore, it is an important way to reduce the mortality of the animals. In order to fulfill this need, RWD launched a new generation of homeothermic monitoring system "Thermostat", the control unit acquired a 4.3 inch color LCD touch screen, allows users to clearly view the subject's real-time body temperature, set the target temperature, view the heating state and set audible alarms when subject's body temperature turns abnormal. A specific control separate circuit for heating pad and rectal probe, respectively to ensure the stability and safety of the system.



heating pad (69023 / 69024 / 69025)



temperature probe (69022)

### Features:

- Real-time display of the working status of the heating pad, and monitor the working state of every channel, show: IDLE, WARNING, STANDBY, and WORKING.
- Heating range: 20-45°C, 0.1°C resolution, Degrees Celsius and Fahrenheit can be switched.
- Power supply: 100-240V, 50/60Hz.
- Control unit size: L205 x W148 x H80 mm.
- Control unit weight: 1.3kg.

### Order Information:

| Type              | Model | Product Description                                | Remarks  |
|-------------------|-------|--|--|
| Complete system   | 69020 | ThermoStar Temperature Controller, 69023 pad incl. | A heating pad (69023) and a temperature probe (69022) incl.        |
|                   | 69026 | ThermoStar Temperature Controller, 69024 pad incl. | A heating pad (69024) and a temperature probe (69022) incl.        |
|                   | 69027 | ThermoStar Temperature Controller, 69025 pad incl. | A heating pad (69025) and a temperature probe (69022) incl.        |
| Temperature probe | 69022 | Temperature Probe                                  | 1.5mm diameter, 30mm rod length, the length of the probe is 150cm. |
| Heating pad       | 69023 | Heating Pad-12.0*20.5cm                            | the length of heating pad is 150cm.                                |
|                   | 69024 | Heating Pad-9.0*17.0cm                             |  |
|                   | 69025 | Heating Pad-7.0*10.0cm                             |  |

## RLT-24 Veterinary Laser Therapy

RLT-24 Veterinary Laser Therapy is non-invasive, deep penetrating, painless and without side effects, and can be used for recalcitrant wounds, chronic diabetic ulcers, neurological injuries. It can be used in the study of diseases such as intractable injuries, chronic diabetic ulcers, neurological damage and inflammatory pain.



### Features:

- Intelligence Solutions: more than 1,300 recommended treatment built in, it can be adjusted according to the specific experiment.
- Multiple modes: support three laser emission modes: CW mode, Pulse mode and ISP mode. The combination of multiple modes can improve the treatment efficiency, reduce thermal burns, and intensify laser treatment.
- High power: the laser power is up to 24W, with remarkable treatment effect, stable light, safe and fast.
- Four wavelengths: arbitrary combination of four wavelength lasers of 650nm, 810nm, 915nm and 980nm, different wavelengths of laser can be selected for different experiments.
- User-friendly design: LCD touch screen; the rear-mounted large optical fibre winding disk, which is convenient for storage and protects the optical fibre and prolongs its service life.
- Features: built-in video tutorial board, "hands-on" teaching, and no difficulty in operation; automatic storage of treatment records for easy information traceability.

### Order Information:

| Cat No. | Description              |
|---------|--------------------------|
| RLT-24  | Veterinary Laser Therapy |

### Included:

| Description              | Amount | Function  | Remarks  |
|--------------------------|--------|---|----------|
| Veterinary Laser Therapy | 1      | Laser therapy   | Included |
| Optical fibre handpiece  | 1      | Laser transmission and control  |          |
| Lens                     | 5      | Change the beam size to realize optimal treatment and diversity   |          |
| Protective goggles       | 6      | 3 for staff and 3 for animals for eyes protection   |          |
| Interlock                | 1      | The interlock should be connected to the device to prevent any personnel entering the room and being exposed to the laser irradiation |          |
| Power adaptor            | 1      | Device power input  |          |
| Toolkit                  | 1      | Store lens and optical fibre handpiece  |          |

## R-CU600 Smart Animal ICU

R-CU600 The Smart Animal ICU, provides animal patients with the most comfortable medical care environment to help animals recover faster. It is suitable for postoperative recovery, emergency care of critical animals, delivery, infant care, respiratory diseases, aging diseases, animal atomization therapy and other animal diseases treatment.



### Features:

- **Smart parameters:** a variety of model animal parameters are built in the system, which can be selected according to the actual needs of the experiment;
- **Touch screen operation:** high definition 7-inch touch screen, intuitive, parameters at a glance;
- **Trebling purification function:** air filter / UV sterilization system /high concentration negative ion environmental, avoid cross infection of animals;
- **Multiple alarm:** the system covers temperature and humidity, CO<sub>2</sub> concentration, sensor failure, heating failure and other abnormal alarms;
- **Digital control:** inner wall wrapped 3D temperature control system, heating balance; Independent external humidifier, stable performance, convenient injection; Sensitive infrared probe, accurate monitoring of CO<sub>2</sub> concentration;
- **Satisfy a variety of experimental conditions:** infusion / oxygen supply / atomization therapy.

### Order Information:

| Model   | Product Description |
|---------|---------------------|
| R-CU600 | Smart Animal ICU    |

## Surgical Instruments & Matrix

RWD surgical instruments and Matrices products cover more than 10 types of instruments and consumables such as surgical scissors, clamps, forceps, auxiliary instruments, sterilizers, molds, sutures. There are thousands of specifications as well offer well-equipped surgical kits which could satisfied with all kind of animal modeling experiment. More details please refer to the catalog from page 66 to 81, or visit our web: www.rwdstco.com.


### The surgical kits included:

Surgical instrument kit for MCAO, spinal cord injury, universal dissecting, cardio-cerebral perfusion surgery, brain micro injection surgery and micro dialysis surgery etc. which could satisfied with the whole experiment operation for rodents from skin dissection to sutured.




### Surgical scissors included:

Micro spring scissors, delicate scissors, surgical scissors, dissecting scissors, boon cutters and ophthalmic scissors etc. which could satisfied with different experiment requirement of tissues dissecting, cutting, and separate.



### Surgical forceps included:

Tissue forceps, hemostatic forceps, needle holder and bone rongeurs etc. which could satisfied with the requirement of tissues clamping, fixation, hemostasis, stitich and bone trimming and cutting.



### Surgical forceps included:

Micro forceps, tissue forceps and dressing forceps and ophthalmic forceps etc. it is mainly used for the clamping and separation of soft tissues such as blood vessels, nerves, organs, skin etc.



### Auxiliary instruments such as:

Scalpel, retractors, hooks, blood vessel clamps, sutures <with needle>, sterilization box and marking instrument etc. to meet the requirement of skin dissection.



### The matrices included:

Brain and hear matrix, stainless material offered solid and durable quality which could be heated, cooled, autoclaved and scrubbed for long time using. Sterilizer mainly used for Instruments' high temperature sterilization with high efficiency. MCAO sutures made of Silica gel wrapped nylon thread material which could easily enter the common carotid artery and embolize the middle cerebral artery without puncturing the inner wall of the vessel, which greatly improves the success rate and repeatability of the modeling.



# Neural Signal Regulation Solution

## Optogenetics System

Optogenetics uses light-responsive proteins to control biological processes. The unmatched spatiotemporal precision, high tunability, reversibility, and few side effects of light makes optogenetics a powerful strategy to study neuroscience and complex biological systems.

RWD new intelligent optogenetics system, on the basis of the last generation machines (laser and waveform generator integrated into one), has made another leap in the performance.



### Features:

- Built-in CW(continuous wave)mode, one-button power test.
- The test/experiment mode can be switched with one key, and the test power value can be directly written into the experiment program.
- Equipped with a remote control, which can control multiple optogenetics devices simultaneously within a range of ten meters.
- Three external trigger modes: Edge/Real-time/Gate, to meet more experimental requirements.
  - Edge mode1:** If the trigger voltage is higher than or equal to the set value of trigger voltage, the program will run; the trigger is only featured by the "on" function and the program's duration depends on its set duration.
  - Real-time mode2:** If the trigger voltage is higher than or equal to the set value of the trigger voltage, the program will run; if it is lower than that value, the program will stop; the trigger is characterized by the "on/off" function. The duration of the program depends on the duration of the trigger.
  - Gate mode3:** The running power relies on the first subprogram selected and other parameters are defined by the external device, and only response to the square wave with frequencies  $\leq 500\text{Hz}$ .
- High resolution signal processing, external signal delay  $< 1\text{ms}$ .
- High stability, long-term operation power fluctuation is less than 1%.
- Support program group custom editing, import and export, experimental programs and records can be exported externally to review.
- With signal input/output interface, it can be compatible with electrophysiology, fiber photometry system, behavioral instruments and other equipment.

### Software Interface:



Parameter Setting Interface



Program Running Interface



Trigger Interface

**Technical Parameters:**

| Parameter item | Specification  |
|----------------|--|
| Power          | 1-100mW, 1mW resolution                              |
| Operation mode | Remote control/touch screen/external trigger         |
| Input          | 2V-10V   |
| Output         | 10mV-5V  |
| Safety         | Password; High temperature alarm; Program run prompt |
| Storage        | Internal memory/USB flash disk                       |

**Order Information:**

| Model   | Product Description                    |
|---------|--|
| IOS-465 | Intelligent Optogenetics System(465nm) |
| IOS-589 | Intelligent Optogenetics System(589nm) |

**Patchcord and Optical Fiber**

- Low insertion loss • Good repeatability • High return loss • Stable temperature
- Good mutual insertion performance • Optional fiber
- Applicable wavelength: 400nm-1100nm • Numerical aperture: 0.22, 0.39,0.50
- Default length is 1m, other length:1.5-5m, step 0.5m (customized)
- Core diameter: 100 $\mu$ m, 200 $\mu$ m, 300 $\mu$ m, 400 $\mu$ m

**Order Information:**

| Model             | Product Description   |
|-------------------|---|
| R-FC-PC-N2-100-L1 | FC/PC Optical Fiber Patch Cables(m) - 100 $\mu$ m,0.22NA,1m                                   |
| R-FC-PC-N2-200-L1 | FC/PC Optical Fiber Patch Cables(m) - 200 $\mu$ m,0.22NA,1m                                   |
| R-FC-PC-N2-300-L1 | FC/PC Optical Fiber Patch Cables(m) - 300 $\mu$ m,0.22NA,1m                                   |
| R-FC-PC-N2-400-L1 | FC/PC Optical Fiber Patch Cables(m) - 400 $\mu$ m,0.22NA,1m                                   |
| R-FC-PC-N3-200-L1 | FC/PC Optical Fiber Patch Cables(m) - 200 $\mu$ m,0.39NA,1m                                   |
| R-FC-PC-N3-300-L1 | FC/PC Optical Fiber Patch Cables(m) - 300 $\mu$ m,0.39NA,1m                                   |
| R-FC-PC-N3-400-L1 | FC/PC Optical Fiber Patch Cables(m) - 400 $\mu$ m,0.39NA,1m                                   |
| R-FC-PC-N5-200-L1 | FC/PC Optical Fiber Patch Cables(m) - 200 $\mu$ m,0.50NA,1m                                   |
| R-FC-PC-N5-400-L1 | FC/PC Optical Fiber Patch Cables(m) - 400 $\mu$ m,0.50NA,1m                                   |
| R-FC-L-N2-100-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 100 $\mu$ m Core,0.22NA,1m |
| R-FC-L-N2-200-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 200 $\mu$ m Core,0.22NA,1m |
| R-FC-L-N2-300-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 300 $\mu$ m Core,0.22NA,1m |
| R-FC-L-N2-400-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 400 $\mu$ m Core,0.22NA,1m |
| R-FC-P-N2-100-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 100 $\mu$ m Core,0.22NA,1m  |
| R-FC-P-N2-200-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 200 $\mu$ m Core,0.22NA,1m  |
| R-FC-P-N2-300-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 300 $\mu$ m Core,0.22NA,1m  |
| R-FC-P-N2-400-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 400 $\mu$ m Core,0.22NA,1m  |
| R-FC-L-N3-200-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 200 $\mu$ m Core,0.39NA,1m |
| R-FC-L-N3-300-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 300 $\mu$ m Core,0.39NA,1m |
| R-FC-L-N3-400-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 400 $\mu$ m Core,0.39NA,1m |
| R-FC-P-N3-200-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 200 $\mu$ m Core,0.39NA,1m  |
| R-FC-P-N3-300-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 300 $\mu$ m Core,0.39NA,1m  |
| R-FC-P-N3-400-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 400 $\mu$ m Core,0.39NA,1m  |
| R-FC-L-N5-200-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 200 $\mu$ m Core,0.50NA,1m |
| R-FC-L-N5-400-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 1.25 mm Ceramic Ferrule, 400 $\mu$ m Core,0.50NA,1m |
| R-FC-P-N5-200-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 200 $\mu$ m Core,0.50NA,1m  |
| R-FC-P-N5-400-L1  | Optic Fiber(m), FC/PC Adapter with $\Phi$ 2.5 mm Ceramic Ferrule, 400 $\mu$ m Core,0.50NA,1m  |



Patchcord



Optical Fiber

**Fiber Optic Rotary Joint**

Applicable to long time stimulation and observation experiments on awake and free moving animals, avoids twining of the fiber optic.


**Order Information:**

| Model    | Product Description            |
|----------|--------------------------------|
| R-FC-1x1 | FC-FC Optical Fiber Commutator |

## Ceramic Ferrule

### Features:

- Material: Ceramic, Optical fiber
- Type: White Series/Black Series
- Applicable wavelength: 400-2200nm
- Core diameter: 100µm, 200µm, 300µm, 400µm
- Numerical aperture: 0.22, 0.39, 0.50
- Outer diameter: 1.25mm (6.4mm long), 2.5mm (10.5mm long)
- Optical fiber length: 2mm-25mm (Ø1.25mm), 2mm-20mm (Ø2.5mm), step 0.5mm (Other lengths can be customized)



### Order Information:

| No. | Model             | Product Description   |
|-----|-------------------|---|
| 1   | R-FOC-L200C-39NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 200um Core, 0.39NA, 20/pkg                              |
| 2   | R-FOC-L300C-39NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 300um Core, 0.39NA, 20/pkg                              |
| 3   | R-FOC-L400C-39NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 400um Core, 0.39NA, 20/pkg                              |
| 4   | R-FOC-F200C-39NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 200um Core, 0.39NA, 20/pkg                               |
| 5   | R-FOC-F300C-39NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 300um Core, 0.39NA, 20/pkg                               |
| 6   | R-FOC-F400C-39NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 400um Core, 0.39NA, 20/pkg                               |
| 7   | R-FOC-L200C-22NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 200um Core, 0.22NA, 20/pkg                              |
| 8   | R-FOC-L300C-22NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 300um Core, 0.22NA, 20/pkg                              |
| 9   | R-FOC-L400C-22NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 400um Core, 0.22NA, 20/pkg                              |
| 10  | R-FOC-F200C-22NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 200um Core, 0.22NA, 20/pkg                               |
| 11  | R-FOC-F300C-22NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 300um Core, 0.22NA, 20/pkg                               |
| 12  | R-FOC-F400C-22NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 400um Core, 0.22NA, 20/pkg                               |
| 13  | R-FOC-L200C-50NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 200um Core, 0.50NA, 20/pkg                              |
| 14  | R-FOC-F200C-50NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 200um Core, 0.50NA, 20/pkg                               |
| 15  | R-FOC-L400C-50NA  | Fiber Optic Cannulae with Φ1.25 mm Ceramic Ferrule, 400um Core, 0.50NA, 20/pkg                              |
| 16  | R-FOC-F400C-50NA  | Fiber Optic Cannulae with Φ2.5 mm Ceramic Ferrule, 400um Core, 0.50NA, 20/pkg                               |
| 17  | R-FOC-BL200C-39NA | Fiber Optic Cannulae(Black) with Φ1.25 mm Ceramic Ferrule, 200um Core, 0.39NA, 20/pkg                       |
| 18  | R-FOC-BL300C-39NA | Fiber Optic Cannulae(Black) with Φ1.25 mm Ceramic Ferrule, 300um Core, 0.39NA, 20/pkg                       |
| 19  | R-FOC-BL400C-39NA | Fiber Optic Cannulae (Black)with Φ1.25 mm Ceramic Ferrule, 400um Core, 0.39NA, 20/pkg                       |
| 20  | R-FOC-BF200C-39NA | Fiber Optic Cannulae (Black)with Φ2.5 mm Ceramic Ferrule, 200um Core, 0.39NA, 20/pkg                        |
| 21  | R-FOC-BF300C-39NA | Fiber Optic Cannulae (Black)with Φ2.5 mm Ceramic Ferrule, 300um Core, 0.39NA, 20/pkg                        |
| 22  | R-FOC-BF400C-39NA | Fiber Optic Cannulae (Black)with Φ2.5 mm Ceramic Ferrule, 400um Core, 0.39NA, 20/pkg                        |
| 23  | R-FOC-BL100C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ1.25 mm Ceramic Ferrule, 100um Core, 0.22NA, 20/pkg |
| 24  | R-FOC-BL200C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ1.25 mm Ceramic Ferrule, 200um Core, 0.22NA, 20/pkg |
| 25  | R-FOC-BL300C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ1.25 mm Ceramic Ferrule, 300um Core, 0.22NA, 20/pkg |
| 26  | R-FOC-BL400C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ1.25 mm Ceramic Ferrule, 400um Core, 0.22NA, 20/pkg |
| 27  | R-FOC-BF100C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ2.5 mm Ceramic Ferrule, 100um Core, 0.22NA, 20/pkg  |
| 28  | R-FOC-BF200C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ2.5 mm Ceramic Ferrule, 200um Core, 0.22NA, 20/pkg  |
| 29  | R-FOC-BF300C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ2.5 mm Ceramic Ferrule, 300um Core, 0.22NA, 20/pkg  |
| 30  | R-FOC-BF400C-22NA | Fiber Optic Cannulae (Black, yellow nylon coating)with Φ2.5 mm Ceramic Ferrule, 400um Core, 0.22NA, 20/pkg  |
| 31  | R-FOC-BL200C-50NA | Fiber Optic Cannulae (Black)with Φ1.25 mm Ceramic Ferrule, 200um Core, 0.50NA, 20/pkg                       |
| 32  | R-FOC-BF200C-50NA | Fiber Optic Cannulae (Black)with Φ2.5 mm Ceramic Ferrule, 200um Core, 0.50NA, 20/pkg                        |
| 33  | R-FOC-BL400C-50NA | Fiber Optic Cannulae (Black)with Φ1.25 mm Ceramic Ferrule, 400um Core, 0.50NA, 20/pkg                       |
| 34  | R-FOC-BF400C-50NA | Fiber Optic Cannulae (Black)with Φ2.5 mm Ceramic Ferrule, 400um Core, 0.50NA, 20/pkg                        |

Notes: The above white series specifications can be customized to black.

## Ceramic Cannula & Ceramic Ferrules Protective Cap



Ceramic Cannula



Ceramic Ferrules Protective Cap

### Order Information:

| Model      | Product Description   |
|------------|---|
| R-MS-1.25  | Ceramic Sleeves, Φ1.25mm (20pcs/ pkg)                       |
| R-MS-2.5   | Ceramic Sleeves, Φ2.5mm (20pcs/ pkg)                        |
| R-BMS-1.25 | Black Ceramic Sleeves, Φ1.25mm (20pcs/ pkg)                 |
| R-BMS-2.5  | Black Ceramic Sleeves, Φ2.5mm (20pcs/ pkg)                  |
| R-DC-2.5   | FC/PC Ceramic Ferrules Protective Cap- Φ2.50mm (200pcs/pkg) |
| R-DC-1.25  | LC/PC Ceramic Ferrules Protective Cap- Φ1.25mm (200pcs/pkg) |

## Fiber Tongs

- Suitable for 100µm-800µm fiber for stripping fiber coating



## Fiber Photometry System ▶

### R810 Dual Color Multichannel Fiber Photometry System

The fiber photometry system records changes in the intensity of fluorescence of neurons in specific brain regions to characterize changes in the activity of neuronal populations. The fiber photometry system can detect population neurons in free-moving animals for long periods of time, which is helpful to explore the correlation between neuron activity and animal behavior. R810 Dual Color Multichannel Fiber Photometry System has two excitation light sources(410nm and 470nm). The 410nm light source can be used to reflect the background noise signal, thus ensuring the acquisition of true fluorescence data. The system can acquire data from no less than 9 channels simultaneously and is suitable for simultaneous recording of multiple brain regions or multiple animals.



#### Features:

##### Convenient and Efficient

- Simply adjust the knob to ensure fiber focus
- The software integrates data acquisition, analysis and data plotting modules
- With a variety of acquisition methods and start & end setting conditions

##### Feature-Rich Extensions

- 4 Input and 4 Output interface for easy connection with other equipment such as optogenetics and electrophysiology for the closed-loop study of stimulation and recording
- Equipped with electrophysiological grounding interface to avoid interference when used in conjunction with electrophysiological experiments
- The software can synchronize and mark multiple special behavioral events or external input signals during the experiment
- The whole system includes the behavior camera, and the software integrates the behavior recording module, which can record the animal movement trajectory and the duration of the ROI area, and supports TTL signal output

##### Data Visualization

- Preview mode facilitates real-time adjustment of parameter settings according to pre-experiment
- Real-time display of signal markers allows direct observation of signal changes
- The analysis results can be exported to CSV, Heat-map,  $df/f$ , Z-Score and other chart formats.

##### High Sensitivity


- Highly sensitive CMOS camera, high frame rate of 300FPS, high quantum conversion efficiency  $\geq 80\%$ , low noise, high gain
- Equipped with photobleaching apparatus to reduce the interference of autofluorescence of accessories

| Parameter item                  | Specification         |
|---------------------------------|-----------------------|
| Frame Rate of CMOS              | 0~300FPS              |
| CMOS resolution                 | 2048*2048px~256*256px |
| Exposure time setting range     | 1~100 ms              |
| Fluorescence band               | 500-535nm             |
| LED power adjustable percentage | 0~100%                |
| 410nm LED Output Power Range    | 0~100 $\mu$ W         |
| 470nm LED Output Power Range    | 0~100 $\mu$ W         |

## Tricolor Multichannel Fiber Photometry System



The system is equipped with laptop(include pre-installed software), low autofluorescence optical fiber, behavior camera, behavior camera bracket.




It can support up to 9 channels for simultaneous collection, which is suitable for simultaneous recording of multiple nerve nuclei or simultaneous recording of multiple animals.



Dual-camera fluorescence acquisition, signal interference-free, exposure time: 1-100ms; gain multiple: 1-100; frame rate of CMOS: 0-300FPS.




The excitation light source consists of 3 wavelengths of light (410nm,470nm,560nm), which can simultaneously record 2 ranges of excitation light (500 - 535nm, 575 - 628nm) — from the excitation of green fluorescent proteins (e.g. GCaMP) and red fluorescent proteins (e.g. RCaMP). 410nm can reflect the background noise signal.




4 Input interfaces are configured to support 4 kinds of external signal Input automatic markers; 4 Output interfaces, supporting output TTL signal to trigger external third-party equipment, to meet the closed-loop research.




The system is equipped with fluorescence data acquisition and analysis software, and the software integrates the behavior record analysis module, which can realize the animal tracking automatic marking.



- The software presets 6 kinds of fluorescence excitation output modes, which can be used in different experimental application scenarios.
- The excitation light power can be directly adjusted on the operating software, the adjustment range is 0~100%, the adjustment precision is 0.1μW.



- The software has three marking functions, external marker, manual marker, automatic output marker. And it supports setting 10 markers at the same time.
- The acquisition stage can output TTL signal to trigger external devices when the fluorescence signal is greater than the set  $\Delta F/F$  threshold or when the animal enters a certain ROI area.



Equipped with special software—contains signal acquisition records, data analysis, behavioral video acquisition and analysis functions. The software includes the following analysis functions:

- ①The data can be smoothed, baseline correction, motion correction, view and output the processing process;
- ②Flexible selection and addition of event markers;
- ③One-click analysis to generate heat maps and mean  $\pm$  sem;
- ④Effective peak statistical analysis, AUC analysis;
- ⑤Behavioral data can be analyzed trajectories and generate heat maps;
- ⑥Analysis results can be exported to CSV, Heat-map, df/f, Z-Score and other formats.

## Micropipette Puller ▶

### MP-500 Micropipette Puller

RWD MP-500 micropipette puller is developed to solve the microelectrode problems in patch clamp and microinjection experiments. It has avant-garde intelligent operation interface, unique fixed position design of the heating filament and other user-friendly design. MP-500 integrates stability, intelligence and humanization, which can meet the requirements of microelectrode-related experiments and enhance user experience at the same time. In summary, MP-500 is an excellent choice for microelectrode related experiments.



#### Four Highlights:

Unique fixed position design of the heating filament,  
The integrated plug-in humidity control chamber,  
Full-color capacitor bilingual operation interface,  
The pulling of micropipettes is stable and efficient

#### Features:

- The unique fixed position design of the heating filament makes it easier to replace various models of filaments.
- The user interface adopts capacitive color touch screen (1024\*600px) with high sensitivity.
- The integrated plug-in humidity control chamber is easy to disassemble, which can keep the air around the filament and the capillary glass dry during the pulling process to avoid the adverse effect of humid air on the pulling result.
- Two microelectrodes with symmetrical tips are stably produced every time. The tip diameter can be less than 0.1  $\mu\text{m}$ , the taper length 3~15 mm, and the impedance range from 1 to above 100 $\Omega$  with high repeatability.
- According to the type of the capillary glass tube and the heating filament and the kinds of glass microelectrodes to be pulled (patch clamp microelectrodes, intracellular recording microelectrodes, etc.), the corresponding pull protocols are provided as a reference and can be modified according to special requirements.
- The device can store up to 99 pulling programs.
- The system diagnostic function can be used to manually check whether all the components of the puller are in normal working condition.
- The preheat mode can reduce the influence of the accumulation of jaw heat on the consistency of the pulling results during continuous pulling.
- The safe heat mode can avoid the damage of the filament from high temperature.
- The front panel cover adopts a roller design, which is convenient for users to open easily.

#### Order Information:

| Model of micropipette puller              | Product Description  |
|---|--|
| MP-500 (HF-3025B)                         | Users can choose one of the three kinds of commonly used heating filament according to actual needs. One has already been installed on the puller before leaving the factory, and the other one is for replacement. In addition, we will add an additional filament that can be used for different types of experiments. |
| MP-500 (HF-3030B)                         |  |
| MP-500 (HF-2545B)                         |  |
| Model of heating filament                 | Platinum-iridium alloy(platinum90%, iridium10%)W× L of square box; HF=Heating Filament, B=Box  |
| HF-3025B                                  | 2.5×3.0mm  |
| HF-3030B                                  | 3.0×3.0mm  |
| HF-2545B                                  | 4.5×2.5mm  |
| Model of borosilicate glass with filament | Borosilicate glass with filament(O.D.×I.D.×L)B=Borosilicate, F=Filament  |
| B-15086-10F                               | 1.5×0.86×100mm   |
| B-10050-10F                               | 1.0×0.5×100mm  |
| B-10078-10F                               | 1.0×0.78×100mm   |
| B-12069-10F                               | 1.2×0.69×100mm   |
| B-12094-10F                               | 1.2×0.94×100mm   |
| B-150110-10F                              | 1.5×1.1×100mm  |

## Motorized Micromanipulator ▶

### MM-500 Motorized Micromanipulator

Electric micromanipulators are mostly used to manipulate recording electrodes and stimulate the movement of electrodes in electrophysiological experiments.

The MM-500 was designed to meet a wide variety needs for the scientific community, and is suitable for patch clamp experiments, microinjection, extracellular recording, intra-cellular recording and precision robotic positioning applications.



#### Features:

- Small and light size, humanized bilingual operation interface, flexible program debugging.
- The TFT screen shows the position and moving speed of the activated manipulator as well as X, Y and Z.
- High-resolution stepper motor meets four-axis drive, the maximum distance of X, Y, Z and diagonal movement is 25mm, with good stability and small drift.
- A rotary photoelectric encoder can control two manipulators and switch freely. Origin/working position function can meet the requirement of rapid return of electrode to original position/working position.
- Compatible with mainstream micromanipulation platform, can be used with any microscope, can hold a variety of probes and microsyringes.

#### Technical parameters:

| Parameter item                | Specification   |
|-------------------------------|---|
| Travel                        | 25 mm on all four axes                                    |
| Maximum Speed                 | 2.9mm/s   |
| Resolution                    | 0.132 $\mu$ m/0.033 $\mu$ m                               |
| Long Term Stability           | 10nm/h at 24°C.   |
| Repeated positioning accuracy | Single axis $\leq$ 1 $\mu$ m, three axis $\leq$ 2 $\mu$ m |

#### Order Information:

| Model of micropipette puller | Product Description   |
|------------------------------|---|
| MM-500-2                     | Motorized micromanipulators(Double-handed setup)  |
| MM-500-L                     | Motorized micromanipulators(Left-handed setup)  |
| MM-500-R                     | Motorized micromanipulators(Right-handed setup)   |
| MM-500L                      | Left-handed mechanical manipulator  |
| MM-500R                      | Right-handed mechanical manipulator   |
| MT-500                       | Microscope moving platform: Contain 3 supporting columns, 1XY mobile platform, 1microscope mounting bracket, 1XY mobile platform mounting bracket |

## In Vivo Imaging Solution

### Laser Speckle Imaging System ▶

Based on the Laser Speckle Contrast Imaging(LSCI) technology, RFLSI III provides real-time and full-field blood perfusion calculation. In non-contact way, reliable quantified data and high definition images could be obtained for various of pre-clinical research. It can be used to record the microcirculation of any exposed organ.



12X Zoom with Autofocus, 120 frames/sec

Non-contact, No contrast agent

#### Product Features

Real time, Full field imaging

Multi-output: Data, Video, Image; on-line state and off-line state

#### Performance Characteristics:

- The 2048×2048 resolution of the images that directly generated is the highest resolution among all the similar products.
- 12X Zoom with autofocus and 120 frames per second ensure the details of the blood perfusion calculated are as much as the eye can see.
- Stable laser output combine with optimized algorithm guarantee the accuracy of the quantification of blood perfusion.
- Changes of the vessel diameter contained in the recorded data along with videos and images could be analysed in on-line state or off-line state to maximum efficiency.
- Solutions of the scientific research with timely after-sale service offered by RWD will make a great assistant in your daily research.

#### References:

- Activation of the XBP1s/O-GlcNAcylation pathway improves functional outcome after cardiac arrest and resuscitation in young and aged mice
- Snake C-Type Lectins Potentially Contribute to the Prey Immobilization in *Protobothrops mucrosquamatus* and *Trimeresurus stejnegeri* Venoms
- Protein O-GlcNAcylation alleviates small intestinal injury induced by ischemia-reperfusion and oxygen-glucose deprivation
- Mechanically Induced Vasospasm Evaluation of Spasmolytic Efficacy of 10 Pharmaceutical Agents Using Laser Speckle Contrast Imaging
- Aging exacerbates impairments of cerebral blood flow autoregulation and cognition in diabetic rats
- Microvascular Injury in Mild Traumatic Brain Injury Accelerates Alzheimer-like Pathogenesis in Mice
- SP6616 as a Kv2.1 inhibitor efficiently ameliorates peripheral neuropathy in diabetic mice
- Arterial Supercharging Is More Beneficial to Flap Survival Due to Quadruple Dilation of Venules


► Related research topics

- Cerebral Cortex Blood Perfusion Monitoring
- MCA-Occlusion Assessment
- Cortical Spreading Depression Observation
- Hindlimb Ischaemia Research
- Skin Burn/Skin Flap Transplantation Research
- Organ Microcirculation Observation

Technical Parameters:


|                            |                                  |
|----------------------------|----------------------------------|
| Type                       | RFLSI III                        |
| Scope of perfusion         | 0 – 5000PU                       |
| Imaging acquisition rate   | up to 120 frames/sec             |
| Working distance           | 10-35cm                          |
| Monitoring Area            | 1.57 × 1.57mm – 90 × 90mm        |
| Image resolution up to     | 2048 × 2048                      |
| Zoom range                 | 12X Optical Zoom with Auto-focus |
| Maxmium spatial resolution | up to 3µm/pixel                  |
| Wavelength                 | 785nm                            |
| Indicated Laser Wavelength | 660nm                            |

Software interface:



Middle Cerebral Artery Occlusion Mouse Model with Realtime Data of Blood Perfusion and Pipe Diameter.

► Case show



In Vivo Imaging Solution

# Animal Administration and Blood Sampling Solution

## Drug Delivery Cannula

Cannula implantation is the international general method for acute and chronic intracranial drug administration in animal experiments. RWD offers customized services for users, the users can customize according to specific experimental requirements.



### Single Cannula

#### Components:

#### Guide cannula



Guide cannula consists of a threaded cylindrical plastic pedestal and a stainless steel tube. It has to be preimplanted into brain and fixed with dental cement. Both dummy cannula and internal cannula can be inserted into guide cannula individually.

#### Dummy cannula




Stainless steel wire and threaded plastic cap make a dummy cannula. It can be inserted into guide cannula and fixed by screw to seal the guide cannula and prevent outflow of tissue fluid after injection. The length of the stainless steel wire can be customized.

#### Metal Cannula



The metal materials avoid biting after the cannula implantation.

#### Internal cannula




Internal cannula can be fixed into guide cannula through fixing screw and its tail end is connected with PE tube for injections. The length can be customized.

#### Fixing screw



Fixing screw fixes internal cannula to guide cannula.


#### Cording Diagram:



Single-Cannula connection

#### Diagrams:

- A - The length of the stainless steel tubing below the plastic pedestal.
- B - The height of the cylindrical plastic base.
- C - The length of the stainless steel tubing below the plastic pedestal, can be customized.
- D - External diameter of stainless steel guide cannulas.
- E - External diameter of cylindrical plastic base.
- F1 - External diameter of Fixing screw.
- F2 - External diameter of Dummy Cannulas.
- H1 - Height of Fixing screw.
- H2 - Height of Dummy Cannulas.
- G1 - The internal cannulas projection length beyond the guide cannulas tip, can be customized.
- G2 - The dummy cannulas projection length beyond the guide cannulas tip, can be customized.

**Ordering Instructions:****Diagram for Customization**

- D-outer diameter of the Guide Cannula (Unit:g) , five sizes 22, 23, 24, 26 and 27g are available.
- C-length (Unit: mm) of the metal tube below the plastic pedestal (B).
- G1-length (Unit: mm) of the protrudent Injector through the Guide Cannula, 0.5 or 1.0mm is recommended.
- G2-length (Unit: mm) of the protrudent Dummy Guide, 0.5 or 1.0mm is recommended.

**Order information (Unit: mm):**

| Matching | Guide Cannula<br>(OD x ID) | Internal Cannula<br>(OD x ID) | Dummy Cannula<br>(plastic OD) | Dummy Cannula<br>(metal OD) | Fixing Screw | PE Tube | Range of customized length |        |        |
|----------|----------------------------|-------------------------------|-------------------------------|-----------------------------|--------------|---------|----------------------------|--------|--------|
|          |                            |                               |                               |                             |              |         | C                          | G1     | G2     |
| Cannula  | 62001(23g)<br>0.64x0.45    | 62201<br>0.41x0.25            | 62101<br>0.4                  | 62106<br>0.4                | 62502        | 62329   | 0.5~15                     | 0.5~15 | 0.5~15 |
|          | 62002(24g)<br>0.56x0.38    | 62202<br>0.36x0.2             | 62102<br>0.3                  | 62107<br>0.3                |              |         |                            |        |        |
|          | 62003(26g)<br>0.48x0.34    | 62203<br>0.3x0.14             | 62102<br>0.3                  | 62107<br>0.3                |              |         |                            |        |        |
|          | 62004(27g)<br>0.41x0.25    | 62204<br>0.21x0.11            | 62104<br>0.2                  | 62108<br>0.2                |              |         |                            |        |        |
|          | 62005(22g)<br>0.71x0.45    | 62201<br>0.41x0.25            | 62101<br>0.4                  | 62106<br>0.4                |              |         |                            |        |        |

**Notes:** The height and external diameter of the cylindrical plastic base are 8mm and 3.5mm respectively. Plastic or metal dummy cannula is optional.

## Double Cannula

**Introduction:**

RWD has launched I, II, III and IV series double cannula systems for acute and chronic intracranial drug administration in animal experiments.

Series I: Recommended for rat or mouse use, 0.48mm (outer diameter) x 0.34mm (internal diameter).

Series II: Recommended for rat use, 0.64mm (outer diameter) x 0.45mm (internal diameter).


Series III: Specifically for mouse use, 0.48mm (outer diameter) x 0.34mm (internal diameter). The height of cylindrical plastic pedestal of Guide cannula is different from Series I.

Series IV: 0.41mm (outer diameter) x 0.25mm (internal diameter).

**Application:**

Double cannula system is used for bilateral implantation, suitable for injecting one or two drugs into two different brain regions.

## Cording Diagram:



## Comments:

- B - Height of cylindrical plastic pedestal which is 7.8mm for Series I, II & IV and 5.5mm for Series III
  - ★C - Length of stainless steel tube below the cylindrical plastic pedestal. The default value is 13.7mm for Series I & III, 20.0mm for Series II and 10.0mm for Series IV. The length can be customized with default value as maximum
  - ★C.C - Center distance between double guide cannulas, double dummy cannulas and double injectors. These 3 values should be the same in one system. It can be customized (range: 0.5mm-3.0mm)
- D - Outer diameter of stainless steel tube and it has 3 sizes: 0.48mm for Series I & III, 0.64mm for Series II and 0.41mm for Series IV
- E - Outer diameter of the cylindrical plastic pedestal and it comes with 3.5mm or 5.0mm. It is related to the outer diameter and centre distance between double guide cannulas
- F1 - Outer diameter of fixing screw and it comes with 5.0mm (when E=3.5mm) or 7.0mm (when E=5.00mm)
- F2 - Outer diameter of guide cannula cap and it comes with 5.5mm (when E=3.5mm) or 7.0mm (when E=5.0mm)
- H1 - Height of the fixing screw and it comes with 7.5mm (when E=3.5mm) and 8.0mm (when E=5.0mm)
- H2 - Height of guide cannula cap and it comes with 7.00mm (when E=3.5mm) and 7.5mm (when E=5.0mm)
  - ★G1 - Length of protruding part after internal cannula is inserted into guide cannula. This length can be customized
  - ★G2 - Length of protruding part after dummy cannula is inserted into guide cannula. This length can be customized

## Ordering Instructions:

Determine the target brain region for drug administration and then determine target values of C, C.C, G1 and G2. These values can not exceed the default ones. Please refer to :

- ★ part of comments of cording diagram for more detail.
- ★ C: Maximal customized length is 13.7mm for Series I & III, 20.0mm for Series II and 10.0mm for Series IV.
- ★ C.C: Customized range is 0.5mm-3.0mm.
- ★ G1: Maximal customized length is 6.0mm for Series I, 4.0mm for Series II, 8.0mm for Series III and 4.0mm for Series IV.
- ★ G2: Maximal customized length is 2.0mm for Series I, 5.0mm for Series II, 4.0mm for Series III and 5.0mm for Series IV.

## Order Information: (Unit: mm)

## Series I: Recommended for rat or mice

| Matching  | C.C | Guide Cannula | Internal Cannula | Dummy Cannula | Cap   | Fixing Screw | PE Tube | Range of customized length |         |         |
|---|-----|---------------|------------------|---------------|-------|--------------|---------|----------------------------|---------|---------|
|   |     |               |                  |               |       |              |         | C                          | G1      | G2      |
| Cannula (26G)<br>O.D.0.48mmxI.D.0.34 mm<br>Dummy Cannula:<br>O.D.0.30mm<br>Internal Cannula:<br>O.D.0.30mm x I.D.0.14mm<br>Plastic Pedestal of Guide Cannula:<br>Height B = 7.8mm<br>Diameter E = 3.5mm | 0.5 | 62021         | 62250            | 62121         | 62523 | 62521        | 62320   | 0.5~13.7                   | 0,0.5~6 | 0,0.5~2 |
|   | 0.6 | 62022         | 62250            | 62122         |       |              |         |                            |         |         |
|   | 0.8 | 62023         | 62223            | 62123         |       |              |         |                            |         |         |
|   | 1.0 | 62024         | 62224            | 62124         |       |              |         |                            |         |         |
|   | 1.2 | 62025         | 62225            | 62125         |       |              |         |                            |         |         |
|   | 1.3 | 62026         | 62226            | 62126         |       |              |         |                            |         |         |
|   | 1.4 | 62027         | 62227            | 62127         |       |              |         |                            |         |         |
|   | 1.5 | 62028         | 62228            | 62128         |       |              |         |                            |         |         |
|   | 1.6 | 62029         | 62229            | 62129         |       |              |         |                            |         |         |
|   | 1.8 | 62030         | 62230            | 62130         |       |              |         |                            |         |         |
| Plastic Pedestal of Guide Cannula:<br>Diameter E = 5 mm   | 2.0 | 62031         | 62231            | 62131         | 62524 | 62522        | 62329   |                            |         |         |
|   | 2.2 | 62032         | 62232            | 62132         |       |              |         |                            |         |         |
|   | 2.4 | 62033         | 62233            | 62133         |       |              |         |                            |         |         |
|   | 2.6 | 62034         | 62234            | 62134         |       |              |         |                            |         |         |
|   | 2.7 | 62035         | 62235            | 62135         |       |              |         |                            |         |         |
|   | 3.0 | 62036         | 62236            | 62136         |       |              |         |                            |         |         |

## Series II: Recommended for rat

| Matching   | C.C | Guide Cannula | Internal Cannula | Dummy Cannula | Cap   | Fixing Screw | PE Tube | Range of customized length |         |         |
|--|-----|---------------|------------------|---------------|-------|--------------|---------|----------------------------|---------|---------|
|  |     |               |                  |               |       |              |         | C                          | G1      | G2      |
| Cannula(23G)<br>O.D.0.64mmx I.D.0.45 mm<br>Dummy Cannula:<br>O.D.0.40mm<br>Internal Cannula:<br>O.D.0.41mm x I.D.0.25mm<br>Plastic Pedestal of Guide Cannula:<br>Height B = 7.8mm<br>Diameter E = 3.5 mm | 0.8 | 62037         | 62237            | 62137         | 62523 | 62521        | 62323   | 0.5~20                     | 0,0.5~4 | 0,0.5~5 |
|  | 1.0 | 62038         | 62238            | 62138         |       |              |         |                            |         |         |
|  | 1.2 | 62039         | 62239            | 62139         |       |              |         |                            |         |         |
|  | 1.3 | 62040         | 62240            | 62140         |       |              |         |                            |         |         |
|  | 1.4 | 62041         | 62241            | 62141         |       |              |         |                            |         |         |
|  | 1.5 | 62042         | 62242            | 62142         |       |              |         |                            |         |         |
|  | 1.8 | 62043         | 62243            | 62143         |       |              |         |                            |         |         |
| Plastic Pedestal of Guide Cannula :<br>Diameter E = 5 mm   | 2.0 | 62044         | 62244            | 62144         | 62524 | 62522        | 62329   |                            |         |         |
|  | 2.2 | 62045         | 62245            | 62145         |       |              |         |                            |         |         |
|  | 2.4 | 62046         | 62246            | 62146         |       |              |         |                            |         |         |
|  | 2.5 | 62047         | 62247            | 62147         |       |              |         |                            |         |         |
|  | 2.8 | 62048         | 62248            | 62148         |       |              |         |                            |         |         |
|  | 3.0 | 62049         | 62249            | 62149         |       |              |         |                            |         |         |

## Series III: Recommended for mice

| Matching   | C.C   | Guide Cannula | Internal Cannula | Dummy Cannula | Cap   | Fixing Screw | PE Tube | Range of customized length |         |         |
|--|-------|---------------|------------------|---------------|-------|--------------|---------|----------------------------|---------|---------|
|  |       |               |                  |               |       |              |         | C                          | G1      | G2      |
| Cannula (26G)<br>O.D.0.48mmxI.D.0.34 mm<br>Dummy Cannula:<br>O.D.0.30mm<br>Internal Cannula :<br>O.D.0.30mm x I.D.0.14mm<br>Plastic Pedestal of Guide Cannula:<br>Height B = 5.5mm<br>Diameter E = 3.5mm | 0.5   | 62050         | 62250            | 62121         | 62523 | 62521        | 62320   | 0.5~13.7                   | 0,0.5~8 | 0,0.5~4 |
|  | 0.8   | 62052         | 62223            | 62123         |       |              |         |                            |         |         |
|  | 1.0   | 62053         | 62224            | 62124         |       |              |         |                            |         |         |
|  | 1.2   | 62054         | 62225            | 62125         |       |              |         |                            |         |         |
|  | 1.3   | 62055         | 62226            | 62126         |       |              |         |                            |         |         |
|  | 1.4   | 62056         | 62227            | 62127         |       |              |         |                            |         |         |
|  | 1.5   | 62057         | 62228            | 62128         |       |              |         |                            |         |         |
|  | 1.6   | 62058         | 62229            | 62129         |       |              |         |                            |         |         |
|  | 1.8   | 62059         | 62230            | 62130         |       |              |         |                            |         |         |
| 2.0  | 62060 | 62231         | 62131            |               |       |              |         |                            |         |         |

**Series IV: Recommended for mice**

| Matching  | C.C   | Guide Cannula | Internal Cannula | Dummy Cannula | Cap   | Fixing Screw | PE Tube | Range of customized length |         |         |
|---|-------|---------------|------------------|---------------|-------|--------------|---------|----------------------------|---------|---------|
|   |       |               |                  |               |       |              |         | C                          | G1      | G2      |
| Cannula (27G)<br>O.D.0.41mm x I.D.0.25 mm<br>Dummy Cannula:<br>O.D.0.20mm<br>Internal Cannula:<br>O.D.0.21mm x I.D.0.11mm<br>Plastic Pedestal of Guide Cannula:<br>Height B = 7.8mm<br>Diameter E = 3.5mm | 0.5   | 62061         | 62251            | 62161         | 62523 | 62521        | 62320   | 0.5~10                     | 0,0.5~4 | 0,0.5~5 |
|   | 0.6   | 62062         | 62251            | 62162         |       |              |         |                            |         |         |
|   | 0.8   | 62063         | 62263            | 62163         |       |              |         |                            |         |         |
|   | 1.0   | 62064         | 62264            | 62164         |       |              |         |                            |         |         |
|   | 1.2   | 62065         | 62265            | 62165         |       |              |         |                            |         |         |
|   | 1.3   | 62066         | 62266            | 62166         |       |              |         |                            |         |         |
|   | 1.4   | 62067         | 62267            | 62167         |       |              |         |                            |         |         |
|   | 1.5   | 62068         | 62268            | 62168         |       |              |         |                            |         |         |
|   | 1.6   | 62069         | 62269            | 62169         |       |              |         |                            |         |         |
|   | 1.8   | 62070         | 62270            | 62170         |       |              |         |                            |         |         |
| 2.0   | 62071 | 62271         | 62171            |               |       |              |         |                            |         |         |

Notes : (1)For cap, plastic or metal is optional. Metal cap 62525 (corresponding plastic cap 62523), Metal cap 62526 (corresponding plastic cap 62524).

(2)The difference between Series I and Series III is the height of the pedestal, the former is 7.8mm and the later is 5.5mm, and Series III is specially designed for mice

## PE & Silicone Tubings

PE(Polyethylene) / HDPE(High Density Polyethylene) tube, with good fluidity, high thermal stability and good flexibility, is widely used in the field of life science. Currently, RWD has 9 types of PE tube and 1 type of silicone tube, including the international standard tubes PE10, PE20 and PE50, which can be used in microdialysis, animal artery intubation, venous cannula, bile duct intubation, subcutaneous implantation of buried pipe, lateral ventricle buried pipe, portal vein intubation and so on. These tubes can be used together with syringe, guide cannula, microdialysis connector and luer stubs.

### Order Information:

| No. | Model | Product Description                      |
|-----|-------|--|
| 1   | 62320 | PE Tubing-0.85*0.42mm, pkg of 10m        |
| 2   | 62321 | PE Tubing (HDPE)-0.85*0.42mm, pkg of 10m |
| 3   | 62322 | PE Tubing-0.50*0.25mm, pkg of 10m        |
| 4   | 62323 | PE Tubing-1.00*0.50mm, pkg of 10m        |
| 5   | 62324 | PE10 Tubing-0.61*0.28mm, pkg of 10m      |
| 6   | 62325 | PE20 Tubing-1.09*0.38mm, pkg of 10m      |
| 7   | 62326 | PE25 Tubing-0.91*0.46mm, pkg of 10m      |
| 8   | 62327 | PE50 Tubing-0.97*0.58mm, pkg of 10m      |
| 9   | 62328 | PE60 Tubing-1.20*0.80mm, pkg of 10m      |
| 10  | 62329 | PE Tubing-1.50*0.50mm, pkg of 10m        |

Note: Except 62322, other PE tube can withstand high-pressure steam sterilization(121°C, 30min).



### Screws



| Model | Product Description |
|-------|---------------------|
| 62513 | Screw-M1.0xL2.0mm   |
| 62514 | Screw-M1.2xL2.0mm   |
| 62515 | Screw-M1.4xL2.5mm   |

# Osmotic Pumps




The osmotic pump has become the best tool for continuous dose of unrestrained laboratory animals. It can be implanted subcutaneously or intraperitoneally in animals. Minimizes side effects and experimental variables, and reduces the risk of stress due to frequent animal handling. And it is used for continuous delivery of a wide range of drugs, including antibodies, growth factors, cytokines, chemotherapeutics, addictive drugs, hormones, steroids and more.



### Features:

- Ensure continuous and uninterrupted administration, maintaining the preset speed for the administration of the animal.
- Storage capsules have good drug compatibility, and protection for short half-life drugs.
- Capsules have good biocompatibility and do not affect the physiology of animals.
- Easy to use, It allows one-time operation in an experiment instead of frequent operation.
- Avoid the stress caused by repeated administration, reducing the impact of non-experimental factors on the experimental results.
- Small size, can be applied to experiments on mice, neonatal rats and other small animal.
- Ensure that the drug is delivered to the target site accurately.

The osmotic pump can be directly implanted in the abdominal cavity or the back of the neck for drug administration, or it can be applied to specific tissue and organ, such as blood vessels, brain, spinal cord, lungs, liver, etc., by connecting the catheter or other accessories.



- Flow regulator: The sole outlet for drug solution, which is used to connect the catheter to target administration area.
- Storage capsules: Fill with solution. It has good corrosion resistance to acid and alkali.
- Shell: Firmness, designed to protect the internal structure of osmotic pump, with fluid passability characteristics.
- Salt interlayer: absorption of tissue fluid by osmosis. After the salt interlayer absorbs the tissue fluid, it can only be extruded inwards to slowly pump the solution out of the capsule.
- Pellicle: prevents the salt layer from diffusing into the tissue, and ensures that the tissue fluid can enter the salt layer to squeeze the storage capsule and pump out the drug solution

### Technical Parameters:

| Capsule<br>(10pcs/box)        | Drug capacity | 100 µL                                |       |       |       | 200 µL      |       |       |       | 2 mL  |       |       |
|-------------------------------|---------------|---------------------------------------|-------|-------|-------|-------------|-------|-------|-------|-------|-------|-------|
|                               | Type          | 1003D                                 | 1001W | 1002W | 1004W | 2003D       | 2001W | 2002W | 2004W | 2ml1W | 2ml2W | 2ml4W |
| Duration                      | 3d            | 1w                                    | 2w    | 4w    | 3d    | 1w          | 2w    | 4w    | 1w    | 2w    | 4w    |       |
| Pumping Rate(µL/hr)           | 1.0           | 0.5                                   | 0.25  | 0.125 | 2.0   | 1.0         | 0.5   | 0.25  | 10.0  | 5.0   | 2.5   |       |
| Flow regulator<br>(10pcs/box) | Cap           | Without                               |       |       |       | PE material |       |       |       |       |       |       |
|                               | Draft tube    | 304 stainless steel material          |       |       |       |             |       |       |       |       |       |       |
| Infusion tube<br>(1pcs/box)   | Syringe       | Flat head for the injection of liquid |       |       |       |             |       |       |       |       |       |       |

### Osmotic pump kit ordering information:

|                       | Model | Description   |
|-----------------------|-------|---|
| Brain Infusion Kit    | Bic-3 | For 1-3mm deep brain area administration (10 Bags/ Package) |
|                       | Bic-5 | For 1-3mm deep brain area administration (10 Bags/ Package) |
| Catheter Delivery Kit | Mic   | For long-distance muscle administration (10 Bags/ Package)  |



Different accessory kits

Animal Administration and Blood Sampling Solution

## Blood Sampling and Administration System

The catheter is surgically inserted into the target location (jugular vein, carotid artery, bile duct, portal vein, etc.), and the Implant Button is then fixed subcutaneously, allowing simultaneous drug administration and blood sampling. It is suitable for manual continuous and repeated blood sampling and drug administration, long-term automatic drug administration and blood sampling by machine.

### Features:

- Suitable for up to 10 sampling sites such as jugular vein, carotid artery, femoral vein and bile duct.
- The design is compact and firm, not easy to be damaged by animals, and the biocompatibility has been verified by experiments.
- It can be matched with syringes, catheters and other sampling methods.
- Effectively reduce the stress reaction caused by repeated injection.



### Components:

#### Implant Button



The Implant Button consists of transparent shell, magnet, felt and capillary tube. Catheter is connected to capillary tube and inserted into a blood vessel. The Button is implanted subcutaneously and fixed by suturing the epidermis.

#### Pharmport



The Pharmport consists of shell (red and white) and capillary tube, sizes: 25g and 22g. The Pharmport is connected with the Implant Button as a simple infusion and blood sampling device.

#### Catheter connector



The Catheter connector is composed of shell and capillary tube, size: 25G, 22G. For connecting Pharmport and PE tubes in vitro.

#### Syringe connector



Syringe connector is composed of transparent shell and capillary tube, with a size of 25G. The front end of the injection needle can be directly inserted into the Pharmport for infusion and blood sampling, the rear end can be connected with the standard syringe.

#### Implant Button Cap



The Implant Button cap is made of high quality duralumin alloy and magnet. When the Implant Button is not in use, The cap can be used to ensure that the injection end is not contaminated. The magnet allows the operator to quickly attach the cap to Implant Button.

#### Handling



The Handling is made of high quality duralumin alloy and magnet. Fix the animal with one end of the Handling to pull out the Implant Button Cap, and then use the other end (including magnet) to continue to fix the animal for infusion or blood sampling operation.

#### Catheter



The size of the catheter matches the capillary of the Implant Button, and the catheter can be inserted into the blood vessels of experimental animals. The node ensures that the catheter is not easy to fall off after implantation.

**Order Information:**

| Product Description                                | Model    |
|--|----------|
| Implant Button (single channel 25G for mice-white) | BMSW25   |
| Implant Button (mouse dual channel 25G/25G)        | BMD25/25 |
| Implant Button (single channel 22G for rat-white)  | BRSW22   |
| Implant Button (single channel 22G for rat-Rad)    | BRSR22   |
| Implant Button (Rat dual channel 22G/22G)          | BRD22/22 |
| Catheter connector (22G)                           | C22      |
| Catheter connector (25G)                           | C25      |
| Syringe connector                                  | SC-1     |
| Pharmport (22G Red)                                | P22-R    |
| Pharmport (22G White)                              | P22-W    |
| Pharmport (25G Red)                                | P25-R    |
| Pharmport (25G White)                              | P25-W    |

| Product Description                              | Model  |
|--|--------|
| Implant Button Cap (mouse single channel-silver) | CMSC-S |
| Implant Button Cap (mouse single channel-gold)   | CMSC-G |
| Implant Button Cap (mouse dual channel-silver)   | CMDC-S |
| Implant Button Cap (mouse dual channel-gold)     | CMDC-G |
| Implant Button Cap (rat-silver)                  | CRSC-S |
| Implant Button Cap (rat-gold)                    | CRSC-G |
| Handling   | BH     |
| Carotid Artery Catheters for Rats                | PURCA  |
| Jugular Vein Catheters for Rats                  | PURJV  |
| Bile Duct Catheters for Rats                     | PURBD  |
| Portal Vein Catheters for Rats                   | PURPV  |
| Femoral Artery Catheters for Rats                | PURFA  |
| Femoral Vein Catheters for Rats                  | PURFV  |
| Jugular Vein Catheters for Guinea pigs           | PUCJV  |
| Carotid Artery Catheters for Guinea pigs         | PUCCA  |
| Carotid Artery Catheters for Mice                | PUMCA  |
| Jugular Vein Catheters for Mice                  | PUMJV  |

| Product Description                  | Model  |
|--------------------------------------|--------|
| Swivels (Single Channel 22G)         | SS1-22 |
| Swivels (Single Channel 25G)         | SS1-25 |
| Swivels Mount for Mice               | MLA-M  |
| Swivels Mount for Rats               | MLA-R  |
| Tether (mouse single channel)        | TM-1   |
| Tether (mouse dual channel)          | TM-2   |
| Tether (drug administration for Rat) | TR-1dr |
| Tether (blood sampling for Rat)      | TR-1bl |
| Tether (rat dual channel)            | TR-2   |
| Bile Loop                            | DZXHQ  |
| PE Tubing 22                         | BC-22  |
| PE Tubing 25                         | BC-25  |

## Syringe Pumps

### 462 Syringe Pump

RWD series syringe micro-injection pumps are widely used in scientific research. The upgrade from the first generation to the current second-generation products includes single-channel and dual-channel options.



#### Features:

- Control the stepper motor for stable and high-precision operation at low speed.
- The communication with the touch screen of the host computer can control the output torque and rotation speed of the stepper motor according to the input, and then control the injection thrust and injection speed.
- Reserved communication interface (BNC, USB-B, and PC communication, RS485 multi-injection pump communication).
- Can input and output, 15-pin D-sub Connector.
- With encoder, feedback speed, alarming when speed is abnormal, equipment operation suspend.
- Safety certification: EMC, ESD, and CE.
- Running indicator: convenient for users to intuitively discover the current status of the pump.

#### Order Information:

| Model | Name                 | Note         |
|-------|----------------------|--------------|
| R462  | RWD 462 Syringe Pump | Dual-channel |

#### Technical Parameters:

|                         |  |
|-------------------------|--|
| Operating mode          | Injection/withdrawal                       |
| Injection accuracy      | $\pm < 0.35\%$                             |
| Repeatability           | $\pm < 0.05\%$                             |
| Syringe range           | 0.5 $\mu$ L-60 mL                          |
| Minimum flow rate       | 0.0001 $\mu$ L / min (0.5 $\mu$ L syringe) |
| Maximum flow rate       | 102 mL/min (60 mL syringe)                 |
| Linear thrust           | 50lbs (23 kg)                              |
| Operating temperature   | 4 °C to 40 °C                              |
| Storage temperature     | -10 °C ~ 70 °C                             |
| Humidity                | 20% to 80%, non-condensing                 |
| Communication interface | USB-B, BNC, RS485 and TTL ports            |
| Advance stroke per step | 0.046 $\mu$ m / step                       |
| Display                 | 5.0-inch LCD touch screen                  |

## R480 Nanoliter Injection Pump

The newly launched R480 Nanoliter Injection Pump is highly precise & airtight, stable & easy to operate, which can be used to accurate microinjection or suction for oocytes, larva, protozoa, animal brain, eyeball, etc.



### Features:

- Stable operation and high injection precision
- High airtightness, ensures no leakage during injection or withdrawal
- The brightness of the screen can be adjusted to meet different experimental situations
- Added the function to set needle retention time. A timer is no longer required if the needle is retained after the injection

### Technical Parameters:

|                                     |                           |
|-------------------------------------|---------------------------|
| Operating mode                      | Injection/withdrawal      |
| Injection Volume                    | 0.6 nL-5000 nL            |
| Minimum flow rate                   | 1 nL/s                    |
| Maximum flow rate                   | 200 nL/s                  |
| Emptying rate                       | 10nL-200nL/s              |
| Withdrawing rate                    | 10 nL-200 nL/s            |
| Glass Capillary size                | OD1.14mm, ID0.53mm        |
| Minimum injection volume resolution | 0.1nL                     |
| Circle                              | 1-8000                    |
| Display                             | 5.0-inch LCD touch screen |

### Order Information:

| Accessory Name                                  | Description   | Quantity |                         |
|---|---|----------|-------------------------|
| R480 Nanoliter Injection Pump                   | Host(including Injector Head )                                | 1        | Standard configurations |
| Power wire                                      | Plug-in power wire for device power input                     | 1        |                         |
| 3.5-inch (8.89cm) glass capillary               | Attached to the Injector head for injection;OD1.14mm,ID0.53mm | 100      |                         |
| 7-inch (17.78cm) glass capillary                | Attached to the Injector head for injection;OD1.14mm,ID0.53mm | 100      |                         |
| Fixed fittings                                  | For conjunction with brain stereotaxer and micromanipulator   | 1        |                         |
| Filling needle                                  | For Filling Mineral Oil: 30 Gauge, 2" length                  | 1        |                         |
| Sealing material                                | For replacing the sealing material inside the Injector Head   | 1        |                         |
| Foot switch                                     | For foot controlled injection                                 | 1        | Optional Accessories    |
| Mineral Oil                                     | 100mL/ bottle   | 1        |                         |
| Special fixing accessories for micromanipulator | With special accessories for RWD Motorized Micromanipulator   | 1        |                         |

# Pathological Section Solution

## Minux® Rotary Microtomes

Rotary microtome is a device for paraffin section, widely used in medical and biological research. It is used for pathological diagnosis in clinical, and used for research at the tissue and cell level.



### ► Features:

#### Visual pointer

- » Designed with the visual pointer on precise positioning system, easy to adjust the sample to the target angle, reducing trimming and sample waste.
- » Adjust multiple devices to the same angle easily and quickly.

#### The side knobs

- » Integrated with the machine, ergonomic design, make operation more comfortable and coordinated.
- » The sample feeding speed can be adjusted freely between 0-1800µm/s, more choice and more efficient.
- » Durable military-grade devices, service life is more than 20 million times.

#### Automatic section

- » It is stable and maintains slice consistency.
- » No extra skills required, easy to complete the operation.
- » Some specimens such as embryos and brains need to be cut for hundreds of pieces. It is hard for you, but easy for automatic section.

#### Easily removable waste tray

- » With magnetic function, easy to disassemble and install.
- » Special coating materials tray, easy to clean up, shorten the cleaning time to a few seconds.

## ► Technical Parameters

| Parameters   | Model | S700 Semi-Automated | S700A Fully-Automated |
|--|-------|---------------------|-----------------------|
| Section thickness setting range: 0.5 - 100 µm        |       | ✓                   | ✓                     |
| Trimming section thickness setting range: 1 - 800 µm |       | ✓                   | ✓                     |
| Electric coarse feed speed: 0-1800µm/s               |       | ✓                   | ✓                     |
| Specimen retraction: 0-250µm                         |       | ✓                   | ✓                     |
| Specimen feed: 28mm, vertical stroke: 70mm           |       | ✓                   | ✓                     |
| Specimen orientation: X/Y axis 8°, Z axis 360°       |       | ✓                   | ✓                     |
| Maximum specimen size: 55*50*30mm                    |       | ✓                   | ✓                     |
| Automated mode and manual mode switch                |       |                     | ✓                     |
| Emergency stop button                                |       |                     | ✓                     |
| Cut window function                                  |       |                     | ✓                     |
| Slice speed selection                                |       |                     | ✓                     |
| Four slice modes: Int.,Single,Cont., Multi.          |       |                     | ✓                     |
| Pedals(optional)                                     |       |                     | ✓                     |

## ► Order Information

| Cat No.  | Product Description  |
|----------|--|
| S700A    | Fully-Auto rotary microtome, include a box of blades(50 pcs) |
| S700     | Semi-Auto rotary microtome, include a box of blades(50 pcs)  |
| S700-F3  | Universal clamp  |
| S700-F4  | Standard clamp   |
| S700-F7  | Blade holder assembly  |
| S700-F9  | Blade holder upper bracket                                   |
| S700-F10 | Easy-clean waste tray  |
| S700-F6  | Antistatic waste tray  |

## ► Disposable Microtome Blades-low profile



### ► Features:

- Sharp and durable, cost-effective
- For rotary microtomes and cryostats use
- Packed in a push-out distribution box for easy access

## Order Information

| Cat No.  | Product Description   |
|----------|---|
| PB-1 Pro | Disposable Microtome Blades, 80mm*8mm*0.25mm, 50 pcs, for paraffin and frozen section |

## Minux® Cryostats

Cryostat is a device for rapidly cooling tissue to a certain hardness under low temperature conditions and then slicing. It is a quick and easy method.

Hospitals mainly use it for rapid pathological diagnosis during surgery. Researchers use it for preparation of biological tissue samples for in situ hybridization, immunohistochemistry, etc.



### ► Features

#### Stable

- **Stable temperature:** advanced temperature control system, small temperature fluctuations. Even if the window is opened for more than 1 hour, the temperature is still stable.
- **Stable slice quality:** high-precision cross rail and screw, five-phase high-resolution sampling motor and horizontal sampling slider made of aviation aluminum, form a precise sampling system to ensure stable slice quality.

#### Safe

- **UVC disinfection:** the latest low-temperature UVC disinfection technology. The irradiation intensity can effectively kill bacteria and viruses, disinfect the inner surface and the air of the cabinet, and reduce the risk of contamination by infectious materials.
- **Waste liquid management system:** a waste liquid volume display on the touch screen. There are reminders of full bottle and not in place to prevent waste liquid from flowing out and contaminating the ground.

#### Efficient

- **Efficient operation:** the color touch screen, with a joystick that integrates coarse sample feeding, trimming and section switching, and slice parameter adjustment, constitutes an efficient operating system.

- **Quick positioning:** precise positioning system with  $\pm 8^\circ$  visual pointer, easy to adjust the specimen clamp to the target angle, reducing trimming and sample waste.
- **Automatic sleep and wake-up:** saving energy, simplify manual operation and reduce waiting time.
- **Sample installation and marking:** Buckle design, convenient installation and firm fixation. The unique label design of specimen discs makes it easy to record and identify each sample to prevent confusion.
- **Efficient and thorough defrosting:** automatic defrosting for cryochamber, specimen head, quick-freeze shelf, more thorough. For manual defrosting, there are different parts to choose to improve efficiency.

## ► Technical Parameters

| Parameters  | Model | FS800 Semi-Automated | FS800A Automated |
|---|-------|----------------------|------------------|
| Section thickness setting range: 0.5 - 100 $\mu\text{m}$                          |       | ✓                    | ✓                |
| Trimming section thickness setting range: 1 - 800 $\mu\text{m}$                   |       | ✓                    | ✓                |
| Electric coarse feed speed: 0-1800 $\mu\text{m/s}$ adjustable                     |       | ✓                    | ✓                |
| Specimen retraction: 0-250 $\mu\text{m}$  |       | ✓                    | ✓                |
| Specimen feed: 28mm   |       | ✓                    | ✓                |
| Vertical stroke: 70mm   |       | ✓                    | ✓                |
| Specimen orientation with visual pointer: X/Y axis $8^\circ$ , Z axis $360^\circ$ |       | ✓                    | ✓                |
| Number of refrigerated sites: 15  |       | ✓                    | ✓                |
| Number of fast cooling sites: 2   |       | ✓                    | ✓                |
| Cryochamber temperature range: $0^\circ\text{C}\sim\sim 35^\circ\text{C}$         |       | ✓                    | ✓                |
| Automatic and manual defrosting   |       | ✓                    | ✓                |
| Automatic and manual UVC disinfection   |       | ✓                    | ✓                |
| Automated mode and manual mode switch   |       |                      | ✓                |
| Emergency stop button   |       |                      | ✓                |
| Cut window function   |       |                      | ✓                |
| Slice speed selection   |       |                      | ✓                |
| Four slice modes: Int., Single, Cont., Multi.                                     |       |                      | ✓                |
| Pedals(optional)  |       |                      | ✓                |

## ► Order Information

| Cat No. | Type     | Product Description   |
|---------|----------|---|
| FS800   | UV       | Semi-Auto,with UVC disinfection function  |
|         | Specimen | Semi-Auto,with UVC disinfection,specimen cooling system                               |
|         | Vacuum   | Semi-Auto,with UVC disinfection,Vacuum,specimen cooling system                        |
| FS800A  | Auto     | Fully-Auto,with UVC disinfection,specimen cooling system,auto section function        |
|         | Max      | Fully-Auto,with UVC disinfection,vacuum,specimen cooling system,auto section function |

# Cellular and Molecular Biology Research Solution

## Single Cell Suspension Dissociator

RWD Single cell suspension dissociator uses tissue processing tubes, tissue dissociation kits, and built-in optimization programs for preparing high-activity single-cell suspension and tissue homogenate. The preparation of single-cell suspension with high cell viability can be completed within 15-30 minutes, which makes the experiments more reproducible. 4 independent operating channels and heating jackets can automate the processes and improve the efficiency of tissue processing. It is widely used in single-cell sequencing, primary cell culture, fluorescence-activated cell sorting, cell therapy, and other research fields.



Single Cell Suspension Dissociator

### Features

- Intelligent**
  - Pre-defined standard processes and customized processes supported.
  - Processes can be imported and synchronized via USB ports.
- Gentle**
  - Stable motor speed ranges from 20 to 4000 rpm.
  - Customized tissue dissociation kits and tissue processing tubes can ensure cell viability and reduce experimental losses.
- Safe**
  - Multiple alarms guarantee the safety of both operators and samples.
  - Fully sealed tube prevents contamination of cells.
- Efficient**
  - 4 independent channels work simultaneously for different tissue processing.
  - Heating jackets are available for fully automated sample preparation.
  - Capable of single-cell suspension and tissue homogenate.

### Technical Parameters:

| Model              | DSC-400        |
|--------------------|----------------|
| Sample volume      | 20-4000mg      |
| Buffer volume      | 0.3-10mL       |
| Speed              | 20-4000 rpm    |
| Temperature        | 37°C           |
| Channels           | 4              |
| Operation method   | Touch screen   |
| Dimensions (W*D*H) | 215*380*190 mm |
| Weight             | ~9.0 kg        |

### Order Information:

| Cat No.   | Product Description  |
|-----------|--|
| DSC-400   | Single cell suspension dissociator                         |
| SCT-25    | Tissue processing tube(25pcs/box)                          |
| SCT-100   | Tissue processing tube(100pcs/box)                         |
| DHTE-5001 | High Activity Tumor Tissue Enzymatic Digestion Kit (Mouse) |

► Consumables:



**Tissue Processing Tube**

- Ethylene oxide(EO) sterilized
- Independent packaging make it convenient for storage and use

**Order Information:**

| Model         | SCT-25    | SCT-100    |
|---------------|-----------|------------|
| Sample volume | 20-4000mg |            |
| Buffer volume | 0.3-10mL  |            |
| Dimensions    | 35*89mm   |            |
| Package       | 25pcs/box | 100pcs/box |



**Enzymatic Digestion Kit**

- Targeted optimization of enzymatic hydrolysis formula can ensure high cell viability
- Preserving cell surface epitopes leads to no effect for downstream experiments

**Order Information:**

| Product Description  | Model      | Number of reactions |
|--|------------|---------------------|
| High Activity Tumor Tissue Enzymatic Digestion Kit (Mouse)           | DHTE-5001  | 50T                 |
| High Activity Neonatal Brain Enzymatic Digestion Kit (Mouse and Rat) | DHNBE-5002 | 50T                 |
| High Activity Adult Brain Enzymatic Digestion Kit (Mouse and Rat)    | DHABE-5003 | 50T                 |
| High Activity Spleen Tissue Enzymatic Digestion Kit (Mouse)          | DHSE-5004  | 50T                 |
| Other Tissue Enzymatic Digestion Kit *                               |            |                     |

**Cell Separation Kit**



HCSC-25 High-capacity separation column



Mouse CD3<sup>+</sup> Cell Separation Kit

► **Features**

- High cell viability, high purity, low stimulation, and low residue.
- Cell culture, flow cytometry, and single-cell sequencing can be carried out after separation.

**Order Information:**

| Cat No.  | Product Description                        |
|----------|--|
| HCSC-25  | High-capacity separation column(25pcs/box) |
| SSA      | Separation stand                           |
| HCS-01   | High-capacity Separator                    |
| K1301-10 | Mouse CD3 <sup>+</sup> Cell Separation Kit |


## CO<sub>2</sub> Incubator

### Provide a more stable cell culture environment

D180-P is a well-designed air-jacketed CO<sub>2</sub> incubator for all kinds of cell lines or tissues culture, including primary cells, stem cells, and so on.

It can provide a more suitable culture environment for cells by controlling temperature, CO<sub>2</sub> concentration precisely and by maintaining a high relative humidity (rH). The built-in 140°C dry-heat sterilization program helps to prevent bacterial contamination. All the data and 7-days performance trend graphing are visible on the 7-inch touch screen, and stored data can be download through USB.

With reliable temperature stability, pure environment, and simple operation, D180-P will be an indispensable device in the fields of cancer research, stem cell research, immunology research, pharmaceutical research, and neurobiology research.



### Product characteristics

#### Stable environment for cell culture

- High temperature resistant infrared (IR) CO<sub>2</sub> sensor enables real-time monitoring the change of CO<sub>2</sub> concentration.
- Six-sided air jacketed heating method with Pt1000 platinum thermistor significantly shorten the recover time to target temperature.
- Standard humidity sensor monitors real-time humidity changes.

#### Contamination control of Inner chamber

- In-chamber HEPA filter combined with an active fan system provides ISO 5 (Class 100) air cleanliness.
- 140°C high temperature sterilization provides overnight cycle operation to eliminate potential contamination.

## Intelligent operation

- 7-inch touch screen provides complete data visibility including three core parameters (Temp. , CO<sub>2</sub> and rH) and alarm records to monitor incubator operation in real time.
- 8000 operation or alarm records ensure the safety of samples and incubator.
- Customized alarm settings allow users to modify upper and lower alarm limits.
- Within 7 days parameters curve display (Temp. , CO<sub>2</sub> and rH) provides real time monitoring.



High Efficiency Particulate Air Filter System



Stainless Steel Waterpan

## ► Technical Parameters

|                               |                               |   |
|-------------------------------|-------------------------------|---|
| Temperature                   | Control & Range               | ±0.1°C, 5°C above ambient to 50°C               |
|                               | Heating Method                | Air Jacketed                                    |
|                               | Uniformity                    | ±0.3°C @ 37°C                                   |
|                               | Alarm Setting                 | ±5°C  |
|                               | Sensor                        | Pt1000 platinum thermistor                      |
|                               | Controller                    | PID   |
|                               | Readability & Stability       | 0.1°C   |
|                               | Sterilization Method          | 140°C dry heat sterilization                    |
| CO <sub>2</sub> Concentration | Control & Range               | Better than ±0.1% ; 0-20%                       |
|                               | Sensor                        | High temperature resistant infrared (IR) sensor |
|                               | Inlet Pressure                | 1.0 bar (recommended)                           |
|                               | Uniformity                    | ±0.3%   |
|                               | Readability & Stability       | 0.10%   |
| Humidity                      | Alarm Setting                 | ±1%   |
|                               | Humidity Pan Volume           | ~3.8 L  |
|                               | Display                       | 1% RH   |
| Shelves                       | Dimensions                    | 47.0×47.0 cm                                    |
|                               | Surface Area (Standard & Max) | 214.6 cm <sup>2</sup> ; 220.7 cm <sup>2</sup>   |
|                               | Number (Standard & Max)       | 3 ; 16  |
| Size & Weight                 | Interior Volume               | 175 L   |
|                               | Interior Dimensions           | 51.0 (W) × 54.7 (D) × 67.5 (H) cm               |
|                               | Exterior Dimensions           | 66.0 (W) × 67.0 (D) × 97.0 (H) cm               |
|                               | Net Weight                    | 110 kg (242.5 lbs)                              |

## ► Order Information

| Cat No. | Product Description  |
|---------|--|
| D180-P  | CO <sub>2</sub> incubator (Equipped with humidity display and high temperature sterilization function), 220V-240V VAC 50/60 Hz |
| D181-P  | CO <sub>2</sub> incubator (Equipped with humidity display and high temperature sterilization function), 120V-120V VAC 50/60 Hz |

## Automated Cell Counters

RWD Automated Cell Counter adopts Trypan Blue staining method, intelligent image recognition and advanced optical imaging technology to achieve accurate and automatic cell analysis. The algorithm can identify both live and dead cells within 9 seconds. Information such as cell counts, viability, diameter can be obtained by just one-click. With fluorescence cube, C100 can display images in BF and FL channels simultaneously and clearly present cell counts and cell morphology. It is an optimal device to analyze cell viability and transfection efficiency in immunology and vaccine development, cell therapy, tumor research, stem cell, metabolism research and others.

BF



C100-SE Automated Cell Counter

FL



C100 Automated Cell Counter



### Product Characteristics

#### High accuracy

- The self-developed algorithm can identify clustered cells.
- Intelligently cell marks contribute to precise cell recognition.
- RBC counting.
- Low autofluorescence and high-transmit filter enable clear imaging.

#### Time-saving

- User-friendly UI and built-in user guide interfaces ensure easy operation.
- Up to 1000 counting results can be stored. No need another device for extra storage.

#### User-friendly

- Autofocus and auto exposure.
- Powerful computing capability ensures high efficiency.

#### Compact design

- Compact size fits into different placements.
- Removable modular FL cubes allow easy switches among various scenes.


### Accessories and consumables



Disposable Cell Counting Slide



FL Cube



- Clustered cell recognition algorithm



- Clear fluorescence imaging



- Pre-dilution calculator



- Cell splitting calculator

## ► Technical Parameters

| Parameters                        | Description  |  |
|-----------------------------------|--|--|
| Model                             | C100-SE  | C100   |
| Optional fluorescence application | N/A  | <ul style="list-style-type: none"> <li>● (Ex :375/28nm ; Em:460/50nm)</li> <li>● (Ex :480/30nm ; Em:535/40nm)</li> <li>● (Ex :540/25nm ; Em:620/60nm)</li> </ul> |
| Amplification Factor              | 2.5×, 5.0 Megapixel  |  |
| Focus Method                      | Manual focusing & Autofocusing                             |  |
| Counting Area                     | 2.15 mm×1.62 mm  |  |
| Cell Type                         | Cell lines, stem cells, primary cells, pollens, beer yeast |  |
| Cell Size Range                   | 4-60 μm (Optimal: 7-60 μm)                                 |  |
| Cell Concentration Range          | 10 <sup>4</sup> -10 <sup>7</sup> cells/mL                  |  |
| Cell counting Time                | Less than 9s   |  |
| Historical Data Storage           | 1000 counting reports and images at most                   |  |
| Languages                         | Chinese and English  |  |
| Dimensions                        | 212 mm(W)* 264 mm(H)* 165 mm (D)                           |  |
| Weight                            | ~3.0kg   | ~3.2kg   |

| FL Cube    | Excitation | Emission | General Dyes   |
|------------|------------|----------|--|
| FL Cube-01 | 375/28nm   | 460/50nm | DAPI, Hoechst, BFP   |
| FL Cube-02 | 480/30nm   | 535/40nm | Acridine Orange (AO), Calcein-AM, Alexa fluor 488, SYBR Green, FITC, GFP |
| FL Cube-03 | 540/25nm   | 620/60nm | Propidium Iodide (PI), DsRed, dTomato, RFP                               |

## Order Information

| Cat No.    | Product Description              |
|------------|----------------------------------|
| C100-SE    | Automated Cell Counter (BF)      |
| C100       | Automated Cell Counter (BF + FL) |
| FL Cube-01 | ● (Ex :375/28nm ; Em:460/50nm)   |
| FL Cube-02 | ● (Ex :480/30nm ; Em:535/40nm)   |
| FL Cube-03 | ● (Ex :540/25nm ; Em:620/60nm)   |
| DS-50      | Disposable Slide (50 pcs/box)    |

## Microcentrifuges

### M1324R Refrigerated Microcentrifuge

RWD M1324R microcentrifuge is a refrigerated model, it can work stably and quietly at the maximum speed of 15000rpm. Even at the highest speed, it will still make sure the stable temperature centrifugation of samples at -10°C~40°C. Optional rotors for different needs (PCR tubes, 1.5/2 mL EP tubes, etc.).

Small footprint and low noise, touch screen interaction, all your user-defined preset programs can be called by one click. All the functions are designed according to your needs, which simplifies centrifugation.



M1324R Microcentrifuge

### Product characteristics

#### Efficient

- Quick precooling. From RT to 4 °C in only 8 minutes.
- The seamlessly fit for inner cavity winding technology enables high-efficient energy transfer.

#### Stable and Safe

- High-density foaming technology benefits energy-saving and heat-preservation.
- Standard aerosol-tight rotor with aerospace-grade aluminum alloy withstands high temperature sterilization and avoids biohazards.
- Qualified ultrasonic non-destructive test, Schenck balancing test and extreme blasting test to protect personnel safety.
- Real-time rotor dynamic balancing monitoring system and multiple alarm functions guarantee safe operation.

#### Simple

- Combination of touch screen and function keys provides quick and convenient setting and parameter use.
- 12 custom programs enable immediate centrifuge operation in multiple scenarios.
- Magnetic condensed water collection pan to eliminate inner water accumulation.
- Mechanical lock design for the centrifuge lid enables easy opening and closing.
- Quick-Lock rotor, closes with only a 1/2 turn.

#### Compact

- One-pipette-width compact design.
- Lid hinge design with suspending hovering optionally.

Self-locking lid

Standard  
Aerosol-tight rotor

Combination of buttons  
and touch interaction

Small, quiet and durable



## ► Technical Parameters

| Parameters          | Description       |
|---------------------|-------------------|
| Max. Speed          | 15000rpm          |
| Max. RCF            | 21130×g           |
| Max. Capacity       | 2mL×24            |
| Display             | Touch screen      |
| Program             | 12                |
| Temperature Control | -10°C to 40°C     |
| Noise Level         | <54dB (A)         |
| Dimension (W*H*D)   | 304mm*304mm*517mm |
| Weight (with rotor) | 33kg              |

## ► Order Information

| Model    | Product Description   |
|----------|---|
| M1324R   | Microcentrifuge, Refrigerated, Rotor with Aerosol-tight lid |
| M-F24G   | Aerosol-tight Lid, 24×1.5/2 mL tube rotor                   |
| M-F24    | 24×1.5/2 mL tube rotor                                      |
| M-F4PCR  | 4×PCR tube rotor  |
| M-F18Kit | 18×Spin column rotor  |
| M-F24QG  | Aerosol-tight Lid, 24×1.5/2.0 mL tube Quick- Lock rotor     |
| M-F24Q   | 24×1.5/2 mL tube Quick- Lock rotor                          |

## M1324 Ventilated Microcentrifuge

RWD M1324 microcentrifuge is a ventilated model, it can satisfy most of common standard applications. And even at the highest speed, it can still work stably and quietly. It is the best choice for laboratory separation and provides great flexibility with its optional rotors for different needs (PCR tubes, 1.5/2 mL EP tubes, etc.).



M1324 Microcentrifuge

- Max. speed 15,800 rpm
- Max. RCF 23,444 × g
- 4.3-inch touch screen control
- Low noise and high-speed accuracy
- Multi-step speed up and down available (10 acceleration and braking ramps)



Quick-Lock rotor

- Quick-Lock rotor, closes with only a 1/8 turn
- Automatic unbalance detection
- Self-locking device, over-speed safety device, and automatic alarm device
- Inching function with brief centrifuging

### ► Technical Parameters

| Parameters          | Description          |
|---------------------|----------------------|
| Max. Speed          | 15800rpm (23444 × g) |
| Max. Capacity       | 2mL × 24             |
| Noise Level         | ≤ 51dB(A)            |
| Display             | Touch screen         |
| Program             | 48                   |
| Dimension (W*H*D)   | 262mm*282mm*385mm    |
| Weight (with rotor) | 22kg                 |

### ► Order Information

| Model    | Product Description                                       |
|----------|---|
| M1324    | Microcentrifuge, Ventilated, Rotor with Aerosol-tight lid |
| M-F24G   | Aerosol-tight Lid, 24 × 1.5/2 mL tube rotor               |
| M-F24    | 24 × 1.5/2 mL tube rotor                                  |
| M-F4PCR  | 4 × PCR tube rotor  |
| M-F18Kit | 18 × Spin column rotor                                    |
| M-F24QG  | Aerosol-tight Lid, 24 × 1.5/2.0 mL tube Quick-lock rotor  |
| M-F24Q   | 24 × 1.5/2 mL tube Quick-Lock rotor                       |

# Animal Surgical Instrument & Matrices

RWD surgical instruments cover surgical scissors, surgical forceps, orthopedic instruments, scalpels and blades, vascular forceps, needle forceps, with more than 1,000 kinds of models.

## Surgery Instrument Kit

### SP0001-G General Surgery Instrument Kit for Rodents



| Cat.No.   | Product Description  | Qty |
|-----------|--|-----|
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm          | 1   |
| S11002-08 | VANNAS Spring Scissors (Triangular)-S/S Cvd/8*1.5mm/8cm          | 1   |
| S12009-11 | IRIS-Fine Scissors (Round Type)-S/S Str/25*5.6mm/11.5cm          | 1   |
| S12010-11 | IRIS-Fine Scissors (Round Type)-S/S Cvd/23.5*6mm/11.5cm          | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm            | 1   |
| S12004-09 | IRIS-Fine Scissors (Round Type)-S/S Cvd/24*4.6mm/9.5cm           | 1   |
| F11001-11 | Micro Forceps-Str, head length 30mm, Tip 0.2*0.12mm, 11cm        | 1   |
| F11004-11 | Micro Forceps-Cvd, head length 26.8mm, Tip 0.2*0.12mm, 11cm      | 1   |
| F12010-10 | Dressing Forceps-Str, Tip width 1.8mm, teeth length 15mm, 10.5cm | 1   |
| F12011-10 | Dressing Forceps-Cvd, Tip width 1.9mm, teeth length 15mm, 10.5cm | 1   |
| F21001-12 | HALSTED Artery Forceps-Str/23.5*1mm/12.5cm                       | 2   |
| F21002-12 | HALSTED Artery Forceps-Cvd/24*1mm/12.5cm                         | 2   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                     | 2   |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                       | 2   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm     | 1   |
| S13029-14 | MAYO-STILLE Dissecting & Operating Scissors-B/B Str/37*9mm/14cm  | 1   |
| S15001-11 | SPENCER Ligature Scissors (Slender Type)-27*6mm/11cm             | 1   |
| S21010-15 | BOHLER Rongeurs (Double)-Light Cvd/2*19mm/15cm                   | 1   |
| S22004-11 | Bone Cutters with Flat Blades (SGL)-15*6.8mm/11cm                | 1   |
| S14014-12 | Operating Scissors (Round Type)-S/S Str/37.6*8.5mm/12.5cm        | 1   |
| S31011-01 | 11# Scalpel Blades (Box of 100pcs)                               | 1   |
| S32003-12 | Scalpel Handles with Ruler 3# Solid-12.5cm                       | 1   |
| S33009-06 | Scalpel Blade Remover-6.5cm                                      | 1   |
| F13019-12 | ADSON 1x2 Teeth Tissue Forceps, Tip width 1.8mm, 12cm            | 1   |
| R21027-12 | STEVENS Hooks, 1 Angled Tooth, 5mm depth, 12.5cm                 | 1   |
| R22009-01 | ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm               | 1   |
| R22005-45 | 3x3 Teeth Retractors-Blunt, 30mm Spread, 5cm                     | 1   |
| R51003-11 | Dual-end Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm      | 1   |
| F35204-70 | NA Nylon Sutures w/Needle-△3/8/2.5×5/30cm/7-0 (50/Box)           | 0.2 |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)         | 0.2 |
| R31006-06 | SS Micro Clamps-Cvd/L*W 6.6*1mm/17mm                             | 5   |
| R31005-06 | SS Micro Clamps-Str/L*W 6*1mm/19mm                               | 5   |
| R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm       | 1   |
| SP0001-P  | Instrument Storage Portfolio                                     | 1   |

### SP0003-R Microsurgery Instrument Kit for Rat, Cat, Rabbit etc.



| Cat.No.   | Product Description  | Qty |
|-----------|--|-----|
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm        | 1   |
| S12010-11 | IRIS-Fine Scissors (Round Type)-S/S Cvd/23.5*6mm/11.5cm      | 1   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm | 1   |
| F11001-11 | Micro Forceps-Str, head length 30mm, Tip 0.2*0.12mm, 11cm    | 1   |
| R31005-06 | SS Micro Clamps-Str/L*W 6*1mm/19mm                           | 5   |
| R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm   | 1   |
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm      | 1   |
| S11002-08 | VANNAS Spring Scissors (Triangular)-S/S Cvd/8*1.5mm/8cm      | 1   |
| F21001-12 | HALSTED Artery Forceps-Str/23.5*1mm/12.5cm                   | 1   |
| F21002-12 | HALSTED Artery Forceps-Cvd/24*1mm/12.5cm                     | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)     | 0.2 |
| S32003-12 | 11# Scalpel Blades (Box of 100pcs)                           | 1   |
| S31011-01 | Scalpel Handles with Ruler 3# Solid-12.5cm                   | 1   |
| SP0001-P  | Instrument Storage Portfolio                                 | 1   |

### SP0003-M Microsurgery Instrument Kit for Mice



| Cat.No.   | Product Description  | Qty |
|-----------|--|-----|
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm        | 1   |
| S12004-09 | IRIS-Fine Scissors (Round Type)-S/S Cvd/24*4.6mm/9.5cm       | 1   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm | 1   |
| F11001-11 | Micro Forceps-Str, head length 30mm, Tip 0.2*0.12mm, 11cm    | 1   |
| R31005-04 | SS Micro Clamps-Str/L*W 4*0.75mm/16mm                        | 5   |
| R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm   | 1   |
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm      | 1   |
| S11002-08 | VANNAS Spring Scissors (Triangular)-S/S Cvd/8*1.5mm/8cm      | 1   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                 | 1   |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                   | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)     | 0.2 |
| S32003-12 | 11# Scalpel Blades (Box of 100pcs)                           | 1   |
| S31011-01 | Scalpel Handles with Ruler 3# Solid-12.5cm                   | 1   |
| SP0001-P  | Instrument Storage Portfolio                                 | 1   |

**SP0005-R**

**Brain Microinjection Surgery Instrument Kit for Rat**



| Cat.No.   | Product Description  | Qty |
|-----------|--|-----|
| S31011-01 | Scalpel Handles with Ruler 3# Solid-12.5cm                       | 1   |
| S32003-12 | 11# Scalpel Blades (Box of 100pcs)                               | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm            | 1   |
| S12006-10 | IRIS-Fine Scissors (Round Type)-S/S Cvd/22*4.5mm/10.5cm          | 1   |
| F12010-10 | Dressing Forceps-Str, Tip width 1.8mm, teeth length 15mm, 10.5cm | 2   |
| F12011-10 | Dressing Forceps-Cvd, Tip width 1.9mm, teeth length 15mm, 10.5cm | 2   |
| F21001-12 | HALSTED Artery Forceps-Str/23.5*1mm/12.5cm                       | 2   |
| F21002-12 | HALSTED Artery Forceps-Cvd/24*1mm/12.5cm                         | 2   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm     | 1   |
| S14014-12 | Operating Scissors (Round Type)-S/S Str/37.6*8.5mm/12.5cm        | 1   |
| F35204-70 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                     | 0.2 |
| R22009-01 | ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm               | 1   |
| R51003-11 | Dual-end Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm      | 1   |
| SP0001-P  | Instrument Storage Portfolio                                     | 1   |

**SP0005-M**

**Brain Microinjection Surgery Instrument Kit for Mice**



| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| S31011-01 | Scalpel Handles with Ruler 3# Solid-12.5cm                                  | 1   |
| S32003-12 | 11# Scalpel Blades (Box of 100pcs)  | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   |
| S12004-09 | IRIS-Fine Scissors (Round Type)-S/S Cvd/24*4.6mm/9.5cm                      | 1   |
| S14014-10 | Operating Scissors (Round Type)-S/S Str/32.7*8.45mm/10.5cm                  | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 2   |
| F12006-10 | IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm | 2   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 2   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm                | 1   |
| F35204-70 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 0.2 |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                                  | 2   |
| R51003-11 | Dual-end Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm                 | 1   |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0007-M**

**Univeral Dissecting Surgery Instrument Kit for Mice**



| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   |
| S12004-09 | IRIS-Fine Scissors (Round Type)-S/S Cvd/24*4.6mm/9.5cm                      | 1   |
| S13003-11 | STEVENS Fine Dissecting Scissors (Flat Type)-B/B Str/25*4mm/11.5cm          | 1   |
| S13004-11 | STEVENS Fine Dissecting Scissors (Flat Type)-B/B Cvd/28*5mm/11.5cm          | 1   |
| S15001-09 | SPENCER Ligature Scissors (Slender Type)-20.3*4mm/9cm                       | 1   |
| S32003-12 | 11# Scalpel Blades (Box of 100pcs)  | 1   |
| S31015-01 | 15# Scalpel Blades (Box of 100pcs)  | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 2   |
| F12006-10 | IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm | 2   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 2   |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                                  | 2   |
| F13029-10 | IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm                  | 1   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm                | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)                    | 0.2 |
| R22005-45 | 3x3 Teeth Retractors-Blunt, 30mm Spread, 5cm                                | 1   |
| R31005-04 | SS Micro Clamps-Str/L*W 4*0.75mm/16mm                                       | 5   |
| R42002-12 | Spinal Cord Hook - Tip Dia. 3mm/12cm  | 1   |
| R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm                  | 1   |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0007-R**

**Univeral Dissecting Surgery Instrument Kit for Rat**



| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm               | 1   |
| S12006-10 | IRIS-Fine Scissors (Round Type)-S/S Cvd/22*4.5mm/10.5cm             | 1   |
| S12010-11 | IRIS-Fine Scissors (Round Type)-S/S Cvd/23.5*6mm/11.5cm             | 1   |
| S13005-14 | METZENBAUM Dissecting Scissors (Round Type)-B/B Str/40*6.3mm/14.5cm | 1   |
| S13006-14 | METZENBAUM Dissecting Scissors (Round Type)-B/B Cvd/44*7.5mm/14.5cm | 1   |
| F21001-12 | HALSTED Artery Forceps-Str/23.5*1mm/12.5cm                          | 2   |
| F21002-12 | HALSTED Artery Forceps-Cvd/24*1mm/12.5cm                            | 2   |
| F12010-10 | Dressing Forceps-Str, Tip width 1.8mm, teeth length 15mm, 10.5cm    | 2   |
| F12011-10 | Dressing Forceps-Cvd, Tip width 1.9mm, teeth length 15mm, 10.5cm    | 2   |
| F13029-10 | IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm          | 1   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm        | 0.2 |
| S15001-11 | SPENCER Ligature Scissors (Slender Type)-27*6mm/11cm                | 1   |
| S32003-12 | 11# Scalpel Blades (Box of 100pcs)                                  | 1   |
| S31015-01 | 15# Scalpel Blades (Box of 100pcs)                                  | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)            | 1   |
| R22009-01 | ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm                  | 1   |
| R31005-06 | SS Micro Clamps-Str/L*W 6*1mm/19mm                                  | 5   |
| R42002-12 | Spinal Cord Hook - Tip Dia. 3mm/12cm                                | 1   |
| R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm          | 1   |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0008-M Cardio-cerebral Perfusion Surgery Instrument Kit for Mice**



| Cat.No.   | Product Description   | Qty | Cat.No.   | Product Description  | Qty |
|-----------|---|-----|-----------|--|-----|
| S14014-11 | Operating Scissors (Round Type)-S/S Str/32.7*8.45mm/10.5cm                  | 1   | F35204-70 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                 | 0.2 |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   | R21027-12 | STEVENS Hooks, 1 Angled Tooth , 5mm depth, 12.5cm            | 1   |
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm                     | 1   | R22005-45 | 3x3 Teeth Retractors-Blunt, 30mm Spread, 5cm                 | 1   |
| F12007-10 | IRIS Dissecting Forceps-Large Cvd, Tip width 0.9mm, teeth length 10mm, 10cm | 1   | F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm | 1   |
| F12010-09 | Dressing Forceps-Str, Tip width 1.8mm, teeth length 13.5mm, 9.5cm           | 2   | R31005-04 | SS Micro Clamps-Str/L*W 4*0.75mm/16mm                        | 5   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 2   | R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm   | 1   |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                                  | 2   | R42002-12 | Spinal Cord Hook - Tip Dia. 3mm/12cm                         | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)                    | 0.2 | SP0001-P  | Instrument Storage Portfolio                                 | 1   |

**SP0008-R**  
**Cardio-cerebral Perfusion Surgery Instrument Kit for Rat**


| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| S14014-12 | Operating Scissors (Round Type)-S/S Str/37.6*8.5mm/12.5cm                   | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   |
| S11002-08 | VANNAS Spring Scissors (Triangular)-S/S Cvd/8*1.5mm/8cm                     | 1   |
| F12007-10 | IRIS Dissecting Forceps-Large Cvd, Tip width 0.9mm, teeth length 10mm, 10cm | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 2   |
| F21001-12 | HALSTED Artery Forceps-Str/23.5*1mm/12.5cm                                  | 2   |
| F21002-12 | HALSTED Artery Forceps-Cvd/24*1mm/12.5cm                                    | 2   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)                    | 0.2 |
| F35204-70 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 0.2 |
| R21027-12 | STEVENS Hooks, 1 Angled Tooth , 5mm depth, 12.5cm                           | 1   |
| R22005-45 | 3x3 Teeth Retractors-Blunt, 30mm Spread, 5cm                                | 1   |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm                | 1   |
| R31005-06 | SS Micro Clamps-Str/L*W 6*1mm/19mm  | 5   |
| R42002-12 | Spinal Cord Hook - Tip Dia. 3mm/12cm  | 1   |
| R34001-14 | Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm                  | 1   |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0011-B MCAO**  
**Surgical Instrument Kit for Rat & Mice**


| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| S32003-12 | Scalpel Handles with Ruler 3# Solid-12.5cm                                  | 1   |
| S31015-01 | 15# Scalpel Blades (Box of 100pcs)  | 1   |
| S14014-15 | Operating Scissors (Round Type)-S/S Str/48*10mm/15.5cm                      | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   |
| S12006-10 | IRIS-Fine Scissors (Round Type)-S/S Cvd/22*4.5mm/10.5cm                     | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 1   |
| F12006-10 | IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm | 1   |
| F13029-10 | IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm                  | 1   |
| F13030-10 | IRIS 1x2 Teeth Tissue Forceps-Light Cvd, Tip width 1.3mm, 10.5cm            | 1   |
| R21013-01 | MANNERFELT Hooks, 2 Teeth (7mm), Blunt, 15.5cm*                             | 4   |
| R31005-06 | SS Micro Clamps-Str/L*W 6*1mm/19mm  | 5   |
| F11028-13 | Micro Forceps-Cvd, head length 38mm, Tip 0.1*0.06mm, 13cm                   | 1   |
| F11027-13 | Micro Forceps-Str, head length 37.5mm, Tip 0.1*0.06mm, 13cm                 | 1   |
| F11003-13 | Micro Forceps-45°Angled, head length 30mm, Tip 0.15*0.1mm, 13cm             | 1   |
| F11002-11 | Micro Forceps-Str, head length 36.5mm, Tip 0.15*0.1mm, 11cm                 | 1   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 1   |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                                  | 1   |
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm                     | 1   |
| F31047-14 | OLSEN-HEGAR Needle Holders with Scissors-Str, 17.3*1.7mm/14cm               | 1   |
| F35204-70 | NA Nylon Sutures w/Needle-△3/8/2.5×5/30cm/7-0 (50/Box)                      | 0.2 |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0009-R General Surgery Instrument Kit for Rat**


| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| F11020-11 | Micro Forceps-Str, head length 36mm, Tip 0.05*0.01mm, 11cm                  | 1   |
| F11021-11 | Micro Forceps-45°Angled, head length 30mm, Tip 0.1*0.06mm, 11cm             | 1   |
| S32003-12 | Scalpel Handles with Ruler 3# Solid-12.5cm                                  | 1   |
| S31011-01 | 11# Scalpel Blades (Box of 100pcs)  | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)                    | 0.2 |
| F35204-70 | NA Nylon Sutures w/Needle-△3/8/2.5×5/30cm/7-0 (50/Box)                      | 0.2 |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm                | 1   |
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm                     | 1   |
| S11002-08 | VANNAS Spring Scissors (Triangular)-S/S Cvd/8*1.5mm/8cm                     | 1   |
| R22009-01 | ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm                          | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 1   |
| F12006-10 | IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm | 1   |
| F12007-10 | IRIS Dissecting Forceps-Large Cvd, Tip width 0.9mm, teeth length 10mm, 10cm | 1   |
| F13029-10 | IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm                  | 1   |
| S12009-11 | IRIS-Fine Scissors (Round Type)-S/S Str/25*5.6mm/11.5cm                     | 1   |
| S12010-11 | IRIS-Fine Scissors (Round Type)-S/S Cvd/23.5*6mm/11.5cm                     | 1   |
| F21001-12 | HALSTED Artery Forceps-Str/23.5*1mm/12.5cm                                  | 1   |
| F21002-12 | HALSTED Artery Forceps-Cvd/24*1mm/12.5cm                                    | 1   |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0009-M**  
**General Surgery Instrument Kit for Mouse**


| Cat.No.   | Product Description   | Qty |
|-----------|---|-----|
| F11020-11 | Micro Forceps-Str, head length 36mm, Tip 0.05*0.01mm, 11cm                  | 1   |
| F11021-11 | Micro Forceps-45°Angled, head length 30mm, Tip 0.1*0.06mm, 11cm             | 1   |
| S32003-12 | Scalpel Handles with Ruler 3# Solid-12.5cm                                  | 1   |
| S31011-01 | 11# Scalpel Blades (Box of 100pcs)  | 1   |
| F35401-50 | NA Terylene Sutures w/Needle-○3/8/3×10/90cm/5-0 (50/Box)                    | 0.2 |
| F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm                | 1   |
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm                     | 1   |
| S11002-08 | VANNAS Spring Scissors (Triangular)-S/S Cvd/8*1.5mm/8cm                     | 1   |
| R22029-04 | Colibri Eye Specula-15mm spread, 4cm  | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 1   |
| F12006-10 | IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm | 1   |
| F12007-10 | IRIS Dissecting Forceps-Large Cvd, Tip width 0.9mm, teeth length 10mm, 10cm | 1   |
| F13029-10 | IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm                  | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   |
| S12004-09 | IRIS-Fine Scissors (Round Type)-S/S Cvd/24*4.6mm/9.5cm                      | 1   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 1   |
| F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                                  | 1   |
| SP0001-P  | Instrument Storage Portfolio  | 1   |

**SP0010-B Basic Surgical Instrument Kit for Rat & Mice**


| Cat.No.   | Product Description   | Qty | Cat.No.   | Product Description  | Qty |
|-----------|---|-----|-----------|--|-----|
| S11001-08 | VANNAS Spring Scissors (Triangular)-S/S Str/8*1.5mm/8cm                     | 1   | F22003-10 | HARTMAN Mosquito Forceps-Cvd/20.8*1mm/10cm                   | 1   |
| S12005-10 | IRIS-Fine Scissors (Round Type)-S/S Str/23*5mm/10.5cm                       | 1   | F31047-12 | OLSEN-HEGAR Needle Holders with Scissors-Str, 10*2.15mm/12cm | 1   |
| F12005-10 | IRIS Dissecting Forceps-Str, Tip width 1mm, teeth length 14mm, 10cm         | 1   | S33007-12 | ZIEGLER Knives-8mm Cutting Edge/1.5mm Width/12cm             | 1   |
| F12006-10 | IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm | 1   | R31005-06 | SS Micro Clamps-Str/L*W 6*1mm/19mm                           | 1   |
| F11020-11 | Micro Forceps-Str, head length 36mm, Tip 0.05*0.01mm, 11cm                  | 1   | R31006-06 | SS Micro Clamps-Cvd/L*W 6.6*1mm/17mm                         | 1   |
| F22002-10 | HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm                                | 1   | SP0001-P  | Instrument Storage Portfolio                                 | 1   |

## Micro Spring Scissors



### VANNAS Spring Scissors

Cat.No.: S11001-08

- Triangular | Straight
- Sharp&Sharp
- Cutting Edge (Length): 8.0mm
- Width: 1.5mm | 8.0cm
- Price: \$99.00



Cat.No.: S11002-08

- Triangular | Curved
- Sharp&Sharp
- Cutting Edge (Length): 8.0mm
- Width: 1.5mm | 8.0cm
- Price: \$102.00



Cat.No.: S11003-08

- Triangular | Angled
- Sharp&Sharp
- Cutting Edge (Length): 6.3mm
- Width: 1.35mm | 8.0cm
- Price: \$106.00



### VANNAS Spring Scissors

Cat.No.: S11035-08

- Triangular | Straight
- Sharp&Sharp
- Cutting Edge (Length): 3.5mm
- Width: 1.5mm | 8.0cm
- Price: \$92.00



Cat.No.: S11036-08

- Flat | Straight
- Sharp&Sharp
- Cutting Edge (Length): 8.2mm
- Width: 1.5mm | 8.0cm
- Price: \$118.00



Cat.No.: S11037-08

- Flat | Curved
- Sharp&Sharp
- Cutting Edge (Length): 8.0mm
- Width: 1.5mm | 8.0cm
- Price: \$118.00



### VANNAS Spring Scissors

Cat.No.: S11007-12

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 19mm
- Width: 4.25mm | 12.5cm
- Price: \$45.5.00



Cat.No.: S11019-14

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 14.5mm
- Width: 3.2mm | 14cm
- Price: \$121.00

Cat.No.: S11022-14

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 14.5mm
- Width: 3.2mm | 14cm
- Price: \$94.00



### VANNAS Spring Scissors

Cat.No.: S11014-12

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 5.5mm
- Width: 1.9mm | 12.5cm
- Price: \$103.00



Cat.No.: S11008-12

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 14mm
- Width: 0.2mm | 12cm
- Price: \$48.00

## Fine Scissors



Cat.No.: S12004-09

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 24mm
- Width: 4.6mm | 9.5cm
- Price: \$25.00

### IRIS-Fine Fine Scissors

Cat.No.: S12003-09

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 24.5mm
- Width: 5mm | 9cm
- Price: \$24.00



Cat.No.: S12005-10

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 23mm
- Width: 5mm | 10.5cm
- Price: \$24.00



Cat.No.: S12006-10

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 22mm
- Width: 4.5mm | 10.5cm
- Price: \$25.00



Cat.No.: S12011-12

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 31mm
- Width: 5.8mm | 12.5cm
- Price: \$24.00



Cat.No.: S18005-10

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 20.5mm
- Width: 5mm | 10.5cm
- Price: \$57.00



Cat.No.: S18005-09

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 20.5mm
- Width: 4.5mm | 9cm
- Price: \$56.00



Cat.No.: S18002-09

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 21.5mm
- Width: 5mm | 9cm
- Price: \$19.00



Cat.No.: S18002-10

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 23mm
- Width: 5mm | 10cm
- Price: \$22.80



Cat.No.: S18004-11

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 22mm
- Width: 5mm | 11.5cm
- Price: \$31.00



Cat.No.: S18004-10

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 20mm
- Width: 5mm | 10.5cm
- Price: \$36.00

### IRIS Fine Scissors

Cat.No.: S12010-11

- Curved
- Sharp&Sharp
- Cutting Edge (Length): 23.5mm
- Width: 6mm | 11.5cm
- Price: \$25.00



Cat.No.: S12031-11

- Full Curved
- Sharp&Sharp
- Cutting Edge (Length): 29mm
- Width: 6mm | 1.5cm
- Price: \$14.50

### Fine Scissors

Cat.No.: S12034-12

- Straight
- Sharp&Sharp
- Cutting Edge (Length): 17mm
- Width: 6mm | 12.5cm
- Price: \$12.50

## Dissecting & Operating Scissors

- Cat.No.: S13020-11**
- Straight
  - Sharp&Sharp
  - Cutting Edge (Length): 20mm
  - Width: 5mm | 12cm
  - Price: \$24.00

- Cat.No.: S13005-14**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 40mm
  - Width: 6.3mm | 14.5cm
  - Price: \$24.00

- Cat.No.: S13035-16**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 45mm
  - Width: 9.5mm | 16.5cm
  - Price: \$31.00



### Delicate Scissors

- Cat.No.: S13001-10**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 25.5mm
  - Width: 6mm | 10cm
  - Price: \$28.00

- Cat.No.: S14001-13**
- Straight
  - Sharp&Blunt
  - Cutting Edge (Length): 40.5mm
  - Width: 9.4mm | 13cm
  - Price: \$17.00

- Cat.No.: S14001-15**
- Straight
  - Sharp&Blunt
  - Cutting Edge (Length): 50mm
  - Width: 10mm | 15.5cm
  - Price: \$19.00




### Operating Scissors

- Cat.No.: S14001-10**
- Straight
  - Sharp&Blunt
  - Cutting Edge (Length): 33mm
  - Width: 7.5mm | 10.5cm
  - Price: \$14.90

- Cat.No.: S14014-12**
- Straight
  - Sharp&Sharp
  - Cutting Edge (Length): 37.6mm
  - Width: 8.5mm | 12.5cm
  - Price: \$17.00

- Cat.No.: S14014-13**
- Straight
  - Sharp&Sharp
  - Cutting Edge (Length): 45mm
  - Width: 9.8mm | 13cm
  - Price: \$18.00

- Cat.No.: S14014-15**
- Straight
  - Sharp&Sharp
  - Cutting Edge (Length): 48mm
  - Width: 10mm | 15.5cm
  - Price: \$21.00



### Operating Scissors

- Cat.No.: S14014-10**
- Straight
  - Sharp&Sharp
  - Cutting Edge (Length): 32.7mm
  - Width: 8.45mm | 10.5cm
  - Price: \$14.30

- Cat.No.: S14015-11**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 36mm
  - Width: 7mm | 11.5cm
  - Price: \$14.30

- Cat.No.: S14015-12**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 37mm
  - Width: 7mm | 12.5cm
  - Price: \$14.30

- Cat.No.: S14015-14**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 46mm
  - Width: 7.5mm | 14.5cm
  - Price: \$17.00



### Operating Scissors

- Cat.No.: S14015-10**
- Straight
  - Blunt&Blunt
  - Cutting Edge (Length): 35mm
  - Width: 7.5mm | 10.5cm
  - Price: \$13.40

## Specialist Scissors & Bone Cutters & Rongeurs



### Stitch Scissors

Cat.No.: S15006-09

- Angled
- Sharp&Blunt
- Cutting Edge (Length): 13.5mm
- Width: 4mm | 9.5cm
- Price: \$30.00



Cat.No.: S15001-09

- Slender Type
- Cutting Edge (Length): 20.3mm
- Width: 4mm | 9cm
- Price: \$21.00



### BEEBE Bone Cutters

Cat.No.: S16001-10

- Straight
- Cutting Edge (Length): 21mm
- Width: 10mm | 10.5cm
- Price: \$10.00



Cat.No.: S22004-11

- Flat Blades
- Cutting Edge (Length): 15mm
- Width: 6.8mm | 11cm
- Price: \$46.00



Cat.No.: S21020-14

- Straight
- Cutting Edge (Length): 0.7mm
- Width: 6mm | 14.5cm
- Price: \$96.50



Cat.No.: S21021-14

- Straight
- Cutting Edge (Length): 1mm
- Width: 6.5mm | 14cm
- Price: \$85.00



Cat.No.: S21023-14

- Curved
- Cutting Edge (Length): 0.7mm
- Width: 7mm | 14cm
- Price: \$200.00



Cat.No.: S21024-14

- Curved
- Cutting Edge (Length): 1mm
- Width: 6.5mm | 14cm
- Price: \$88.00



### Rongeurs

Cat.No.: S21019-14

- Straight
- Cutting Edge (Length): 0.5mm
- Width: 7.5mm | 14.5cm
- Price: \$101.50

### BEEBE Bone Cutters

Cat.No.: S22005-12

- Flat Blades
- Cutting Edge (Length): 15.2mm
- Width: 7.4mm | 12.5cm
- Price: \$23.50

## Scalpel&Micro Forceps



### Scalpel

Cat.No.: S33006-13

- Straight
- Cutting Edge (Length): 25mm
- Width: 3mm | 13cm
- Price: \$33.00



### Scalpel Handles(3#)

Cat.No.: S32003-12

- 12.5cm
- Price: \$12.30



### Scalpel Handles(4#)

Cat.No.: S32004-13

- 14cm
- Price: \$8.90



### Scalpel Blades(10#)

Cat.No.: S31010-01

- Box of 100pcs
- Price: \$25.00



### Scalpel Blades(11#)

Cat.No.: S31011-01

- Box of 100pcs
- Price: \$25.00



### Scalpel Blades(12#)

Cat.No.: S31012-01

- Box of 100pcs
- Price: \$25.00



### Scalpel Blades(15#)

Cat.No.: S31015-01

- Box of 100pcs
- Price: \$25.00



### Micro Forceps

Cat.No.: F11029-11

- Straight
- Tips: 0.14mmx0.09mm | 11cm
- Price: \$44.00



### Micro Forceps

Cat.No.: F11030-11

- 45° Angled
- Tips: 0.1mmx0.06mm | 11cm
- Price: \$55.00



### Micro Forceps

Cat.No.: F11031-11

- Curved
- Tips: 0.07mmx0.03mm | 11cm
- Price: \$53.00



### Micro Forceps

Cat.No.: F11001-11

- Straight
- Tips: 0.2mmx0.12mm | 11cm
- Price: \$24.00



Cat.No.: F11002-11

- Straight
- Tips: 0.15mmx0.1mm | 11cm
- Price: \$25.00



Cat.No.: F11020-11

- Straight
- Tips: 0.05mmx0.01mm | 11cm
- Price: \$37.00



- Cat.No.: F11003-13
- 45°Angled
- Tips: 0.15mmx0.1mm | 13cm
- Price: \$33.00

- Cat.No.: F11021-11
- 45°Angled
- Tips: 0.1mmx0.06mm | 11cm
- Price: \$44.00

**Micro Forceps**

- Cat.No.: F11003-11
- 45°Angled
- Tips: 0.3mmx0.15mm | 11cm
- Price: \$29.00



- Cat.No.: F13004-11
- Curved
- Tips: 0.3mmx0.15mm | 11.5cm
- Price: \$29.00

- Cat.No.: F11004-11
- Curved
- Tips: 0.2mmx0.12mm | 11cm
- Price: \$32.00

**Micro Forceps**

- Cat.No.: F11022-11
- Curved
- Tips: 0.1mmx0.06mm | 11.5cm
- Price: \$42.00



- Cat.No.: F11027-13
- Curved
- Tips: 0.1mmx0.06mm | 13cm
- Price: \$90.00

**Micro Forceps**

- Cat.No.: F11028-13
- Straight
- Tips: 0.1mmx0.06mm | 13cm
- Price: \$90.00

## Surgical Forceps



- Cat.No.: F12003-12
- Straight
- Tips width: 1.5mm | 12cm
- Price: \$22.00

- Cat.No.: F13029-10
- Straight
- Tips width: 1.3mm | 10.5cm
- Price: \$28.00

- Cat.No.: F13030-10
- Curved
- Tips width: 1.3mm | 10.5cm
- Price: \$30.00

- Cat.No.: F13009-11
- Straight
- Tips width: 1.6mm | 11cm
- Price: \$65.00

**ADSON Dressing Forceps**

- Cat.No.: F12002-12
- Straight
- Tips width: 0.6mm | 11cm
- Price: \$22.00



- Cat.No.: F12006-10
- Curved
- Tips width: 1.0mm | 10cm
- Price: \$25.00

- Cat.No.: F12007-10
- Large Curved
- Tips width: 0.9mm | 10cm
- Price: \$25.00

- Cat.No.: F14012-10
- Large Curved
- Tips width: 0.5mm | 10cm
- Price: \$25.00

**IRIS Dissecting Forceps**

- Cat.No.: F12005-10
- Straight
- Tips width: 1.0mm | 10cm
- Price: \$23.00



- Cat.No.: F18011-09
- Straight | 5mm Platform
- Tips : 0.4x0.9mm | 9cm
- Price: \$53.00

**McPHERSON Forceps**

## Artery Forceps & Needle Holders



Cat.No.: F21002-12

- Curved
- Cutting Edge (Length): 24mm
- Width: 1mm | 12.5cm
- Price: \$14.90

### Artery Forceps

Cat.No.: F21001-12

- Straight
- Cutting Edge (Length): 23.5mm
- Width: 1mm | 12.5cm
- Price: \$14.30



Cat.No.: F18001-10

- Angled
- Head Width : 5.4mm | 10.5cm
- Price: \$19.00

### Cover Glass Forceps

Cat.No.: F18015-10

- Angled
- Head Width: 6mm | 10.6cm
- Price: \$14.90



Cat.No.: F22003-10

- Curved
- Cutting Edge (Length): 20.8mm
- Width: 1mm | 10cm
- Price: \$18.00

Cat.No.: F22010-12

- Straight
- Head Width : 0.6mm | 12.5cm
- Price: \$14.00

Cat.No.: F22002-10

- Straight
- Cutting Edge (Length): 20mm
- Width: 0.8mm | 10.5cm
- Price: \$17.00

### Artery Forceps

Cat.No.: F22011-12

- Curved
- Head Width : 0.6mm | 12.5cm
- Price: \$30.00



Cat.No.: F31047-14

- Straight | with Scissors
- Cutting Edge (Length): 17.3mm
- Width: 1.7mm | 14cm
- Price: \$49.00

### Needle Holders

Cat.No.: F31047-12

- Straight | with Scissors
- Cutting Edge (Length): 10mm
- Width: 2.15mm | 12cm
- Price: \$46.00



### BARRAQUER Needle Holders

Cat.No.: F31017-12

- Straight | saws & lock
- Cutting Edge (Length): 10mm
- Width: 0.9mm | 12cm
- Price: \$96.00



Cat.No.: F31034-12

- Straight
- Cutting Edge (Length): 15.2mm
- Width: 1.75mm | 13cm
- Price: \$23.00

Cat.No.: F31034-14

- Straight
- Cutting Edge (Length): 17.5mm
- Width: 1.75mm | 14cm
- Price: \$24.00

### Needle Holders

Cat.No.: F31025-13

- Straight
- Cutting Edge (Length): 17.5mm
- Width: 1.75mm | 13cm
- Price: \$23.00

## Ear Punch & Micro Clamps



Cat.No.: C22002-10

- Straight
- Punch Dia. : 2mm | 10.5cm
- Price: \$35.00

### Ear Punch

Cat.No.: C22001-10

- Straight
- Punch Dia. : 1mm | 10.5cm
- Price: \$35.00



### Micro Clamps

Cat.No.: R34001-14

- Straight
- Cutting Edge (Length): 8.4mm
- Width: 5mm | 14cm
- Price: \$41.00



Cat.No.: R31005-04

- Straight
- Cutting Edge (Length): 4mm
- Width: 0.75mm | 16mm
- Price: \$36.00



Cat.No.: R31005-06

- Straight
- Cutting Edge (Length): 6mm
- Width: 1.0mm | 19mm
- Price: \$36.00



Cat.No.: R31005-08

- Straight
- Cutting Edge (Length): 8mm
- Width: 2.0mm | 18mm
- Price: \$36.00



Cat.No.: R31005-10

- Straight
- Cutting Edge (Length): 10mm
- Width: 2.0mm | 20mm
- Price: \$36.00



Cat.No.: R31006-04

- Curved
- Cutting Edge (Length): 6.3mm
- Width: 1.0mm | 17mm
- Price: \$45.00



Cat.No.: R31006-06

- Curved
- Cutting Edge (Length): 6.6mm
- Width: 1.0mm | 17mm
- Price: \$45.00

The vascular clip of the R31005 and R31006 series must be used with the R34001-14 vascular clip holder

## Hooks & Retractors



### Hooks

Cat.No.: R21027-12

- 1 Angled Tooth
- 5mm depth | 12.5cm
- Price: \$19.00



### Retractors

Cat.No.: R22029-03

- 15mm spread | 3cm
- Price: \$19.00



### Retractors

Cat.No.: R22005-45

- 30mm Spread | 5cm
- Price: \$36.00



### Retractors

Cat.No.: R22011-50


- 30mm Spread | 5.5cm
- Price: \$45.00



### Retractors

Cat.No.: R22012-45

- 22mm Spread | 4.5cm
- Price: \$17.00



### Retractors

Cat.No.: R22006-55

- 36mm Spread | 5.5cm
- Price: \$38.00



### Retractors

Cat.No.: R22030-03

- 20.5mm Spread | 3cm
- Price: \$130.00

## Tying Forceps



### Tying Forceps

Cat.No.: FC11001R-11

- Straight(SS)
- Tips: 0.3x0.5mm | 11cm
- Price: \$82.00



### Tying Forceps

Cat.No.: FC11002R-08

- Straight(SS)
- Tips: 0.3x0.5mm | 8.5cm
- Price: \$70.00



### Tying Forceps

Cat.No.: FC11001T-11

- Straight(TA)
- Tips: 0.3x0.5mm | 11cm
- Price: \$89.00



### Tying Forceps

Cat.No.: FC11002T-08

- Straight(TA)
- Tips: 0.3x0.5mm | 8.5cm
- Price: \$77.00

## Corneal Forceps



### Corneal Forceps

Cat.No.: FC12001T-11

- Straight(TA)
- Tips: 28mmx0.3mm | 10.7cm
- Price: \$122.00



### Corneal Forceps

Cat.No.: FC12002R-08

- Straight(SS)
- Tips: 28mmx0.3mm | 8.5cm
- Price: \$81.00



### Corneal Forceps

Cat.No.: FC12001R-11

- Straight(SS)
- Tips: 28mmx0.3mm | 10.7cm
- Price: \$108.00



### Corneal Forceps

Cat.No.: FC12002T-08

- Curved(TA)
- Tips: 28mmx0.3mm | 8.5cm
- Price: \$96.00

## Tissue Forceps



### Tissue Forceps

Cat.No.:FC13002R-08

- Straight(SS)
- Tips:0.3x0.5mm | 8.5cm
- Price:\$104.00



### Tissue Forceps

Cat.No.:FC13001T-11

- Straight(TA)
- Tips:0.3x0.5mm | 10.7cm
- Price:\$135.00



### Tissue Forceps

Cat.No.:FC13001R-11

- Straight(SS)
- Tips:0.3x0.5mm | 10.7cm
- Price:\$122.00



### Tissue Forceps

Cat.No.:FC13002T-08

- Straight(TA)
- Tips:0.3x0.5mm | 8.5cm
- Price:\$113.00

## Capsulorhexis Forceps & Cataract Scissors



### Capsulorhexis Forceps

Cat.No.:FC14001T-11

- Curved(TA)
- Tips:0.6x0.6mm | 11cm
- Price:\$162.00



### Capsulorhexis Forceps

Cat.No.:FC14001R-11

- Curved(SS)
- Tips:0.6x0.6mm | 11cm
- Price:\$162.00



### Cataract Scissors

Cat.No.:SC12001R-13

- Straight(SS)
- Tips:19x1mm | 13cm
- Price:\$139.00



### Cataract Scissors

Cat.No.:SC12002R-13

- Curved(SS)
- Tips:19x1mm | 13cm
- Price:\$139.00

## Corneal/Conjunctival Scissors



Cat.No.: SC11001R-12

- Curved(SS)
- Sharp&Sharp
- Cutting Edge (Length): 15mm
- Width: 3mm | 12cm
- Price: \$138.00

### Corneal/Conjunctival Scissors

Cat.No.: SC11002R-11

- Curved(SS)
- Sharp&Sharp
- Cutting Edge (Length): 11mm
- Width: 3mm | 11cm
- Price: \$122.00



### Corneal/Conjunctival Scissors

Cat.No.: SC11001T-12

- Curved(TA)
- Sharp&Sharp
- Cutting Edge (Length): 15mm
- Width: 3mm | 12cm
- Price: \$156.00



### Corneal/Conjunctival Scissors

Cat.No.: SC11002T-11

- Curved(TA)
- Sharp&Sharp
- Cutting Edge (Length): 11mm
- Width: 3mm | 11cm
- Price: \$137.00



### Ear Tag applicator

Cat.No.: C23001-14

- Straight
- Punch Dia. : 2.55mm | 14.5cm
- Width: 3mm | 11cm
- Price: \$41.00

### Ear Tag

Cat.No.: 23001-5

- Black numbers A1-A500
- Price: \$112.00

## Glass Bead Sterilizer

The sterilizer provides a safe, convenient and affective method for sterilizing your surgical instruments before and during surgery, ensuring the tools free of pathogens and microbial contaminants. The highest temperature can reach 300 °C. The surgical instrument could sterilized completed Within 5 seconds.

### Ordering Information:

| Model  | Products Description                    | Price    |
|--------|---|----------|
| RS1501 | 150g Glass Bead Sterilizer-4x8cm(110V)  | \$429.00 |
| RS1502 | 150g Glass Bead Sterilizer-4x8cm(220V)  | \$429.00 |
| RS3001 | 300g Glass Bead Sterilizer-4x14cm(110V) | \$514.00 |
| RS3002 | 300g Glass Bead Sterilizer-4x14cm(220V) | \$514.00 |
| RS1500 | Glass Beads for Sterilizer, 1500g       | \$40.00  |



## Sterilization Trays

### Sterilization Trays

Cat.No.: C12006-01

- Plastic
- Width: 270x160x22mm
- Price: \$70.00

### Sterilization Trays

Cat.No.: C12006-02

- Plastic
- Width: 215x135x22mm
- Price: \$58.00



## Brain & Heart Matrices

RWD Brain and Heart matrices is made of stainless steel, electroless plating, through precision machining and highly polished, strong durability, which can be heated, chilled, autoclaved scrubbed and use in variety of harsh conditions. The surface will not be damaged by thermal pressure, and will not be damaged by mechanical impact, keeping the product in use for a long time.

Brain matrices have coronal sections and sagittal sections, each section can be as 0.5 mm or 1 mm. Widely use in preparing pathological tissue section and cut out TTC dyeing organization.






### Order Information:

| Model | Product Description   | A (mm) | B (mm) | C (mm) | D (mm) | E (mm) | Depth (mm) | Diagram |
|-------|---|--------|--------|--------|--------|--------|------------|---------|
| 68707 | Brain Matrix, Mouse 40-75g, Coronal, 1mm, Stainless steel     | 3.18   | 11.1   | 8.73   | 19.1   | 12.5   | 7.40       |         |
| 68708 | Brain Matrix, Mouse 40-75g, Sagittal, 1mm, Stainless steel    | 3.18   | 11.1   | 8.73   | 19.1   | 12.5   | 7.40       |         |
| 68709 | Brain Matrix, Rat, 175-300g, Coronal, 1mm, Stainless steel    | 4.76   | 17.9   | 14.7   | 36.6   | 27.35  | 8.94       |         |
| 68710 | Brain Matrix, Rat, 175-300g, Sagittal, 1mm, Stainless steel   | 4.76   | 17.9   | 14.7   | 36.6   | 27.35  | 8.94       |         |
| 68711 | Brain Matrix, Rat, 300-600g, Coronal, 1mm, Stainless steel    | 4.76   | 19.8   | 14.7   | 36.6   | 25.7   | 10.91      |         |
| 68712 | Brain Matrix, Rat, 300-600g, Sagittal, 1mm, Stainless steel   | 4.76   | 19.8   | 14.7   | 36.6   | 25.7   | 10.91      |         |
| 68713 | Brain Matrix, Mouse, 40-75g, Coronal, 0.5mm, Stainless steel  | 3.18   | 11.1   | 8.73   | 19.1   | 12.5   | 7.40       |         |
| 68714 | Brain Matrix, Mouse, 40-75g, Sagittal, 0.5mm, Stainless steel | 3.18   | 11.1   | 8.73   | 19.1   | 12.5   | 7.40       |         |
| 68715 | Brain Matrix, Rat, 175-300g, Coronal, 0.5mm, Stainless steel  | 4.76   | 17.9   | 14.7   | 36.6   | 27.35  | 8.94       |         |
| 68716 | Brain Matrix, Rat, 175-300g, Sagittal, 0.5mm, Stainless steel | 4.76   | 17.9   | 14.7   | 36.6   | 27.35  | 8.94       |         |
| 68717 | Brain Matrix, Rat, 300-600g, Coronal, 0.5mm, Stainless steel  | 4.76   | 19.8   | 14.7   | 36.6   | 25.7   | 10.91      |         |
| 68718 | Brain Matrix, Rat, 300-600g, Sagittal, 0.5mm, Stainless steel | 4.76   | 19.8   | 14.7   | 36.6   | 25.7   | 10.91      |         |
| 68719 | Heart Matrix, Mouse 20-75g, Coronal, 1mm, Stainless steel     | 8.3    | 12.7   | —      | —      | —      | 4.8        |         |
| 68720 | Heart Matrix, Rat 175-400g, Sagittal, 1mm, Stainless steel    | 12.7   | 19.9   | —      | —      | —      | 9.6        |         |
| 68721 | Heart Matrix, Mouse 20-75g, Coronal, 0.5mm, Stainless steel   | 8.3    | 12.7   | —      | —      | —      | 4.8        |         |
| 68722 | Heart Matrix, Rat 175-400g, Sagittal, 0.5mm, Stainless steel  | 12.7   | 19.9   | —      | —      | —      | 9.6        |         |

## MCAO Sutures

RWD MCAO suture utilizes high-quality nylon thread to induce ischemic stroke model. The within-group variation of stroke models can be maximally reduced by the super smooth surface and evenness of RWD suture. And without rupturing the vesicle, RWD suture makes its insertion into the animal brain much easier. The stability is also greatly increased for the cerebral ischemia modeling and more consistent data. RWD sutures are available in nylon monofilament material, diameters of silicone-rubber coated part range from 0.19mm to 0.44mm; Length of suture from 2-3cm.



### ► Features

- It adopts high-quality nylon thread with very good flexibility, not easy to puncture blood vessels, which increase the success of stroke model development.
  - Ensure high uniformity and stability of the MCAO model.
  - Silicone-rubber coated in tip, greatly reduce the risk of puncturing blood vessels.
  - The line is marked with black dots, increased the observability during plugging into the blood vessel.
  - The new manufacturing process makes the sutures lighter and easier to use. The success rate of model making is higher and it is more suitable for long-term experiment.
- L1 in mouse suture: 3~4mm
  - L1 in rat suture: 5~6mm
  - L1 could be customized: 2~8mm

### Order Information:

| Model          | Diameter of Head end (mm) | Silicone head length (mm) | Diameter of suture body (mm) | Total Length (mm) | Animal's weight (g) | Packing (pcs) |
|----------------|---------------------------|---------------------------|------------------------------|-------------------|---------------------|---------------|
| MSMC19B007PK50 | 0.17~0.19                 | 3~4                       | 0.07                         | 30                | 15~20g              | 50            |
| MSMC19B008PK50 | 0.17~0.19                 | 3~4                       | 0.08                         | 30                | 15~20g              | 50            |
| MSMC19B100PK50 | 0.17~0.19                 | 3~4                       | 0.1                          | 30                | 15~20g              | 50            |
| MSMC21B007PK50 | 0.20~0.21                 | 3~4                       | 0.07                         | 30                | 21~25g              | 50            |
| MSMC21B008PK50 | 0.20~0.21                 | 3~4                       | 0.08                         | 30                | 21~25g              | 50            |
| MSMC21B100PK50 | 0.20~0.21                 | 3~4                       | 0.1                          | 30                | 21~25g              | 50            |
| MSMC21B120PK50 | 0.20~0.21                 | 3~4                       | 0.12                         | 30                | 21~25g              | 50            |
| MSMC23B007PK50 | 0.22~0.23                 | 3~4                       | 0.07                         | 30                | 26~30g              | 50            |
| MSMC23B008PK50 | 0.22~0.23                 | 3~4                       | 0.08                         | 30                | 26~30g              | 50            |
| MSMC23B100PK50 | 0.22~0.23                 | 3~4                       | 0.1                          | 30                | 26~30g              | 50            |
| MSMC23B120PK50 | 0.22~0.23                 | 3~4                       | 0.12                         | 30                | 26~30g              | 50            |
| MSMC25B120PK50 | 0.24~0.25                 | 3~4                       | 0.12                         | 30                | 31~35g              | 50            |
| MSMC25B150PK50 | 0.24~0.25                 | 3~4                       | 0.15                         | 30                | 31~35g              | 50            |
| MSMC26B150PK50 | 0.25~0.26                 | 3~4                       | 0.15                         | 30                | >35g                | 50            |
| MSRC32B200PK50 | 0.31~0.32                 | 5~6                       | 0.2                          | 50                | <200g               | 50            |
| MSRC35B200PK50 | 0.33~0.35                 | 5~6                       | 0.2                          | 50                | 200~250g            | 50            |
| MSRC37B200PK50 | 0.36~0.37                 | 5~6                       | 0.2                          | 50                | 251~280g            | 50            |
| MSRC40B200PK50 | 0.38~0.40                 | 5~6                       | 0.2                          | 50                | 281~330g            | 50            |
| MSRC42B200PK50 | 0.41~0.42                 | 5~6                       | 0.2                          | 50                | 331~400g            | 50            |
| MSRC45B200PK50 | 0.43~0.45                 | 5~6                       | 0.2                          | 50                | >400g               | 50            |

**000 «Диаэм»**

**Москва**

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

**www.dia-m.ru**

**С.-Петербург**  
+7 (812) 372-6040  
spb@dia-m.ru

**Новосибирск**  
+7 (383) 328-0048  
nsk@dia-m.ru

**Воронеж**  
+7 (473) 232-4412  
vrn@dia-m.ru

**Йошкар-Ола**  
+7 (927) 880-3676  
nba@dia-m.ru

**Красноярск**  
+7 (923) 303-0152  
krsk@dia-m.ru

**Казань**  
+7 (843) 210-2080  
kazan@dia-m.ru

**Ростов-на-Дону**  
+7 (863) 303-5500  
rnd@dia-m.ru

**Екатеринбург**  
+7 (912) 658-7606  
ekb@dia-m.ru

**Кемерово**  
+7 (923) 158-6753  
kemerovo@dia-m.ru

**Армения**  
+7 (094) 01-0173  
armenia@dia-m.ru

