

Pipette Tips Compatibility Chart

Pipette			Pipette Tips														
			10µL			20µL		100µL	200µL			1000µL		Alternative RAININ LTS pipette tips			
Brand	Model	Range (µl)	AI10H/ AI10FH AI10HH	AI10LH/ AI10LFH	AI10DR1	AI20H	AI20FH	AI100FH	AI200H/ AI200FH AI200HH	AI200LH/ AI200LFH	AI200DR1	AI1000H/ AI1000FH	AI1250H/ AI1250FH AI1250HH	AI20HRN/ AI20FHRN	AI200HRN/ AI200FHRN	AI1000HRN/ AI1000FHRN	
Eppendorf		0.1-2.5	✓	✓	✓												
		0.5-10	✓	✓	✓	✓											
		2-20				✓					✓						
		10-100					✓	✓	✓	✓	✓						
		20-200					✓	✓	✓	✓	✓						
		100-1000											✓	✓			
Thermo	F2	1-10	✓	✓		✓											
		2-20							✓	✓	✓						
		20-200					✓	✓	✓	✓	✓						
	F3	1-10	✓	✓		✓											
		5-50					✓	✓	✓	✓	✓						
		10-100					✓	✓	✓	✓	✓						
		20-200					✓	✓	✓	✓	✓						
Gilson		0.5-10	✓	✓	✓												
		2-20			✓	✓											
		20-200					✓	✓	✓	✓	✓						
		100-1200											✓	✓			
RAININ	LTS System	10												✓			
		100													✓		
		200													✓		
		1000														✓	
	Universal RAININ Pipette	2	✓	✓	✓	✓											
		10	✓	✓	✓	✓											
		20				✓											
		100					✓	✓	✓	✓	✓						
		200					✓	✓	✓	✓	✓						
1000											✓	✓					
DLAB		0.1-2.5	✓	✓	✓	✓											
		0.5-10	✓	✓	✓	✓											
		2-20	✓	✓	✓												
		5-50					✓	✓	✓	✓	✓						
		10-100					✓	✓	✓	✓	✓						
		20-200					✓	✓	✓	✓	✓						
		50-200					✓	✓	✓	✓	✓						
		100-1000											✓	✓			