

Pipette accessories

Pipette stands

We offer a variety of pipette stands to store multiple or single pipettes. To ensure safety and convenience, pipettes must be protected and stored safely, preferably in an upright position.

The Finnpiquette F-stand accommodates six pipettes. The sturdy crossbar structure with a slot for each pipette ensures that the pipettes will sit snugly in the stand. The Finnpiquette F-stand is suitable for both single channel and multichannel pipettes.

The new Finnpiquette multichannel stand is a multi-position stand for a multichannel pipette. The pipette

can be stored in an upright position or in a horizontal position. The stand has a small footprint and can accommodate all multichannel Thermo Scientific Finnpiquettes. The Finnpiquette Novus stand is designed for convenient storage of an electronic Finnpiquette Novus pipette. Both single channel and multichannel Finnpiquette Novus models can be used with ease.

Autoclavable or disposable reagent reservoirs for a variety of applications

The autoclavable reagent reservoirs are suitable for many purposes. They come in three different models. The V-shaped bottom of the 60

ml reservoir ensures the recovery of liquids, whereas the divided reservoirs reduce reagent costs. These non-autoclavable disposable reagent reservoirs are made of polystyrene and have three designs ranging from 25 to 100 ml. They assure safety, not only due to the extra-wide base that offers stability but also due to the pour-off spouts that are located in all four corners. All the disposable reservoirs have a V-shaped bottom that maximizes recovery of liquids. In addition, the divided reservoir saves time and reagent costs by enabling handling of two liquids at the same time, also in smaller volumes.



Pipette stands and hangers

| Cat. no. | Description |
|----------|--|
| 9420400 | NEW! Finnpiquette F-Stand (for 6 pipettes) |
| 9420390 | NEW! Finnpiquette Multichannel Stand (for one pipette) |
| 9420 340 | Finnpiquette Carousel Stand (6 pipettes or 3 Novus) |
| 9420 290 | Pipette stand* (for 6 pipettes) |
| 9420 320 | Pipette mini stand (for 3 pipettes) |
| 2206 040 | Shelf hanger* |
| 9420 360 | Novus stand (for one pipette) |

* Fits all Finnpiquette models with a grippy finger rest.

Reagent basin and reservoirs

| Cat. no. | Description | Quantity |
|----------|--------------------------------|----------|
| 9510 027 | Reagent basin, 60 ml | 5/box |
| 9510 037 | Reservoir, 8-well | 10/box |
| 9510 047 | Reservoir, 1-well | 10/box |
| 95128093 | Reagent Reservoir 25ml | 100/box |
| 95128095 | Reagent Reservoir 25ml divided | 100/box |
| 95128085 | Reagent Reservoir 100ml | 100/box |
| 95128080 | Reagent Reservoir demo package | 15/box |