

# MX-RL-Pro


## LCD Digital Tube Rotator, Horizontal

### Features

- LCD displays speed and time
- Adjustable speed range of 10 to 70 rpm
- Provides gentle but effective rolling and mixing
- Can be used in cold rooms or incubators
- Wide variety of Disk accessories and Rotisserie accessories available

### Applications

Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



Specifications	MX-RL-Pro
Motor type	Brushless DC motor
Speed range	10–70rpm
Speed display	LCD
Timer display	LCD
Timer	1–1199min
Operation mode	Continuous / timer
Voltage,Frequency	100–240V, 50/60Hz
Power	40W
Dimension [D × W × H ]	220 × 510 × 260mm
Permissible ambient temperature and humidity	5–40°C,80%RH
Protection class	IP21

**Disk Accessories Type-1 Cat. No.18900160**

Used for MX-RD-Pro  
for 1.5mLx60 centrifuge tubes holder

**Disk Accessories Type-2 Cat. No.18900161**

Used for MX-RD-Pro  
for 15mLx16 centrifuge tubes holder

**Disk Accessories Type-3 Cat. No.18900162**

Used for MX-RD-Pro  
for 50mLx8 centrifuge tubes holder

**Disk support rods Cat. No.18900140**

4pcs in a group  
use with disk accessories

**Rotisserie Accessory Type-3 Cat. No.18900144**

Used for MX-RL-Pro, for 50mLx24  
centrifuge tubes held horizontally

**Rotisserie Accessory Type-4 Cat. No.18900145**

Used for MX-RL-Pro, for 1.5mLx32  
centrifuge tubes held vertically

**Rotisserie Accessory Type-1 Cat. No.18900142**

Used for MX-RL-Pro, for 1.5mLx48  
centrifuge tubes held horizontally

**Rotisserie Accessory Type-5 Cat. No.18900146**

Used for MX-RL-Pro, for 15mLx16  
centrifuge tubes held vertically

**ООО «Диаэм»**

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

**www.dia-m.ru**

**С.-Петербург**  
+7 (812) 372-6040  
spb@dia-m.ru

**Новосибирск**  
+7 (383) 328-0048  
nsk@dia-m.ru

**Воронеж**  
+7 (473) 232-4412  
vrn@dia-m.ru

**Йошкар-Ола**  
+7 (927) 880-3676  
nba@dia-m.ru

**Красноярск**  
+7 (923) 303-0152  
krsk@dia-m.ru

**Казань**  
+7 (843) 210-2080  
kazan@dia-m.ru

**Ростов-на-Дону**  
+7 (863) 303-5500  
rnd@dia-m.ru

**Екатеринбург**  
+7 (912) 658-7606  
ekb@dia-m.ru

**Кемерово**  
+7 (923) 158-6753  
kemeroovo@dia-m.ru

**Армения**  
+7 (094) 01-0173  
armenia@dia-m.ru

