



IKA

designed for scientists

EN



BEAUTIFUL STIRRING | IKA myPlate

IKA myPlate

/// The magnetic stirrer with an individual set-up plate

IKA myPlate combines the excellent technical features of the popular and proven IKA magnetic stirrers colour squid and big squid with the option of creative individualization.



Personal interchangeable frame

Design the set-up plate of the IKA myPlate by yourself. You receive the device with a decorative image template, that can be exchanged easily through the convenient interchangeable frame. Some people have family photos on their desk – you convert one of the best magnetic stirrers from your day-to-day lab work into a canvas for your favorite pictures.

Proven high functionality

Apart from stirring your samples with IKA myPlate on a creative surface, you get all good features of this strong stirrer, just as you are used to from IKA.

- > High chemical resistance
- > Electronically controlled motor
- > Digital speed indicator
- > Speed range between 0 and 2,500 rpm

TECHNICAL DATA IKA MYPLATE

Stirring quantity max. per stirring position (H ₂ O)	1.5 l
Direction of rotation	Right
Speed	min. 0 rpm max. 2,000 rpm
Set-up plate material	Glass
Set-up plate dimensions	Ø 135 mm
Protection class according to DIN EN 60529	IP 54
Voltage	100 – 240 V
Weight	0.5 kg
Dimensions (W × H × D)	150 × 30 × 160 mm
Power input	4 W Standby: 0.1 W

CREATE THE SURFACE DESIGN ON IKA.COM:

Pre-punched templates are included in the scope of delivery. To obtain perfect prints of your own pictures, upload them on the IKA website and create the exact matching print to individualize your magnetic stirrer.

IKA myPlate

Ident. No. 0025005156

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

