





Your needs, our innovations

Wash bottles



Wash bottles

Made a combination of high and low density polyethylene. Translucent bottle, with HDPE cap (blue), delivery tube and seal cap (natural).
Bottle caps are ribbed for a better handling.
High resistance to all liquids, acids and bases.
Include a large neck for safe and easy filling.



code	capac. ml	neck mm	base mm	height mm	case quantity	case weight	case volume
44330000	250	30	65	155	90	5.30	0.078
44300000	500	30	75	185	60	4.85	0.080
44310000	1,000	30	90	222	45	5.02	0.110

Wide mouth wash bottles

Translucent low density polyethylene bottles, with screw cap and delivery tube available in 3 different colours (natural, blue, yellow and red) for easy identification. Feature a one piece cap and delivery tube assembly, thus eliminating air or liquid absorption through the cap. Include a wide neck for safe and easy filling.
Ribbed caps for a better handling while opening and closing.
Moulded-in graduations.
250 ml model is graduated in 25 ml increments, 500 ml model, in 100 ml increments.



code	capac. ml	cap colour	int. neck Ø mm	base Ø mm	height mm	case quantity	case weight	case volume
191637	250	natural	38	60	140	10	0.64	0.018
191637.04	250	blue	38	60	140	10	0.64	0.018
191637.06	250	yellow	38	60	140	10	0.64	0.018
191637.10	250	red	38	60	140	10	0.64	0.018
191638	500	natural	38	75	165	10	1.01	0.019
191638.04	500	blue	38	75	165	10	1.01	0.019
191638.06	500	yellow	38	75	165	10	1.01	0.019
191638.10	500	red	38	75	165	10	1.01	0.019

Integral wash bottles

Made of translucent polyethylene. Feature a one piece cap and delivery tube assembly, which eliminates air or liquid absorption through the cap and ensures a positive leaktight seal. To increase liquid flow, cut the tip of the delivery tube; close the tip by means of the cap attached to the delivery tube.
Caps are ribbed for a better handling while opening and closing.
Tagging points on neck and cap for a tamper proof closure. Bottles include a fill in line.



code	capac. ml	height mm	neck d/n std.	case quantity	case weight	case volume
191633	250	143	GL32	10	0.68	0.008
191634	500	181	GL32	10	0.84	0.015

000 «Диаэм» Москва ул. Магаданская, д. 7, к. 3 • тел./факс: (495) 745-0508 • sales@dia-m.ru www.dia-m.ru

С.-Петербург +7 (812) 372-8040 spb@dia-m.ru	Новосибирск +7 (383) 328-8048 nsk@dia-m.ru	Воронеж +7 (473) 232-4412 vno@dia-m.ru	Йошкар-Ола +7 (927) 880-3676 ok@dia-m.ru	Красноярск +7 (923) 303-0152 krk@dia-m.ru
Казань +7 (843) 210-2080 kazan@dia-m.ru	Ростов-на-Дону +7 (863) 303-5500 rnd@dia-m.ru	Екатеринбург +7 (912) 658-7606 ekb@dia-m.ru	Кемерово +7 (923) 158-6753 kemerovo@dia-m.ru	Армения +7 (094) 01-0173 armenia@dia-m.ru

