


## BagLight® PolySilk®

Blender bag without filter, transparent

**Ref. 131 025** (100 mL) – **Ref. 132 025** (400 mL) – **Ref. 132 050** (400 mL) – **Ref. 132 200** (2000 mL)  
**Ref. 133 025** (3500 mL)



Standards:



Manufactured under:



BagLight® PolySilk® is a blender bag without filter. It is ideal for sample preparation.

### TECHNICAL SPECIFICATIONS

- Non-filter bag
- **PolySilk®**, polyolefin complex
- Flexible and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Shelf-life: 10 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 100/400/2000/3500 mL
- In compliance with ISO 7218, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)
- Approved for food contact - Regulation (EC) N° 1935/2004
- Designed and made in France

### BagLight® PolySilk® 100

**Ref. 131 025**

- Max blending volume: 100 mL
- Optimal blending volume: 5 - 50 mL
- Bag dimensions: 110 x 200 mm
- Box dimensions: 31 x 22 x 7 cm, weight: 1.80 kg
- Pack of 25
- Box of 500

### BagLight® PolySilk® 400

**Ref. 132 025**

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 25
- Box of 500

## BagLight® PolySilk® 400

Ref. 132 050

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 50
- Box of 500

## BagLight® PolySilk® 2000

Ref. 132 200

- Max blending volume: 2000 mL
- Optimal blending volume: 400 - 1500 mL
- Bag dimensions: 250 x 380 mm
- Box dimensions: 48 x 29 x 14 cm, weight: 10.2 kg
- Pack of 25
- Box of 500

## BagLight® PolySilk® 3500

Ref. 133 025

- Max blending volume: 3750 mL
- Optimal blending volume: 400 - 3750 mL
- Bag dimensions: 380 x 560 mm
- Box dimensions: 41 x 34 x 12 cm, weight: 7.70 kg
- Pack of 25
- Box of 250

**000 «Диаэм»**

**Москва**

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)

[www.dia-m.ru](http://www.dia-m.ru)

**С.-Петербург**

+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)

**Новосибирск**

+7 (383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)

**Воронеж**

+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)

**Йошкар-Ола**

+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)

**Красноярск**

+7 (923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)

**Казань**

+7 (843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)

**Ростов-на-Дону**

+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)

**Екатеринбург**

+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)

**Кемерово**

+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)

**Армения**

+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)

