

EPA/TOC Vials



Headspace Vials



ALWSCI
Technologies

ДИАМ

www.dia-m.ru

Be Simple,
Be Professional

Standard + Pre-Cleaned + Certified

SKU	Product Description	Packing
C0000061	40mL Clear Glass Sample Vial 27.5x95mm 24-414 Screw Thread.	100pcs/pk.
C0000062	40mL Amber Glass Sample Vial 27.5x95mm 24-414 Screw Thread.	100pcs/pk.
C0000211	24-414 White Open Top PP Screw Cap with 22mm Natural PTFE/Natural Silicone Septa 3mm Thick.	100pcs/pk.
C0001414	24-414 White Open Top PP Screw Cap Bonded with 22mm Natural PTFE/White Silicone Septa 3mm Thick.	100pcs/pk.
C0000342	40mL 27.5x95mm Clear Glass EPA/TOC Vial 24-414 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality) and Dust Cover. Vial + Cap + Septa are pre-assembled together.	72pcs/pk.
C0000343	40mL 27.5x95mm Amber Glass EPA/TOC Vial 24-414 White Open Top PP Screw Cap with 22mm Natural PTFE/White Silicone 3.0mm thick Septa (EPA Quality) and Dust Cover. Vial + Cap + Septa are pre-assembled together. 72pcs/pk.	72pcs/pk.
C0002296	Pre-cleaned and certified (<10 ppb TOC) 40mL 27.5x95mm clear glass vial 24-414 Purple Open Top PP Screw Cap Bonded with 22mm Natural PTFE/White Silicone 3.0mm thick Septa and Dust Cover. Vial + Cap + Septa are pre-assembled together.	72pcs/pk.

To be continued...

SPME + Gas Analysis + GC/MS

SKU	Product Description	Packing
C0000028	10mL Clear Glass Screw Headspace Vial 22.5x46mm.	100pcs/pk.
C0000030	20mL Clear Glass Screw Headspace Vial 22.5x75mm.	100pcs/pk.
C0000842	18mm Silver Color Open Top Metal Cap with Blue PTFE/White Silicone Septa 1.3mm Thick.	100pcs/pk.
C0000208	18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm Blue PTFE/White Silicone Septa 1.5mm Thick.	100pcs/pk.
C0000974	18mm Silver Color Open Top Metal Cap (8mm hole) with 17.5mm White PTFE/Blue Transparent Silicone Septa 1.5mm Thick (UltraClean).	100pcs/pk.
C0000033	10mL Clear Glass 22.5x46mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial.	100pcs/pk.
C0000036	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Flat Bottom. Crimp Headspace Vial.	100pcs/pk.
C0000039	20mL Clear Glass 22.5x75mm. 20mm Beveled Edge. Rounded-Flat Bottom. Crimp Headspace Vial. Short Neck.	100pcs/pk.
C0000196	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/White Silicone Septa 3mm Thick.	100pcs/pk.
C0000973	20mm Open Top Aluminum Crimp Cap (10mm hole) with 20mm Natural PTFE/Blue Transparent Silicone Septa 3mm Thick (UltraClean).	100pcs/pk.

To be continued...



000 «Диазм»

Москва
ул. Магаданская, д. 7, к. 3 • тел./факс: (495) 745-0508 • sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vnn@dia-m.ru

Иошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru



Autosampler Vials



Screw



Crimp



Snap

Storage Vials



Barcoded



1-1000ML



Customized

Lab Filtration



Syringe
Filter



Syringe



Membrane
Filter

SKU	Product Description	Packing
C0000002	2mL Clear Glass 12x32mm Flat Base 8-425 Screw Thread Vial with Label.	100pcs/pk.
C0000004	2mL Amber Glass 12x32mm Flat Base 8-425 Screw Thread Vial with Label.	100pcs/pk.
C0000133	Black 8-425 Open Top Screw Cap with 8mm Red PTFE/White Silicone Septa 1.5mm Thick.	100pcs/pk.
C0000410	2mL 9-425 Wide Opening PP Screw Vial with Graduations Transparent .	100pcs/pk.
C0000008	2mL Clear Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label.	100pcs/pk.
C0000011	2mL Amber Glass 12x32mm Flat Base 9-425 Screw Thread Vial with Label.	100pcs/pk.
C0000143	Blue 9-425 Open Top Ribbed Screw Cap with 9mm Red PTFE/White Silicone Septa 1mm Thick.	100pcs/pk.
C0000148	Blue 9-425 Open Top Ribbed Screw Cap with Red PTFE/White Silicone Septa 1mm Thick PRE-SLIT.	100pcs/pk.
C0000144	Blue 9-425 Open Top Ribbed Screw Cap Bonded with Red PTFE/White Silicone Septa 1mm Thick.	100pcs/pk.
C0000015	2mL Clear Glass Flat Base 11mm Crimp Vial Wide Opening with Label.	100pcs/pk.
C0000018	2mL Amber Glass Flat Base 11mm Crimp Vial Wide Opening with Label.	100pcs/pk.
C0000169	SILVER 11mm Open Top Crimp Cap Transparent TEF/Natural Rubber red-orange Septa 1mm Thick.	100pcs/pk.
C0000021	2mL Amber Glass Flat Base 11mm Snap Vial Wide Opening with Label.	100pcs/pk.
C0000023	2mL Amber Glass Flat Base 9-425 Screw Thread Vial. Pre-assembled together.	100pcs/pk.
C0000170	CLEAR 11mm Open Top Snap Cap with Red PTFE/White Silicone Septa 1mm Thick.	100pcs/pk.

SKU	Product Description	Packing
C0000024	4mL Clear Glass 15x45mm Flat Base 13-425 Screw Thread Vial.	100pcs/pk.
C0000044	4mL Amber Glass Sample Vial 15x45mm 13-425 Screw Thread.	100pcs/pk.
C0000183	13-425 Black Closed Top PP Cap with Red PTFE/White Silicone Septa 1mm Thick.	100pcs/pk.
C0000051	8mL Clear Glass Sample Vial 17x60mm 15-425 Screw Thread.	100pcs/pk.
C0000052	8mL Amber Glass Sample Vial 17x60mm 15-425 Screw Thread.	100pcs/pk.
C0001370	12mL Clear Glass Sample Vial 19x65mm 15-425 Screw Thread.	100pcs/pk.
C0001371	12mL Amber Glass Sample Vial 19x65mm 15-425 Screw Thread.	100pcs/pk.
C0000191	15-425 Black Closed Top PP Cap with Natural PTFE/White Silicone Septa 2mm Thick.	100pcs/pk.
C0001410	15mL Clear Glass Sample Vial 21x70mm 18-400 Screw Thread.	100pcs/pk.
C0001638	15mL Amber Glass Sample Vial 21x70mm 18-400 Screw Thread.	100pcs/pk.
C0000220	18-400 Black Closed Top PP Cap with Natural PTFE/White Silicone Septa 1.5mm Thick.	100pcs/pk.
C0000057	20mL Clear Glass Sample Vial 27.5x57mm 24-414 Screw Thread.	100pcs/pk.
C0000058	20mL Amber Glass Sample Vial 27.5x57mm 24-414 Screw Thread.	100pcs/pk.
C0000432	22mL Clear Glass Sample Vial 23x85mm. 20-400 Screw Thread.	100pcs/pk.
C0000227	20-400 Black Closed Top PP Cap with Natural PTFE/White Silicone Septa 2mm Thick.	100pcs/pk.

SKU	Product Description	Packing
C0000639	13mm Nylon Syringe Filter 0.45um with Outer Ring and Printing.	100pcs/pk.
C0000640	25mm Nylon Syringe Filter 0.22um with Outer Ring and Printing.	100pcs/pk.
C0000641	25mm Nylon Syringe Filter 0.45um with Outer Ring and Printing.	100pcs/pk.
C0000646	13mm PTFE Hydrophobic Syringe Filter 0.22um with Outer Ring and Printing.	100pcs/pk.
C0000648	25mm PTFE Hydrophobic Syringe Filter 0.22um with Outer Ring and Printing.	100pcs/pk.
C0000649	25mm PTFE Hydrophobic Syringe Filter 0.45um with Outer Ring and Printing.	100pcs/pk.
C0000644	25mm PTFE Hydrophilic Syringe Filter 0.22um with Outer Ring and Printing.	100pcs/pk.
C0000656	25mm PVDF Hydrophobic Syringe Filter 0.22um with Outer Ring and Printing.	100pcs/pk.
C0000657	25mm PVDF Hydrophobic Syringe Filter 0.45um with Outer Ring and Printing.	100pcs/pk.
C0000653	25mm PVDF Hydrophilic Syringe Filter 0.45um with Outer Ring and Printing.	100pcs/pk.
C0000305	47mm Nylon Membrane 0.45um.	100pcs/pk.
C0000425	47mm PTFE Hydrophilic Membrane 0.45um.	100pcs/pk.
C0000916	47mm PVDF Hydrophobic Membrane 0.22um.	100pcs/pk.
C0001243	47mm Sterile MCE Grid Membrane White Membrane/Black 0.45um.	100pcs/pk.

To be continued...