


HYBRIDIZATION INCUBATOR (NB-202/ NB-202R)

Designed to provide only Hybridization or with Rocking in one incubator. Fast Heat-Up and Precise Temperature Control helps efficient sample mixture, Hybridization, incubation.



Accessories



Bottle rack



Bottles



Rocker platform

Features

- **Safety hybridization procedures**
Stable and safe hybridization provides very efficient results.
- **Microprocessor control**
Microprocessor control allows the precise control of temperature.
- **Stabilized rotation**
Keep stable working and low noise by stable rotation.
- **Outer glass door**
Outer glass door allows the observation without opening the door.
- **Safety switch**
Safety switch is equipped to prevent overheating.
- **Easy-to-change bottle rack**
Bottle rack can be replaced very easily to modify bottle holders.
- **Simultaneous use of rotator and rocker (NB-202R)**
Users can rotate and rock the samples simultaneously.
Rocker is installable in incubator(cat.no NB-202R) like the photo.



Specification



Items	Unit	NB-202	NB-202R
Temperature			
range	°C	Ambient +5°C to 80°C	Ambient +5°C to 80°C
accuracy	°C	±0.3°C at +37°C	±0.3°C at +37°C
safety		Over temperature protector S/W	Over temperature protector S/W
control		Microprocessor digital PID	Microprocessor digital PID
Display		LED Display	LED Display
Speed range	rpm	Variable 2 to 35rpm	Variable 2 to 35rpm
Rack capacity	mm, ea	40(Ø)x200(D)mmx8ea 40(Ø)x120(D)mmx8ea	40(Ø)x200(D)mmx8ea 40(Ø)x120(D)mmx8ea
Rocker platform	mm		300x250(mm)
Chamber dimension	mm	320(W)x350(D)x370(H)mm	320(W)x350(D)x370(H)mm
Overall dimension	mm	415(W)x465(D)x540(H)mm	415(W)x465(D)x540(H)mm
Power	V/ Hz	110/220V, 50/60Hz, 400W	110/220V, 50/60Hz, 400W

ООО «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

