

BagOpen®

Bag opener/holder

Ref 211 010 (80/100 mL) - **Ref 211 040** (400 mL) - **Ref 211 350** (3500 mL)



Manufactured under:



BagOpen® easily opens the sample bag without any contact with the sample.

TECHNICAL SPECIFICATIONS

- Stable and robust base
- Polycarbonate sides
- Resistant up to 65°C
- Available for 80/100/400/3500 mL bags
- Designed and made in France

BagOpen® 80/100

Ref: 211 010

- Available for 80 and 100 mL bags
- Dimensions (w x d x h): 9 x 4.6 x 17 cm
- Box (w x d x h): 23 x 12 x 5 cm, gross weight: 200 g

BagOpen® 400

Ref: 211 040

- Available for 400 and 2000 mL bags
- Dimensions (w x d x h): 15 x 7.5 x 25 cm
- Box (w x d x h): 26 x 16 x 6.2 cm, gross weight: 540 g

BagOpen® 3500

Ref: 211 350

- Available for 3500 mL bags
- Dimensions (w x d x h): 32 x 14.4 x 43.6 cm
- Box (w x d x h): 45 x 20 x 5 cm, gross weight: 2 kg

CAUTION

- Do not put in an autoclave

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

