

## BagClip®

Closing clip

**Ref 231 010** (80/100 mL) - **Ref 231 040** (400 mL) - **Ref 231 350** (3500 mL)



Manufactured under:



BagClip® allows a temporary closing in 3 points of any type of bag.

### TECHNICAL SPECIFICATIONS

- Gas and liquid leakproof
- No contact with the sample
- Available for 80/100/400/2000/3500 mL blender bags
- Designed and made in France

### BagClip® 80/100

**Ref: 231 010**

- Available for 80 and 100 mL bags
- Length (white clip): 130 mm
- Length (green tube): 160 mm
- Box (w x d x h): 21.5 x 15 x 5.5 cm, gross weight: 550 g
- Quantity per pack: 50

### BagClip® 400

**Ref: 231 040**

- Available for 400 mL bags
- Length (white clip): 195 mm
- Length (green tube): 225 mm
- Box (w x d x h): 27 x 19 x 6.5 cm, gross weight: 750 g
- Quantity per pack: 50

### BagClip® 3500

**Ref: 231 350**

- Available for 2000 and 3500 mL bags
- Length (white clip): 395 mm
- Length (green tube): 450 mm
- Box (w x d x h): 51 x 10 x 12 cm, gross weight: 1.6 kg
- Quantity per pack: 50

# interscience

## CAUTION

- Do not put in an autoclave
- Do not put in a washing machine
- Can be disinfected with ethanol 70 %, hydrogen peroxyde 6 % or bleach 10 % max

**000 «Диаэм»**

**Москва**

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)

**[www.dia-m.ru](http://www.dia-m.ru)**

**С.-Петербург**  
+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)

**Новосибирск**  
+7 (383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)

**Воронеж**  
+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)

**Йошкар-Ола**  
+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)

**Красноярск**  
+7 (923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)

**Казань**  
+7 (843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)

**Ростов-на-Дону**  
+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)

**Екатеринбург**  
+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)

**Кемерово**  
+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)

**Армения**  
+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)

