


ДИА•М
современная лаборатория

www.dia-m.ru
заказ on-line

DURAN®
LABORATORY
GLASSWARE
CATALOG
Washing Bottle




**DURAN
WHEATON
KIMBLE**

Excellence in your hands


DURAN® Gas Washing Bottle

with fused-in filter disk, with standard ground joint and cap




DURAN® Gas Washing Bottle

Drehsel type head



DURAN® Gas Washing Bottle

head with filter disk, with standard ground joint



Typical application: cleaning ("washing") of gases with solvents.

Cat. No.	Capacity (mL)	d (OD) (mm)	h (mm)	h ₁ (mm)	Neck	Porosity	Hose connection d (OD) (mm)	Pack Unit
257030109	350	60	250	180	29/32	1	10	1

With screw-connection system. The insertion height of the head is adjustable. Individual parts can also be ordered separately.

Typical applications: cleaning ("washing") of gases with solvents.

Cat. No.	Capacity (mL)	DIN Thread (GL)	Porosity	Hose connection d (OD) (mm)	Bowl (mm)	Pack Unit
without filter disk						
247130008	500	45		9		1
with filter disk						
257040101	500	45	1	9	25	1

Cat. No.	Description	Pack Unit
Individual parts		
247130205	Drehsel-head, without filter disc	1
257540109	Drehsel-head, with filter disc (Por. 1)	1
218014401	DURAN® Laboratory bottle 500 mL, with DIN thread GL 45	10
292550603	Plastic hose connection, straight	10
292270508	Screw-caps with aperture, PBT, GL 14	10
292271007	Screw-caps with aperture, PBT, GL 45	10
292282501	VQM rubber ring (26 x 42 x 5 mm)	10

Typical applications: cleaning ("washing") of gases with solvents.

Cat. No.	Description	Capacity (mL)	h (mm)	h ₁ (mm)	Neck	Porosity	Hose connection d (OD) (mm)	Bowl (mm)	Pack Unit
257010107		100	250	180	34/35	1	10	25	1
257020108		250	250	160	45/40	1	10	34	1
Individual parts									
257520107	Gas washing bottle, head with fritted disc (100 mL)					1			1
257520108	Gas washing bottle, head with fritted disc (250 mL)					1			1

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

