

## 5-1. Plant Culture

### Incu Tissue

Incu Tissues are ideal for plant tissue cultures.

As the sprout grows, another Incu Tissue can be connected to the top as an optional frame.

- Embossing on the closing edge of cover and jar
- Optional frame (Cat. No. 310074) is only available for 310070 and 310071
- Cat. No. 310070 can be assembled with Insect Breeding Box (Cat. No. 310075, 310076, 310077)



#### Incu Tissue

Type	Cat. No.	Material	Style	External Dimensions w x l x h (mm)	Internal Dimensions w x l x h (mm)	Sterile	Packaging
	310070	PS	Square	72.00 x 72.00 x 100.00	65.40 x 65.40 x 98.20	+	4 / 120
	310071	PC	Square	72.00 x 72.00 x 100.00	65.40 x 65.40 x 98.20	-	4 / 120
	310072	PP	Square	72.00 x 72.00 x 100.00	65.40 x 65.40 x 98.20	-	4 / 120
	310074	PP	Frame	80.80 x 80.80 x 21.00	63.80 x 63.80 x 21.00	-	30 / 60

### Plant Culture Dish

Plant Culture Dishes are used for sprout culture.

Pre-sterilized, flat dishes provide excellent culture results.

- Excellent air circulation design
- Excellent flatness
- Stable stacking
- Designed for reduced contamination
- Autoclavable (Cat. No. 310103, 310101)
- Gridded bottom (Cat. No. 310200)



#### Plant Culture Dish

Type	Cat. No.	Material	External Dimensions d x h (mm)	Internal Dimensions d x h (mm)	Autoclavable	Sterile	Packaging
	310103	PP	100.00 x 50.00	91.40 x 48.50	+	-	5 / 200
	310100	PS	100.00 x 40.00	91.30 x 38.20	-	+	5 / 200
	310101	PP	100.00 x 40.00	91.30 x 38.20	+	-	5 / 200
	310200	PS	100.00 x 20.00	94.60 x 18.60	-	+	10 / 200
	10090	PS	90.00 x 15.00	85.90 x 12.60	-	+	10 / 500
	10091	PS	90.00 x 20.00	86.20 x 17.70	-	+	10 / 200
	10100	PS	100.00 x 15.00	96.40 x 13.75	-	+	10 / 500

# Plant Culture



## Plant Culture Bottle I

Plant Culture Bottles I are used for sprout or sapling culture.

- Stackable feature for space saving
- Designed for reduced contamination
- Screw cap



01\_Cell Culture

02\_Molecular Analysis

03\_Microbiology


04\_Handling & Storage

05\_Plant & Insect Culture

06\_Clinical Labware

07\_Appendix

### Plant Culture Bottle I

Type	Cat. No.	Material (Bottle / Cap)	Dimensions w x l x h (mm)	Total Vol. (ml)	Sterile	Packaging
	310500	PP / PP	100.00 x 100.00 x 110.00	500.00	-	1 / 80

## Plant Culture Bottle II

Plant Culture Bottles II are used for culturing sapling or callus.

- Durable plastic bottle
- Designed for reduced contamination
- Screw cap
- Autoclavable



### Plant Culture Bottle II

Type	Cat. No.	Material (Bottle / Cap)	Dimensions w x l x h (mm)	Total Vol. (ml)	Sterile	Packaging
	310501	PC / PP	115.00 x 115.00 x 140.00	500.00	-	1 / 50

# Plant Culture



## Phytohealth

Phytohealth is useful for the culture of individual plantlet.

- Stackable feature for space saving
- Embossing on the closing edge of cover and the jar
- Phytohealth bodies and caps are configured separately
- Autoclavable



### Phytohealth

Type	Cat. No.	Material	External Dimensions d x h (mm)	Internal Dimensions-Body d x h (mm)	Internal Dimensions-Cap d x h (mm)	Sterile	Packaging
	310120	PP	120.00 x 80.00	103.00 x 78.60	101.00 x 10.20	-	400
	310121	PP	120.00 x 110.00	103.00 x 78.60	101.00 x 32.50	-	400







## AraHarvest

AraHarvest is a convenient system for culturing and harvesting seeds from individual sapling.

- AraHarvest set: tube, base, basket, top tray, bottom tray (Cat. No. 310060)
- Aracon Tube (Cat. No. 310061)
- Aracon Base (Cat. No. 310062)
- Aracon Basket (Cat. No. 310063)
- Aracon Top Tray (Cat. No. 310064), space saving 50 hole tray
- Aracon Bottom Tray (Cat. No. 310065)



### AraHarvest

Type	Cat. No.	Material	Tube	Base	Basket	Top Tray	Bottom Tray	Packaging
	310060	PS / PP	200 ea	200 ea	200 ea	4 ea	4 ea	
	310061	PS	100 ea					100
	310062	PS		100 ea				100
	310063	PP			100 ea			100
	310064	PP				4 ea		4
	310065	PP					4 ea	4

000 «Диаэм»

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург  
+7 (812) 372-6040  
spb@dia-m.ru

Новосибирск  
+7 (383) 328-0048  
nsk@dia-m.ru

Воронеж  
+7 (473) 232-4412  
vrn@dia-m.ru

Йошкар-Ола  
+7 (927) 880-3676  
nba@dia-m.ru

Красноярск  
+7 (923) 303-0152  
krsk@dia-m.ru

Казань  
+7 (843) 210-2080  
kazan@dia-m.ru

Ростов-на-Дону  
+7 (863) 303-5500  
rnd@dia-m.ru

Екатеринбург  
+7 (912) 658-7606  
ekb@dia-m.ru

Кемерово  
+7 (923) 158-6753  
kemerovo@dia-m.ru

Армения  
+7 (094) 01-0173  
armenia@dia-m.ru

