

# JACKETED GLASS REACTOR


S-2L




Double-layer Glass



Glass Discharge Valve



## Product Components



## SIX Product Advantages

+



Frequency Control Motor




Stable Operation



High Vacuum Degree



Constant Speed Stirring



Precise Temperature Control



Efficient Distillation

## Product Parameter

Model	S-1L	S-2L	S-3L	S-5L
Stirring Power(W )	60			
Stirring Speed(rpm)	0~600(Max. 1300)			
Stirring Shaft Diameter(mm)	12			
Reaction Kettle Capacity(L)	1	2	3	5
Power Supply	220V 50/60 Hz			
Motor Torque(g/cm)	2000	2400	4800	
Note	Floor type			
Explosion-proof	1. The stirring power (W): 180, 250, 370			
	2. Power supply: 220V / 50Hz or 110V / 60Hz			
110 V/60Hz can be customized				

## Product Details



01 /

**Constant Pressure Funnel**

Feed rapidly, no residual, Acid and alkali resistance, Temp Rang:-100°C~400°C.

02 /

**Gear Motor**

Smooth running, high torque, automatic booster, no spark, Explosion proof is optional.



03 /

**Kettle Cover**

It is convenient to clean, disassemble, and install. PT100 temperature probe.



04 /

**Double-layer Galss Kettle Body**

High qualit borosilicate 3.3 glassware.



## Support Equipment



Cooling Chiller



Circulating Water/Oil Bath



GDX Heater-chiller(200°C)




DFY Heater-chiller(98°C)

## All Models Display



## Intelligent Controller



- More convenient,intuitive and accurate
- Temperaure&speed digital display

## Factory Display

