

ДИАМ
сервисная лаборатория

www.dia-m.ru
lab-on-line



NB-201M MINI INCUBATOR

FEATURES

- 15.2 liter chamber
- 6.8kg light weight to carry
- Portable use with the carrying handle
- Economical price & compact design
- Available size to use in a workstation
- Forced air circulation by the internal fan
- Excellent temperature uniformity
- Cooling with very low noise
- 2 stainless steel shelves/3 maximum
- Portable by car with the car charger
- CO₂ control available (Optional)
- Access port installation available (Optional)



000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru



NB-201M MINI INCUBATOR



NB-201M	MINI INCUBATOR
Chamber capacity	15.2 Liter
Temperature range	15C~45C at ambient 25C
Temperature accuracy	± 0.25C at 37C
Temperature control	Micom
Display	LED display
Cooling and heating	By peltier (thermoelectric element)
Internal fan	Built in internal fan
Shelf	2 (Stainless steel) / maximum 3
In&out material	ABS resin
Inside dimension	224(W) X 200(D) X 340(H) mm
Outside dimension	292(W) X 333(D) X 433(H) mm
Weight	6.8kg
Power frequency	DC 12V, AC 110V~220V, 50~60Hz
Power consumption	DC-COLD : 46W, HOT : 48W AC-COLD : 63W, HOT : 63W