


# Whitley Modified Atmosphere Systems - Porthole Comparison Chart



Workstation	Instant Access Porthole System	Manual Porthole System/1	Manual Porthole System/2	Workstation
DG250	-	●	-	DG250
A25	●	-	-	A25
A35	○	-	○	A35
A45	▼ ○	-	▼ ○	A45
A55	○	-	○	A55
A85	-	●	-	A85
A95	-	●	-	A95
A135	-	-	●	A135
A155	-	-	●	A155
M35	○	-	○	M35
M45	▼ ○	-	▼ ○	M45
M85	-	●	-	M85
M95	-	●	-	M95
H35	-	-	●	H35
H45	▼	-	●	H45
H85	-	●	-	H85
H95	-	●	-	H95
H135	-	-	●	H135
i2	-	-	●	i2

- Eliminates sleeves and cuffs
- Instant access
- International Patents
- Ideal for more often / shorter duration use

- 10 Sample transfer capacity
- Sleeve and cuff attachment
- Oval, ergonomic shape
- Integrated shelf for sample transfer
- International Patents

- Sleeve and cuff attachment
- Oval, ergonomic shape
- International Patents
- Ideal for less often / longer duration use

Key: ● Fitted as standard ○ Can be specified with order - Not available

000 «Диаэм» Москва ул. Магаданская, д. 7, к. 3 тел./факс: (495) 745-0508 sales@dia-m.ru www.dia-m.ru

С.-Петербург  
+7 (812) 372-6040  
spb@dia-m.ru

Новосибирск  
+7(383) 328-0048  
nsk@dia-m.ru

Воронеж  
+7 (473) 232-4412  
vrn@dia-m.ru

Йошкар-Ола  
+7 (927) 880-3676  
nba@dia-m.ru

Красноярск  
+7(923) 303-0152  
krsk@dia-m.ru

Казань  
+7(843) 210-2080  
kazan@dia-m.ru

Ростов-на-Дону  
+7 (863) 303-5500  
rnd@dia-m.ru

Екатеринбург  
+7 (912) 658-7606  
ekb@dia-m.ru

Кемерово  
+7 (923) 158-6753  
kemerovo@dia-m.ru

Армения  
+7 (094) 01-0173  
armenia@dia-m.ru

